

Catalogue 2017/2018

- // Clasps and Bars
- // Orthodontics
- // CA® CLEAR ALIGNER
- // Digital Orthodontics
- // Dental Sleep Medicine
- // Pressure Moulding Technique
- // Laboratory Equipment

Contents

// Clasps and Bars

Jaw Fracture Splint	19	N BI-FORM Clasp	16	PROFI Bar	19
J Clasp	18	N Clasp	16	SCHEU-Anchor	18
J Clasp oblique	18	NK Clasp	17	T Clasp	17
JM Clasp	18	O Clasp	16	TK Clasp	17
Lingual bar curved	19	OK Clasp	16		

// Orthodontics

3D Model insulation NEW	28	ECHARRI templates	48	Quad Helix	45
Activation key	37	Face bows	42-43	RESIMIX®	28
Activator tubes	44	Frontal guiding bars		Retainer boxes	29
Adams Clasp	22	acc. to R. Hinz	23	SCHEU-Anchor	22
Arrow Clasp	22	Hinge 44	41	Security key	
BIOPERM-trimmer	39	Hinge system acc. to Herbst	41	for expansion screws	37
BLUE-BLOKKER®	40	Key for expansion screws	37	Set-Up technique	38-39
Bottle Spout	26	Lingual/palatal sheaths	45	Silicone tube	28
CARAT® wire bending pliers	46-47	MENZANIUM®	24	Spray bottle / spray nozzles	29
CETRON®	49	Mica	27	STEADY-RESIN	26-27
CHROMIUM	25	Ortho arch bender	47	Symmetry measurement gauge	44
Coffin spring	23	Ortho designs	29	Triangle Clasp	22
Diamond separating disc	40	Orthodontic insulating agent	28	U-Bar acc. to Karwetzky	23
Distance rings	41	POLYFIX	39	VECTOR® assortment NEW	30-36

// CA® CLEAR ALIGNER

CA® Arch CHROMIUM NEW	58	CA® Foil A2 pd NEW	56	CA® Splint case for patients NEW	57
CA® Books NEW	62	CA® Gauge	63	CA® Stripping tools	63
CA® Carbide bur	62	CA® Grinding set	62	CA® TIP pliers	59
CA® CHECKER	54	CA® Photo support 4.0		CA® Travel Case NEW	57
CA® Coloured pencils	61	and accessories	55	ECHARRI PST Set /	
CA® Cyano Veneer Fast	61	CA® Pipettes	61	PST Set Standard	63
CA® Distalizer Kit acc.		CA® Plastic buttons	61	SD box	62
to Dr. Echarrri NEW	57	CA® POWER GRIP	60	VECTOR® 30 NEW	58
CA® FLASH	61	CA® Retention pd	56	VECTOR® 40 acc. to Prof. Hinz	58
CA® Foil NEW	56	CA® SMART 2D	54		

// Digital Orthodontics **NEW**

Asiga MAX NEW	66	OnyxCeph3™ NEW	
IMPRIMO® LC IMPRESSION	67	CA® SMART 3D – laboratory	68-69
IMPRIMO® LC Model	67	OnyxCeph3™ NEW	
IMPRIMO® LC Splint	67	CA® SMART 3D – clinic	68-69

// Dental Sleep Medicine

AM Aligner NEW	79	Positioning aid		TAP®-T	76
CLEAR-BLOKKER®	74	and other accessories	73	TAP® -T Matrix	78
GEORGE GAUGE™		SD Box	75	TAP® -T Reverse	77
Bite Registration	73	SD Cyano Veneer Fast	74	TAP®-T / TAP®-T	
IST® Kit	72	SD Metal Primer	74	Reverse starter set	77
myTAP™	79	SD Pipettes	74	TAP® T XL NEW	
Occlusion trainer	75	TAP® PAP NEW	79	Lower guiding device	78
OPM	75	TAP® Splint case for patients	78	TELESCOPES HF 5 / HF 9 / HF 15	72
		TAP® Standard	76		

// Pressure Moulding Technique


Material		IMPRELON® S pd	89	Indirect bonding technology	92
BIOCRYL® C	82	ISOFOLAN®	91	LC Burs NEW	99
BIOCRYL® Form X NEW	82	OSAMU-Retainer® Kit	89	LC Insulation	98
BIOCRYL® M	82	Skin Packaging	91	LC Primer	98
BIOPLAST®	88	Spacer Foil	91	OSAMU-Bond	98
BIOPLAST® bleach	88	Accessories		OSAMU Polisher	103
BIOPLAST® Color	87	BIOCRYL®-RESIN	96	Pellet brush	101
BIOPLAST® Multicolor	87	BIOPERM Trimmer	103	Pellet Cover	101
BIOPLAST® XTREME	86	Blocking-out putty	95	POLYFIX	103
BIOPLAST® XTREME Deco	87	BIOPLAST® Insulation agent	95	SD Bond	92
BIOPLAST® XTREME Kit	86	BLUE-BLOKKER®	94	SD Bracket pincers	93
BIOPLAST® XTREME PRO Kit	86	CA® Carbide bur	102	SD Bracket remover	93
BRUX CHECKER® NEW		CLEAR-BLOKKER®	97	SD Composite primer	98
acc. to Prof. Sato	88	Cutting bur	102	SD Cutter	93
COPYPLAST®	90	DIMO®	103	SD Felt wheel	99
COPYPLAST® C	90	DIMO®PRO / DIMO®PRO SLIM	103	SD Foil scissors A	100
DURAN® NEW	83	DURASPLINT®	96	SD Foil scissors B	101
DURAN®+ NEW	83	DURASPLINT® LC	97	SD High-gloss polishing mop	99
DURAN®+ A2 pd	84	Finishing bur	102	SD Polishing brush	99
DURAN®+ white pd NEW	84	Finishing Set	102	SD Polishing paste	99
DURASOFT® blue pd	85	Foam discs	100	SD Scalpel	92
DURASOFT® pd	84	HM Carbide cutter	102	SD Separating disc, serrated	93
DURASOFT® seal	85	HM Tricutter fine	102	SD Silicone	92
HARDCAST®	90	Impression tray handles	100	SIL-KITT	94
IMPRELON®	89			Stainless steel pellets	101

// Laboratory Equipment

Asiga MAX NEW	106	IMPRIMO® LC Splint	107	Parallelometer	113
BIOSTAR®	108	LC 6 Light Oven	111	Sampler of pressure	
Heating Support	111	LC Mini Light	111	moulding material	109
Fixators and mounting splints	112	MINISTAR S® NEW	109	TWINSTAR® P	110
IMPRIMO® LC Impression	107	Model holder	113	TWINSTAR® V	110
IMPRIMO® LC Model	107	Model Platform	111		

Appendix:

Indication chart, material lists, opening and expansion parameters of VECTOR® expansion screws, measurement units, numeric table of contents by REF number



Products are changing, values remain.

In 1926, company founder and dentist Rudolf Scheu laid the foundation for SCHEU-DENTAL. We are proud of 91 years of company history, made in Germany. As family firm now in its third generation, we develop dental products and services, characterised by maximum precision. Our attention to detail has helped us to achieve market leadership in the fields of dental pressure

moulding, sleep medicine and 3D printing and to play a key role in shaping the face of dental technology throughout the world. These achievements established worldwide confidence in our innovative ideas for new products.

Our key priority has remained the same: to deliver top quality on a daily basis with the help of a



1980



2000



2010

2016

Watch 91 years of company history:



1990



team of well-trained and highly motivated staff. With the acquisition of smile dental some years ago SCHEU-DENTAL added another orthodontic specialist to its network, as well as CAD/CAM expert in the guise of CA DIGITAL.

our field service. We are looking forward to your feedback – for in particular the insight in daily laboratory routines helps us to focus on the essential and to develop products that make your work easier.

See it for yourself and discover the wide range of products available from SCHEU-DENTAL in our new catalogue. Please feel free to contact our team in Iserlohn or your contact person in

Yours sincerely,
Christian Scheu





50 years BIOSTAR®: a success story

There is hardly a product so closely linked to the history of SCHEU-DENTAL as the BIOSTAR®. From the very beginning in 1967, it has been a benchmark in the pressure moulding technique and up this very day, its functional principle demonstrates its superiority in daily use.

The advancement of the following machine gen-

erations has been consistently pursued – in that same pioneering spirit: Since 1988, a processor unit and dash panel allow to manually program temperatures and heating and cooling times as well as to precisely supervise and control them. The implementation of scanning technology in 2009 was yet another logical step in order to facili-



1967



since 2015

tate and accelerate operation. Thanks to our short-wave, thermostatically-controlled infra-red heater that has been introduced in 2006, the working temperature is reached within only one second.

The BIOSTAR® is a synonym for our innovative strength and has been a part of SCHEU-DENTAL for 50 years now. What makes us particularly proud: many of our clients got along very well


with the BIOSTAR® since the very beginning and are still working with models of the first or second generation. A clear proof of the top quality made in Germany we feel obliged to until today. Therefore our pressure moulding units are delivered with a three-year warranty and a ten-year warranty on availability of spare parts. A guarantee for the success of further generations – at SCHEU-DENTAL and in the practices and laboratories worldwide.

Perfection is a love for detail.

The realisation of complex precision constructions requires dexterity and experience. Therefore, many parts at SCHEU-DENTAL are still manufactured by hand, as this is the only way to ensure the greatest possible quality of the precision mechanical parts. The staff at our site and regular product controls according to the latest certification requirements guarantee that the company never deviates from this quality standard.








Perfection means bringing innovative strength and quality to life.

At our Research and Development department we come up with clever ideas designed to lighten the daily workload of the dental laboratory and practice. Our varied product family combines the latest discoveries with trusty therapeutic products that still demonstrate the innovative strength that inspired them.

From the CA® CLEAR ALIGNER splint system to our BIOSTAR® pressure moulding unit and our 3D printer family IMPRIMO®: we are always on the lookout for improvement and innovative ideas. Whether investigating the enhancement of resins or developing new products for treating snoring and sleep apnoea – our close cooperation with recognised



research institutes and universities gives rise to sophisticated quality products based on cutting-edge technology.

Today, in conjunction with our subsidiary company CA DIGITAL, we can offer you the complete range of digital services and products – with individual solu-

tions for the CA® CLEAR ALIGNER treatment or the digitally planned lingual retainer MEMOTAIN®. As a founding member of Ortho Alliance®, an orthodontic alliance of four leading manufacturers in Germany, we are working on the specific development of digital processes and products in modern orthodontics to the benefit of both, clinician and patient.



Perfection you can rely on.

Whether a query about setting up your SCHEU-DENTAL equipment such as the BIOSTAR® or a question relating to specific applications – our staff in the Sales and Customer Service Departments are on hand with expert advice at all times. We want you to take advantage of their

sound practical guidance from the outset – a service that our customers can rely on whatever the problem. And this goes beyond your actual purchase: for some years now our SCHEU-ACADEMY offers training courses and seminars for dental clinicians and technicians at a regular basis,



from courses on CA®CLEAR ALIGNER and TAP® certification to our pressure moulding workshops and courses on digital orthodontics. Why not get acquainted with the latest developments in the practice and learn how to achieve even more professional results in the future? Discover our wide

range of training courses and seminars for dental practices and laboratories at www.scheu-dental.com/scheu-academy. Or simply agree an appointment with our field sales team right away. We look forward to hearing from you!

Clasps and Bars



For us at SCHEU-DENTAL
continuity and quality
go hand in hand.

We have been manufacturing clasps for more than 90 years. Based on that experience, our range of clasps and bars has become a benchmark. Made of high-grade stainless steel alloys with excellent resistance to oral conditions, fabricated in a cutting and cold forming process, from one piece without soldering connections and with highly-polished surface. Our clasps are of excellent elasticity and stability – your patients can trust you.





O Clasp

- O Clasp (Ø 0.9 mm)*
- O Clasp (Ø 1.0 mm)**

#1002
#1003

* 10/100 pieces. ** 10 pieces

CE




OK Clasp

- OK Clasp (Ø 0.8 mm)

#1007

10 pieces

CE




N Clasp

- N Clasp (Ø 0.9 mm)*
- N Clasp (Ø 1.0 mm)**
- N Clasp (Ø 1.2 mm)**

#1009
#1010
#1011

10/100 pieces, ** 10 pieces

CE




N BI-FORM Clasp

- N BI-FORM-Clasp (Ø 0.9/1.2 mm)

#1012

10 pieces

CE




NK Clasp

○ NK Clasp (Ø 0.8 mm)

#1013

10/100 pieces

€



T Clasp

○ T Clasp (Ø 0.9 mm)


#1026

○ T Clasp (Ø 1.0 mm)

#1027

10 pieces

€



TK Clasp

○ TK Clasp (Ø 0.9 mm)

#1029

○ TK Clasp (Ø 1.0 mm)

#1030


10/100 pieces

€



SCHEU-Anchor

Unique retention and spring force! SCHEU-anchors consist of stainless steel alloys with excellent resistance to corrosion and oral conditions.



- SCHEU-anchor (Ø 0.6 mm)* #2048
- SCHEU-anchor (Ø 0.7 mm) #2049
- SCHEU-anchor (Ø 0.8 mm) #2050
- SCHEU-anchor (Ø 0.9 mm) #2051
- SCHEU-anchor (Ø 1.0 mm)* #2052
- SCHEU-anchor (Ø 1.2 mm)* #2053

10/100/1000 pieces* 10/100 pieces

CE



J Clasp

- J Clasp (Ø 0.8 mm) #1017
- J Clasp (Ø 0.9 mm) #1018
- J Clasp (Ø 1.0 mm) #1019
- J Clasp (Ø 1.1 mm) #1047

10/100 pieces

CE



J Clasp oblique

- J Clasp oblique (Ø 0.8 mm) #1020
- J Clasp oblique (Ø 0.9 mm) #1021
- J Clasp oblique (Ø 1.0 mm) #1022
- J Clasp oblique (Ø 1.1 mm) #1048

10/100 pieces (in pairs 5/50 left and right 5/50)

CE



JM Clasp

- JM Clasp (Ø 1.2 mm) #1023

10 pieces

CE

PROFI Bar

- PROFI bar straight (length 80 mm) #1040
- PROFI bar curved (length 80 mm) #1039

10 pieces

CE





Jaw Fracture Splint

With 20 buttons. MRI incompatible

- Jaw fracture splint, hard (Ø 1.0 mm) #1038
- Jaw fracture splint, soft (Ø 1.0 mm) #1049

10 pieces

CE




Lingual bar curved

Elliptical, hard stainless steel profile, anatomically curved with stamped retentions.

- Lingual bar, curved (length 60 mm) #1051
- Lingual bar, curved (length 70 mm) #1052
- Lingual bar, curved (length 80 mm) #1053

10 pieces

CE




Our products are still unsurpassed!

Quality is not a question of size. End of the 1950ies, our first orthodontic product was the SCHEU-anchor – a retention element, still unsurpassed in terms of quality, workmanship and spring force. Today our orthodontic line is supplemented and enhanced by innovative products such as VECTOR® expansion screws and the CA® CLEAR ALIGNER splint therapy.



Orthodontics




SCHEU-Anchor

Unique retention and spring force! SCHEU-anchors consist of stainless steel alloys with excellent resistance to corrosion and oral conditions.

- SCHEU-anchor (Ø 0.6 mm)* #2048
- SCHEU-anchor (Ø 0.7 mm)* #2049
- SCHEU-anchor (Ø 0.8 mm) #2050
- SCHEU-anchor (Ø 0.9 mm) #2051
- SCHEU-anchor (Ø 1.0 mm)* #2052
- SCHEU-anchor (Ø 1.2 mm)* #2053

10/100/1000 pieces* 10/100 pieces

CE




Arrow Clasp

Made of the nickel-free alloy MENZANIUM® in cold forming, with high-gloss polished surface. Excellent elasticity and breaking strength. Nickel trace elements from 0 to 0.2% are procedurally possible.

- Arrow clasp (Ø 0.7 mm) #2055
- Arrow clasp (Ø 0.8 mm) #2056

10/100 pieces

CE





Triangle Clasp

Made of the nickel-free alloy MENZANIUM® in cold forming, with high-gloss polished surface. Excellent elasticity and breaking strength. Nickel trace elements from 0 to 0.2% are procedurally possible.

- Triangle clasp (Ø 0.7 mm) #2288
- Triangle clasp (Ø 0.8 mm) #2322

100/1000 pieces

CE



Adams Clasp

Modified form for optimised retention and easy adjustment. Made of the nickel-free MENZANIUM® alloy – nickel trace elements of 0 to 0.2% are procedurally possible.

- Adams clasp (Ø 0.7 mm, width 6 mm) #2104
- Adams clasp (Ø 0.7 mm, width 7 mm) #2105
- Adams clasp (Ø 0.7 mm, width 8 mm) #2106
- Adams clasp (Ø 0.7 mm, width 9 mm) #2107
- Adams clasp (Ø 0.7 mm, width 10 mm) #2108
- Adams clasp (Ø 0.7 mm, width 11 mm) #2109

10/100 pieces

CE



Coffin spring

Made of the nickel-free MENZANIUM® alloy – nickel trace elements of 0 to 0.2% are procedurally possible.

○ Coffin spring (Ø 1.2 mm, length 11 mm)

#2419

10/100 pieces

CE



Frontal guiding bars acc. to R. Hinz

For treatment of Angle Class II cases (lower jaw retrusion/ distal occlusion) by means of bimaxillary appliances as well as for bilateral upper jaw expansion in growing children using an expansion screw. Made of the nickel-free alloy MENZANIUM® – nickel trace elements from 0 to 0.2 % are procedurally possible.

○ Guiding bars * (Ø 1.2 mm, length 14 mm)

#2329

○ Guiding bars * (Ø 1.2 mm, length 18 mm)

#2330

* 5/50 pairs

CE

○ Positioning support, 5 pieces

#2331

For insertion and alignment of guiding bars and expansion screws

○ Angle Tester acc. to Adamik

#2332

For fixing and controlling the insertion angle



U-Bar acc. to Karwetzky

Functional orthodontic appliance enabling the lower jaw to move in the required position depending on the type (1-3). U-bars acc. to Karwetzky consist of stainless steel alloys with excellent resistance to corrosion and oral conditions.

○ U-Bar acc. to Karwetzky

#2118

* 5/50 pairs

CE



MENZANIUM®

MENZANIUM® is a nickel-free* stainless steel alloy, fabricated in a patented high pressure melting process replacing allergenic components of nickel by manganese and nitrogen. Compared to conventional stainless steel wires, corrosion resistance and breaking strength are increased. Wire is available in two versions hard (1400-1600 N/mm²) and spring hard (1800-2000 N/mm²).

Nickel trace elements from 0 to 0.2 % are procedurally possible.



MENZANIUM® Coil wire

Our MENZANIUM® coil wires are diamond drawn and high-gloss polished. Processing by joining techniques such as laser-welding and welding. Measurements, tolerances and quality in accordance with DIN EN ISO 15841.

Versions			hard	spring hard
○ 0.5 mm	.020 inch	(50 m)	-	#8459*
○ 0.6 mm	.024 inch	(40 m)	-	#8460
○ 0.7 mm	.028 inch	(30 m)	#8453	#8461
○ 0.8 mm	.032 inch	(20 m)	#8454	#8462
○ 0.9 mm	.036 inch	(10 m)	#8455	#8463
○ 1.0 mm	.040 inch	(10 m)	#8456*	#8464
○ 1.2 mm	.048 inch	(10 m)	-	#8466

*only as laboratory coil

Laboratory coil: 10-50 m wire length clinic coil: 58-330 m (500 g coil) CE


MENZANIUM® Wire rods

Our MENZANIUM® wire rods, elongated and aligned.



diameter		length	hard	spring hard
0.6 mm	.024 inch	75 mm	-	#8478
0.7 mm	.028 inch	75 mm	#8480	#8467
0.7 mm	.028 inch	150 mm	#8482	#8468
0.8 mm	.032 inch	150 mm	-	#8469
0.8 mm	.032 inch	200 mm	-	#8470
0.9 mm	.036 inch	150 mm	-	#8471
0.9 mm	.036 inch	200 mm	-	#8472
1.0 mm	.040 inch	200 mm	-	

all packs 500 g each CE



Breakdown potential for wires from the group of iron alloys

Electrochemical study of corrosion properties of orthodontic auxiliaries in artificial saliva (reference electrode Ag / AgCl, feed rate: 1m V/sec).

M. Neuhofer: "Experimental Dentistry", dissertation 1993, University of Freiburg, head: Prof. Dr. H. F. Kappert

CHROMIUM Coil wire

CHROMIUM coil wires are diamond-drawn and high-gloss polished. Processing by joining techniques such as laser-welding and welding. Measurements, tolerances and quality in accordance with DIN EN ISO 15841. Wire is available in two versions hard (1400-1600 N/mm²) and spring hard (1800-2000 N/mm²).

Versions			hard	spring hard
○ 0.5 mm	.020 inch	50 m	#8323**	#8330
○ 0.6 mm	.024 inch	40 m	#8324**	#8331
○ 0.7 mm	.028 inch	30 m	#8325	#8332
○ 0.8 mm	.032 inch	20 m	#8326	#8333
○ 0.9 mm	.036 inch	10 m	#8327	#8334
○ 1.0 mm	.040 inch	(10 m)	#8328	#8335
○ 1.1 mm	.044 inch	(10 m)	#8329**	#8336
○ 1.2 mm	.048 inch	(10 m)	#8338**	#8339
○ 1.5 mm	.059 inch	(36 m)	#8382*	#8383*
<i>only as clinic coil</i>				
<i>only as laboratory coil</i>				

Laboratory coil: 10-50 m wire length Clinic coil: 26-324 m (500 g coil)

CE



Chromium Wire rods:

Our Chromium wire rods, elongated and aligned.



diameter		length	hard	spring hard
0.7 mm	.028 inch	75 mm	-	#8306
0.7 mm	.028 inch	100 mm	#8293	#8307
0.7 mm	.028 inch	150 mm	#8294	#8308
0.8 mm	.032 inch	100 mm	-	#8310
0.8 mm	.032 inch	150 mm	#8297	#8311
0.8 mm	.032 inch	200 mm	-	#8312
0.9 mm	.036 inch	150 mm	#8299	#8313
0.9 mm	.036 inch	200 mm	-	#8314

all packs 500g each

CE

STEADY-RESIN S

Cold-curing resin for the "salt and pepper" technique. This polymer meets the requirements of the "salt and pepper" technique, such as high stability and optimum wettability.

- *STEADY-RESIN Neon monomers* #8134
- 1 kg €€

STEADY-RESIN M

Cold-curing resin for the "premix" technique. Characterised by a wide processing range and excellent modelling properties.

- *STEADY-RESIN M Polymer* #8135
- 1 kg €€

STEADY-RESIN VARIO

Cold-curing resin for orthodontics with prolonged processing time, can be used both for "salt and pepper" and for the modelling techniques.

- *STEADY-RESIN VARIO Polymer* #8140
- 1 kg €€



STEADY-RESIN S+M

Universal monomer, can be processed with all STEADY-RESIN polymer versions. The clear monomer can be mixed with STEADY-RESIN coloured monomers.

- *STEADY-RESIN S+M (monomer) clear* #8136
 - *STEADY-RESIN S+M (monomer) pink* #8137
- 500 ml €€

Bottle Spout

Reusable

- *Bottle Spout* #8149
- 1 piece



STEADY-RESIN Neon monomers

Ready-to-use coloured neon monomers for use with STEADY-RESIN S, M and VARIO as well as DURASPLINT® polymer. For streak-free and translucent results.

- STEADY-RESIN Neon monomers neon pink #8350
- STEADY-RESIN Neon monomers neon red #8357
- STEADY-RESIN Neon monomers neon yellow #8358
- STEADY-RESIN Neon monomers neon green #8359
- STEADY-RESIN Neon monomers neon blue #8360
- STEADY-RESIN Neon monomers neon orange #8361

100/250 ml

CE

STEADY-RESIN Coloured monomers

Ready-to-use coloured STEADY-RESIN monomers for use with STEADY-RESIN S, M and VARIO as well as DURASPLINT® polymer.

- STEADY-RESIN Coloured monomers blue #8351
- STEADY-RESIN Coloured monomers yellow #8352
- STEADY-RESIN Coloured monomers red #8353
- STEADY-RESIN Coloured monomers black #8354
- STEADY-RESIN Coloured monomers green #8355
- STEADY-RESIN Coloured monomers violet #8356

100/250 ml

CE



STEADY-RESIN Polymer

Pre-mixed polymer in two colours. Innovative pigments provide contrasting effects.

- STEADY-RESIN Polymer black #8130
- STEADY-RESIN Polymer white #8131

200 g

CE



Mica

For individual use with STEADY-RESIN and DURASPLINT® polymer. Blends rainbow, rainbow green and rainbow violet are completely metal-free and ensure iridescent effects in delicate pastel shades.

- Mica silver #8144
- Mica gold #8145
- Mica blue #8146
- Mica red #8147
- Mica green #8148
- Mica rainbow #8150
- Mica rainbow green #8151
- Mica rainbow violet #8152

50 g

CE





RESIMIX®

Original mixing cups of clear silicone, stable by strong suction base. Very practical pouring spout for medium and large version.

- RESIMIX® 5 ml, 3 pieces #7243
- RESIMIX® 30 ml, 2 pieces #7284
- RESIMIX® 70 ml, 1 piece #7244

Orthodontic insulating agent

Alginate insulation for plaster models. Single apply provides optimum insulation.

- Orthodontic insulating agent #8364

500 ml



NEW!



3D Model insulation

Insulation for 3D models when processing cold-curing resin, such as STEADY-RESIN S. Application with enclosed brush provides optimum insulation. 3D Model insulation is suitable for 3D models generated with 3D printers.

- 3D Model insulation #6510

65 ml


Silicone tube

For covering clasp wires.

- Silicone tube #3211
- inner: Ø 1.0 mm, outer: Ø 2.0 mm

1m

CE





Spray bottle / spray nozzles

Spray bottle with dosing tip and cap is suitable for both monomer and polymer. Perfect dosing by piercing the tip with a needle.

Metal spray nozzles available optionally for easier dosing of polymer and monomer. Nozzles can be mounted on spray bottles.

- Spray bottle 100 ml, 1 piece #8141
- Spray nozzles (monomer), 5 pieces #8138
- Spray nozzles (polymer), 5 pieces #8139

Ortho designs

Monomer-resistant film images to insert in orthodontic plates or pressure moulding material.

- Assortment of film images, 400 pieces #8142
- Assortment of film images, 20 pieces #8143

*please indicate desired motif

CE





Retainer boxes

For orthodontic appliances or BIOPLAST® sports mouth guards. Available in blue, yellow, green, red and violet.

- Retainer boxes (5 colours, 2 each) #5435

10 pieces

NEW!**VECTOR® 30**

Special miniaturised sector screw with bright yellow directional arrow for single tooth movements.

To ensure safe anchorage, retention wires (CHROMIUM Ø 0.8 mm, not included) can be individually inserted and adjusted.

Can also be used in the CA® CLEAR ALIGNER technology for oblique-lateral segment expansion acc. to Dr. Pablo Echarri in conjunction with the CA® Arch CHROMIUM and two VECTOR® 30 screws.




○ VECTOR® 30

#2484

W x H x L: 4.8 x 2.6 x 9.5 mm, expansion 3.0 mm

2/10 pieces, incl. 2/10 keys for expansions screws

CE

**VECTOR® 40 acc. to Prof. Hinz**

A special miniaturised screw for fabrication of CA® splints, allowing activation and rapid transversal expansions (up to 3 mm). The retention arms allow safe anchorage on the CA® splints using cold curing resin. With laser-welded directional arrow.



○ VECTOR® 40

#2486

W x H x L: 3.6 x 2.5 x 19 mm, expansion 3.0 mm

10/50 pieces

CE

**VECTOR® 50**

Skeletal miniaturised expansion screw with bright yellow directional arrow for transversal expansion in the lower jaw, also for narrow upper jaws.

○ VECTOR® 50

#2489

W x H x L: 4.8 x 2.5 x 10.0 mm, expansion 3.0 mm

10/50 pieces

CE



VECTOR® 80

Skeletal expansion screw with bright yellow directional arrow for transversal expansion of the lower jaw. Additional bore hole ensures excellent retention in acrylics.

- VECTOR® 80 #2490
W x H x L: 7.0 x 3.1 x 14.3 mm, expansion 7.5 mm

10/50 pieces

CE



**VECTOR® 90 / VECTOR® 100**

Skeletal expansion screw with bright yellow directional arrow for transversal expansion of the lower jaw. For treatment of Angle Class II cases we recommend using VECTOR® 100 in combination with VECTOR® 900.

- VECTOR® 90 #2491
W x H x L: 6.2 x 3.0 x 10.0 mm, expansion 4.0 mm
- VECTOR® 100 #2492
W x H x L: 6.2 x 3.0 x 12.0 mm, expansion 5.0 mm

10/50 pieces

CE


**VECTOR® 140 / VECTOR® 160**

Skeletal expansion screw with bright yellow directional arrow for transversal and sagittal expansion in the upper and lower jaw.

- VECTOR® 140 #2494
W x H x L: 7.5 x 3.3 x 12.0 mm, expansion 5.0 mm
- VECTOR® 160 #2496
W x H x L: 7.5 x 3.3 x 14.3 mm, expansion 6.0 mm

10/50 pieces

CE



**VECTOR® 190 / VECTOR® 200**



Stable, skeletal expansion screw with bright yellow directional arrow for transversal expansion of the upper jaw and bimaxillary appliances. Additional bore hole ensures excellent retention in acrylics.

- VECTOR® 190 #2497
W x H x L: 9.8 x 3.1 x 10.0 mm, expansion 5.0 mm
- VECTOR® 200 #2498
W x H x L: 9.8 x 3.1 x 14.3 mm, expansion 7.5 mm

10/50 pieces

CE







VECTOR® 350

Skeletal miniaturised sector expansion screw with bright yellow directional arrow for single tooth and segment movement in the upper and lower jaw.

○ VECTOR® 350 #2509
W x H x L: 4.8 x 2.5 x 10.0 mm, expansion 4.0 mm

10/50 pieces

CE





VECTOR® 400

Delicate skeletal sector expansion screw with bright yellow directional arrow. Can also be used for single tooth movements.

○ VECTOR® 400 #2510
W x H x L: 6.2 x 3.0 x 11.0 mm, expansion 3.0 mm

10/50 pieces

CE





VECTOR® 420

Delicate skeletal sector expansion screw with bright yellow directional arrow. Can be used for distalisation and for single tooth movements.

○ VECTOR® 420 #2512
W x H x L: 6.0 x 3.0 x 13.0 mm, expansion 5.0 mm

10/50 pieces

CE



VECTOR® 440

Delicate skeletal sector expansion screw with bright yellow directional arrow, allowing distalisation.

○ VECTOR® 440 #2514
W x H x L: 6.0 x 3.0 x 13.3 mm, expansion 5.0 mm

10 pieces

CE



VECTOR® 440 R

Delicate skeletal sector expansion screw with bright yellow directional arrow for single tooth or segment movement.

- VECTOR® 440 R #2513
W x H x L: 6.0 x 3.0 x 13.3 mm, expansion 5.0 mm

10 pieces

CE



**VECTOR® 450**

Skeletal sector expansion screw with limit stop. Prevents the screw from breaking after maximum expansion, reactivation is possible.

- VECTOR® 450 #2515
W x H x L: 7.0 x 2.8 x 13.0 mm, expansion 4.0 mm

10/50 pieces

CE



**VECTOR® 460**

Skeletal sector expansion screw in stronger version with bright yellow directional arrow.

- VECTOR® 460 #2516
W x H x L: 7.5 x 3.3 x 13.0 mm, expansion 5.0 mm

10 pieces

CE



**VECTOR® 500 / VECTOR® 520**



Traction screw for space closure with black directional arrow on the housing.

- VECTOR® 500 #2520
W x H x L: 6.2 x 3.0 x 13 mm, traction 3.0 mm
- VECTOR® 520 #2522
W x H x L: 6.2 x 3.0 x 18.0 mm, traction 5.0 mm

10 pieces

CE








VECTOR® 600* / VECTOR® 620*

Special screw for rapid palatal expansion. The anatomically bent retention arms are soldered or welded onto molar bands. The specially designed housing profile fits perfectly to the palatal vault. The integrated polyamide bush prevents accidental resetting of the screw. Laser-welded marking for direction and expansion on the housing. Comes with security expansion key.

- VECTOR® 600 #2524
W x H x L: 10.4 x 4.0 x 14.0 mm, expansion 9.0 mm
- VECTOR® 620 #2526
W x H x L: 10.4 x 4.0 x 16.0 mm, expansion 11.0 mm

1/10 pieces

CE 0044






VECTOR® 720 acc. to Bertoni

Special Bertoni screw for insertion in multi-sector plates, enabling separate protrusion and symmetric transversal expansion by two independently operating spindles.

- VECTOR® 720 #2530
W x H x L: 17.0 x 3.6 x 15.7 mm
frontal expansion: 5.5 mm
transversal expansion: 7.0 mm

2 pieces

CE



VECTOR® 730 acc. to Bertoni

Special Bertoni screw for insertion in multi-sector plates, enabling separate protrusion and symmetric or asymmetric transversal expansion by three independently operating spindles.

- VECTOR® 730 #2532
W x H x L: 17.0 x 3.6 x 15.7 mm
frontal expansion: 5.5 mm
transversal expansion: 2 x 3.5 mm

2 pieces

CE

*Not suitable for surgically assisted rapid palatal expansion.

Protect screw body from excessive heat that might damage the integrated polyamide bushing and cause jamming.



VECTOR® 740 acc. to Bertoni

Special Bertoni screw for insertion in multi-sector plates, enabling separate protrusion and symmetric or asymmetric transversal expansion by three independently operating spindles.

- VECTOR® 740 #2540
 $W \times H \times L$: 15.0 × 5.0 × 15.0 mm
 frontal expansion: 2.5 mm
 transversal expansion: 5 mm

2 pieces

CE



**VECTOR® 750 acc. to Bertoni**

Special Bertoni screw in modified, angled and delicate version for insertion in multi-sector plates, enabling separate protrusion and symmetric or asymmetric transversal expansion by three independently operating spindles. For bilateral expansion.

- VECTOR® 750 #2542
 $W \times H \times L$: 15.0 × 5.0 × 15.0 mm
 frontal expansion: 2.5 mm
 transversal expansion: 2 x 2.5 mm

2 pieces

CE


**VECTOR® 800* / VECTOR® 820***

Special miniaturised screw for distraction of the lower jaw using stable retention arms of 60 mm length. With laser marking on housing.


- VECTOR® 800 #2533
 $W \times H \times L$: 7.0 × 3.3 × 16.5 mm, expansion 9.0 mm
- VECTOR® 820 #2534
 $W \times H \times L$: 7.0 × 3.3 × 18.5 mm, expansion 11.0 mm

1/5 pieces with security key for expansion screws

CE 0044



*Not suitable for surgically assisted rapid palatal expansion.




VECTOR® 900

Screw for protrusive double plate for treatment of Angle Class II cases (lower jaw retrusion/distal occlusion) as well as for bilateral expansion of the upper jaw in growing children. Compact design with integrated expansion screw, comes with integrated key for expansion screws and SD box. For lower jaw we recommend the combined use with VECTOR® 100.

- VECTOR® 900 #2536
W x H x L: 14.0 x 20.0 x 13.0 mm, expansion 5.0 mm
- Angle tester acc. to Adamik #2332
For fixing and controlling the insertion angle

1 piece

€




VECTOR® Piston screw assortment

Stainless steel screw for single tooth movement during orthodontic treatment with removable appliances. Three different screw lengths allow large protrusion by a constant and precisely controlled amount of pressure on the tooth.

- Piston screw 4 mm, 5 pieces #2477
- Piston screw 6 mm, 5 pieces #2478
- Piston screw 8 mm, 5 pieces #2479
- Serrated housing nut, 5 pieces #2481

€



Accessories

- Mounting screw, 5 pieces #2480
- Reamer, 1 piece #2482
- Screwdriver F, 1 piece #2483
- Screwdriver M, 1 piece #2083



Key for expansion screws

For precision adjustment of VECTOR® expansion screws.

- Key for expansion screws #2101
- Key for expansion screws, angled * #2130

10/100pieces, * 10 pieces

€€



Security key for expansion screws

Suitable for HF 5, HF 9 and HF 15 telescopes of IST® devices and for all types of rapid palatal screws.

- Security key for expansion screws #2100

10 pieces

€€

Activation key:

Activation key for TAP®-T and TAP®-T Reverse protrusion adjustment device. Can also be used for IST® hexagon screws.

- Activation key #9035

1 piece



Screwdriver F*

To activate VECTOR® piston screw.

- Screwdriver F #2483

1 piece



Screwdriver M*

For VECTOR® mounting screw.

- Screwdriver M #2083

1 piece

* The screwdriver is not a medical device and is not suitable for use in the mouth.



Set-Up plaster

White formaldehyde-free plaster offers best conditions for use in general orthodontics and the CA® CLEAR ALIGNER technology.

Technical data according to DIN EN ISO 6873:
Compressive strength after 24 hours 65 MPa
Expansion after 2 hours < 0.09%

○ Set-Up plaster #5266
2.5/5.0 kgs

Set-Up separating disc

With stable mesh layer, comes with mandrel with reinforced shank (Ø 40 mm).

Application

Hard plaster, ceramics, precious metal, NEM



○ Set-Up separating disc with mandrel, 0.5 mm #5293
1 piece



Set-Up trimmer

Special carbide bur for drilling and vertical cuts

○ Set-Up trimmer #5371





Set-Up manual saw

The angled bar allows precise and controlled separation of the individual tooth segments.

- *Set-Up manual saw* #5372
- *Set-Up saw blades** #5373
(W x H x L: 0.15 x 2.0 x 75.0 mm)
- *Set-Up saw blades** #5379
(W x H x L: 0.20 x 2.0 x 75.0 mm)

10 pieces

Set-Up wax blocks

Fast-melting wax for transpositions in orthodontics.

- *Set-Up wax blocks, red* #5368

3 blocks of 50 g



Set-Up wax rods

Soft Set-Up wax with good adhesion properties.

- *Set-Up wax rods, opaque* #5268

25 rods of approx. 1.2 g



BIOPERM-trimmer / POLYFIX

For grinding and finishing soft-elastic and sandwich material. BIOPERM-trimmer in blue, POLYFIX in brown.

- *BIOPERM trimmer, 2 pieces* #3226
- *POLYFIX on mandrel, 1 piece* #3371
- *POLYFIX without mandrel, 6 pieces* #3372



Diamond separating disc

Diamond coated on both sides (thickness 0.28 mm) for very precise saw cuts on set-up models and for finishing splints.

- Diamond separating disc mounted on mandrel, Ø 22 mm

#5288



BLUE-BLOKKER®

Light-curing, blue-transparent blocking-out and modelling putty with optimum flow properties. For fabrication of space holders for bleaching splints and medical trays as well as for sealing and blocking Set-Up models. Burns without leaving any residue, can therefore also be used in prosthetics for the fabrication of fixed and removable prostheses. Telescopic models are stable and can be withdrawn without tension. Curing in the wavelength range of 380-470 nm.

- BLUE-BLOKKER® 4 x 1.5 g, 2 cannula (Ø 0.6) and 8 cannula (Ø 0.9) #5276
- cannula* (Ø 0.4 mm) #5277
- cannula* (Ø 0.6 mm) #5264
- cannula* (Ø 0.9 mm) #5278

*5 pieces

Hinge system acc. to Herbst*

For treatment of class II cases, for soldering or welding on cast or ready-made bands. Without limit stop.



- Set: 1 pair of hinges and 4 pairs of fastening elements #5479
incl. hexagon socket screws
- 4 pairs of fastening elements #5455
incl. hexagon socket screws
- Set 1: pair of hinges and 4 pairs of fastening elements #5480
incl. slotted screws
- 4 pairs of fastening elements #5456
incl. slotted screws
- 1 pair of hinges #5478

© 0044

Distance rings

When fabricating Herbst* appliances, distance rings can be slid onto connecting rods, enabling progressive re-activation.



- 1 mm #5481
- 2 mm #5482
- 3 mm #5483
- 4 mm #5484
- 5 mm #5485

10 pieces, inner diameter 2.0 mm

© 0044

* Herbst is a registered trademark of Dentaurem GmbH & Co KG, Ispringen.



Hinge 44 with slotted screw

For treatment of class II cases with removable or cemented splints. Laser-welded retention (44 mm) ensures secure anchorage in the splints. Without limit stop.

- 1 pair of hinges and 4 pairs of fastening elements #5490
- 4 pairs of fastening elements #5489

© 0044





ORTHOFIT® Face bows

Face bow made of stainless steel. Inner and outer bow (Ø 1.8 mm) double-sided laser welded. Permanent centre marking and coloured rings to identify the size. Inner bow made of one piece. Outer bow available in short and standard with bent outer loops and standard inner bow or U-loops (Ø 1.15 mm, .045") in five sizes. The bent outer loops with 16 degree angle provide optimal power transmission.

Wire ends of outer loops are laser deburred.

Size	Rod length*	Height inner bow
○ Standard	120 mm	
○ Size 1	83 mm	26.0 mm
○ Size 2	90 mm	28.0 mm
○ Size 3	97 mm	30.5 mm
○ Size 4	104 mm	35.0 mm
○ Size 5	111 mm	39.5 mm


* rod length = stretched length of the inner bow with U-loops.

ORTHOFIT®-Line

Inner bow with pre-bent U-loops and outer bow in five sizes. Same length of outer and inner bow for optimised power projection. Marked with coloured ring.

○ Size 1	red	#2453
○ Size 2	green	#2454
○ Size 3	yellow	#2455
○ Size 4	blue	#2456
○ Size 5	white	#2457

1 pieces CE



ORTHOFIT®-Anatomic

Standard inner bow with two outer bows of different lengths. Marked with coloured ring.

○ Anatomic short	red	#2450
○ Anatomic standard	blue	#2451



1 piece CE



Face bows

Face bow made of high-gloss polished stainless steel. Inner and outer bow (Ø 2.0 mm) double-sided laser welded. Permanent centre marking and coloured caps to identify the size. Inner bow made of one piece. Available with regular or short outer bow with bent outer loops as well as standard inner bow or U loops in five sizes

Wire ends of outer loops are laser deburred.



Size	Rod length*	Height inner bow
○ Standard	120 mm	
○ Size 1	84 mm	26.0 mm
○ Size 2	89 mm	28.0 mm
○ Size 3	94 mm	30.5 mm
○ Size 4	103 mm	35.0 mm
○ Size 5	112 mm	39.5 mm

* rod length = stretched length of the inner bow with U-loops.



Face Bow, high-gloss polished

Stainless steel version with high-gloss polish, stable double-sided laser welding. One-piece inner bow and regular or short outer bow with coloured protection caps. Sizes 1-5 are fitted with loops of 15 ° torque.

○ Standard, short	white	#2439
○ Standard, regular	white	#2119
○ Size 1, regular	yellow	#2334
○ Size 2, regular	blue	#2335
○ Size 3, regular	black	#2336
○ Size 4, regular	green	#2337
○ Size 5, regular	red	#2338

1 piece

CE



Activator tubes acc. to Teuscher

Polymerisable fixation for face-bows.

- Activator tubes #2436
Inner diameter: 1.15 mm/.045"
- Activator tubes #2434
Inner diameter: 1.30 mm/.051"

10 pieces

€




Activator tubes

For insertion of face bows in individually curved Adams clasps, double-sided stable laser welding.

- Activator tubes (length 6.0 mm) #2472
Inner diameter: 1.15 mm/.045"
Wire: 1.0 mm h / 125 mm length
- Activator tubes (length 6.0 mm) #2471
Inner diameter: 1.15 mm/.045"
Wire: 0.8 mm sh/125 mm length
- Activator tubes (length 4.5 mm) #2470
Inner diameter: 1.15 mm/.045"
Wire: 0.7 mm sh/125 mm length

10 pieces

€




Symmetry measurement gauge

6 mm acrylic glass plate with mounting pins for exact positioning on the raphe palati. 2 mm grid allows precise determination of sagittal and transversal asymmetries of dental arc length and width.

- Symmetry measurement gauge #8431

1 piece



Quad Helix

Ready-made expansion arch from CHROMIUM wire (Ø 0.9 mm) for fixation at the upper molars in palatal sheaths.

Available in straight and curved.


- Quad Helix Size 1, straight sheath #8010
- Quad Helix size 1, curved sheath #8011
- Quad Helix size 2, straight sheath #8020
- Quad Helix size 2, curved sheath #8021
- Quad Helix size 3, straight sheath #8030
- Quad Helix size 3, curved sheath #8031

5 pieces

© 0044

	straight			curved		
	8010 Size 1	8020 Size 2	8030 Size 3	8011 Size 1	8021 Size 2	8031 Size 3
Length 1	13	16	17.5	13	16	17.5
Length 2	25	27	29	25	27	29
Width	36	39	45	36	39	45
Height	9	9	9	10	10	10
S	5.5-5					

S = dimension sheath



Lingual/palatal sheaths

Stainless steel sheaths for soldering or laser-welding on molar bands. Available in straight and curved.

- Lingual/palatal sheath, straight #8080
- Lingual/palatal sheath, curved #8081

5 pieces

© 0044



CARAT® Wire bending pliers

Pliers made of high quality stainless steel alloy with exactly machined tool tips and precisely guided sheaths. Ergonomic design with rounded edges for fatigue-free working.


CARAT® Crampon pliers

Pin and wire cutting and bending pliers with chequered and crenate tips, suitable for hard (\varnothing up to 1.2 mm) and soft wires (\varnothing up to 1.5 mm).

○ CARAT® Crampon pliers

#2040

1 piece


**CARAT® Flat Pliers**

Flat needle-nose pliers according to Marburg model with delicate chequered guiding tips.

○ CARAT® Flat pliers

#2041

1 piece


**CARAT® Young pliers**

Loop and arch bending pliers with multiple crenate tips for labial arches and hard wire (\varnothing up to 0.9 mm).

○ CARAT® Young pliers

#2042

1 piece



CARAT® Angle pliers (HM)

Carbide wire bending pliers with a square and a round tip (tungsten carbide), for wires up to 0.7 mm diameter.

○ CARAT® Angle pliers (HM)

#2043

1 piece




CARAT® Aderer pliers large

Aderer wire bending pliers, three prong with reinforced plier tips, for wires up to 0.9 mm diameter.

○ CARAT® Aderer pliers large

#2044

1 piece




CARAT® Aderer pliers medium

Aderer wire bending pliers, three prong with reinforced plier tips, for wires up to 0.7 mm diameter.

○ CARAT® Aderer pliers medium

#2046

1 piece



Ortho arch bender


Bending aid of aluminium with 11 different radii for bending round wire (Ø 0.5 -1.5 mm) and square wires (□ up to 0.8 mm / 0.31").

○ Ortho arch bender

#8450

1 piece





ECHARRI Ceph template

Blue template, used for cephalometric measurements, plaster cast study, clinical examination in the CA® CLEAR ALIGNER treatment as well as in other orthodontic treatments.


○ ECHARRI Ceph template #8432
1 piece



ECHARRI Arch template

Red template, used to measure the length of the dental arch, visualising the deviation from the ideal arc by overlay. Can also be used in the fixed orthodontic treatment for adapting arches.


○ ECHARRI Arch template #8433
1 piece



ECHARRI VTO template

Green template, used for orthodontic or surgical VTO (Visualised Treatment Objective).

○ ECHARRI VTO template #8434
1 piece



ECHARRI Template kit

Comes with ECHARRI Ceph template, ECHARRI Arch template and ECHARRI VTO templates well as detailed instructions.

○ ECHARRI template kit #8435
1 piece

CETRON® Cleansing powder

Mild powder for long lasting care and cleansing of orthodontic plates, splints, dentures and sports mouth guards. To dissolve in 150 ml of water.

○ CETRON® cleansing powder #3239

5/25/250 bags



CETRON® Cleansing spray

Care and cleansing of splints, dentures, orthodontic plates, mouth guard and snoring appliances. Fresh lemon taste reduces typical monomer taste, ensuring highly polished appliances and high wearing comfort at the same time.

○ CETRON® Cleansing spray, 75 ml #3238

75 ml

CE



CETRON® Care set

Care and cleansing of splints, dentures, orthodontic plates, mouth guard and snoring appliances.

○ CETRON® care set #3236
incl. cleansing box, brush and 5 bags CETRON®
cleansing powder of 15 g

1 piece


CA® CLEAR ALIGNER

The right concept
for every demand.

Many things are different with CA® CLEAR ALIGNER, the form of the splint, the material thickness... But the most essential difference: It's up to you to decide – holding the complete value creation in your hands or calling for the services of CA DIGITAL. CA® CLEAR ALIGNER leaves you all the options.

More than 200,000 patients have been successfully treated with CA® CLEAR ALIGNER – a significant number and above all a wide diversity of experience.

www.ca-clear-aligner.com
www.ca-digit.com




One system – many options.





The complete value creation in your hands. The individual components of the CA® CLEAR ALIGNER system give you the opportunity to carry out the treatment from diagnosis to therapy planning, fabrication of set-up models until the final splint fabrication and activation in your practice or own laboratory. The necessary equipment consists of the pressure moulding devices MINISTAR S® or BIOSTAR®, the CA® foils in soft, medium and hard, the software CA® SMART 2D including web cam, CA® photo support and CA® CHECKER.



Interested in getting started with digital splint therapy? You'll find your ideal partner with CA DIGITAL, a service provider offering you the complete range of digital treatment solutions. Whether you prefer receiving printed models of digitally planned set-ups in order to carry out the fabrication of the splints yourself or if you want to order digital planning including splint fabrication: CA DIGITAL is your individual partner for all your requirements.



The 3-Splint-System




Three material thicknesses worn during four weeks for a unique and beautiful smile: during each treatment step, the patient wears three splints of various thicknesses. Thanks to the steady increase in pressure, teeth are efficiently, yet smoothly aligned – with extraordinary comfort. The special design with additional covering of the gingival margin (2-3 mm) increases power transmission and ensures even more efficiency in treatment.

CA® CLEAR ALIGNER: a global success story.

The CA® CLEAR ALIGNER therapy system is mentored by an international board of experienced specialists from science and technology, contributing their expertise in orthodontics, dental technology and material science. This expert panel chaired by Dr. Pablo

Echarri is in constant contact to share new findings in seminars, workshops and technical articles with CA® users.

- more than 2,500 certified orthodontists
- more than 200 certified laboratories
- 17 years of clinical experience
- more than 200,000 CA® treatment cases



Extensive support in every respect.

Professional


Actual case studies of our international specialists are regularly published as CA® Clinical Paper at www.ca-clear-aligner.com/en/b2b/indikationen. New findings in treatment planning are recorded in clinical protocols and shared with our certified CA® partners.

Personal

Our CA® partners are consistently encouraged by direct phone or personal support as well as by practical training in our SCHEU-ACADEMY.

Communicative

We provide free images and videoclips for your website as well as patient information for your practice.





New version!

CA® SMART 2D

The enhancements made to the CA® SMART 2D software in version 4.0 offer new possibilities for tooth movements in the CA® set-up technique. Thanks to the modified presentation of the different views it is now possible to take lateral and vestibular photos in addition to frontal and occlusal photos. Single or multi-user installation, can be carried out by your local system administrator. Comes with CA® web cam.

System requirements

From Windows 7 Service Pack 1, USB 2.0 port, OpenGL Version ≥ 2.1 and DirectX 9.0c, monitor with resolution of min. 1024 px \times 768 px, mouse with scroll wheel

○ CA® SMART 2D

#5311

1 piece



CA® CHECKER

Set consisting of digital precision measuring device with USB port, solid base plate made of polished granite, model holder, spherical bur and level. Precision measurement device equipped with capacitive inclination sensors for determination of torque and angulation values before and after transpositions on set-up models.

System requirements

From Windows 7 Service Pack 1, USB 2.0 port, OpenGL Version ≥ 2.0 and DirectX 9.0c, monitor with resolution of min. 1024 px \times 768 px

○ CA® CHECKER

#5274

1 piece



CA® Photo support 4.0

The stable height adjustable support for the CA® web cam enables reproducible documentation of occlusal and frontal photos. Technical precondition and separate accessory for the CA® SMART 2D software. CA® photo support comes with calibration plate. CA® web cam is not part of the CA® photo support, but part of the CA® SMART 2D software.

Dimensions (W × H × D): 180 × 210 × 195 mm

○ CA® photo support 4.0

#5315

1 piece

Extension set (4.0) for CA® photo support

The extension set for the CA® SMART 2D 4.0 software allows taking lateral and vestibular photos in addition to occlusal and frontal ones.

Contents: 1 CA® model holder, 1 hexagon key, 2 height adjustment blocks, 1 CA® SMART Matrix, pins, sleeves, labels identifying the height adjustment blocks

○ Extension set (4.0) for CA® photo support

#5313

○ CA® model holder for photo support 4.0

#5361

For photo support with CA® SMART 2D software 4.0
For occlusal, frontal and lateral photos

1 piece



Height adjustment block for CA® photo support

For fixing the exact height of a model, used with the CA® photo support, removable; remains at the model during processing, enables parallel work on several patient cases.

Contents: Height adjustment block, labels identifying the height adjustment blocks

○ Height adjustment block for CA® photo support

#5314

1 piece



NEW: CA® Foil extra soft

CA® Foil extra soft / soft / medium / hard*

High-quality, transparent, hard-elastic material for gentle tooth alignment with aesthetic splints in the CA® system. Material with original CA® labelling stating the thickness: extra soft (0.4 mm), soft (0.5 mm), medium (0.625 mm) and hard (0.75 mm).

- CA® foil extra soft* #3441
- CA® foil soft #3405
- CA® foil medium #3406
- CA® foil hard #3407

10/100 pieces *10 pieces

CE



CA® Retention pd*

CA® Retention foils for long-term use, pre-dried and individually vacuum sealed. Material with original CA® labelling stating the thickness (0.75 and 1.0 mm).

- CA® Retention pd (0.75 mm) #3438
- CA® Retention pd (1.0 mm) #3440

10 pieces

CE



NEW!

CA® Foil A2 soft / medium / hard pd*

High-quality, hard-elastic, tooth-coloured (A2) material for gentle tooth alignment with aesthetic splints in the CA® system. Material pre-dried and individually vacuum-sealed.

- CA® Foil A2 soft pd #3453
- CA® Foil A2 medium pd #3454
- CA® Foil A2 hard pd #3455

10 pieces

CE

*Delivery exclusively to certified CA® CLEAR ALIGNER partners.

NEW: Splint remover



CA® Splint case for patients*

Handy zipped case for storing and carrying CA® splints. Re-sealable CA® mini grip bags for passing on the CA® splints to the patient, can be labelled. Splint remover for easy removal of pressure moulded splints.

- CA® Set splint case for patients* #5282
comes with CA® Travel Case, 3 CETRON® cleansing powder bags of 15 g, 3 mini grip bags, Splint remover
- CA® mini grip bags* #5283
- Splint remover, 5 pieces #5359

NEW!



CA® Travel Case*

Retainer box for up to two CA® CLEAR ALIGNER splints. Specific shape to fit easily into purse or luggage. Detachable lid with additional space for expansion key, two ventilation holes.

- CA® Travel Case #5303

NEW!



CA® Distalizer Kit acc. to Dr. Echarri

Kit for the fabrication of a CA® Distalizer acc. to Dr. Echarri for one and double-sided distalisation of molars.




All components available separately, see instructions for article numbers.

- CA® Distalizer Kit acc. to Dr. Echarri #5319
1 CA® pressure spring, 1 CA® guiding rod, 1 CA® telescope sleeve with silicone tube, 1 VECTOR® 350 with key for expansion screws

CE



NEW!



VECTOR® 30

Special miniaturised sector screw with bright yellow directional arrow for single tooth movements. To ensure safe anchorage, retention wires (CHROMIUM Ø 0.8 mm, not included) can be individually inserted and adjusted.

For use in the CA® CLEAR ALIGNER technology for oblique-lateral segment expansion acc. to Dr. Pablo Echarri in conjunction with the CA® Arch CHROMIUM and two VECTOR® 30 screws.

○ VECTOR® 30 #2484
W x H x L: 4.8 x 2.6 x 9.5 mm, expansion 3.0 mm

2/10 pieces, incl. 2/10 keys for expansions screws.

CE

NEW!



CA® Arch CHROMIUM

Pre-bent arch made of CHROMIUM wire as auxiliary element. Tapered ends facilitate sliding-on the screws.

Application

For the fabrication of special CA® CLEAR ALIGNER splints with two built-in VECTOR® 30 sector screws for lateral-oblique segment expansion acc. to Dr. Pablo Echarri.

CHROMIUM wire diameter: 0.8 mm

○ CA® Arch CHROMIUM #5317

5 pieces

CE



VECTOR® 40 acc. to Prof. Hinz

Special miniaturised screw for fabrication of CA® splints, allowing activation and rapid transversal expansions (up to 3 mm). The retention arms allow safe anchorage on the CA® splints using cold curing resin. With laser-welded directional arrow.

○ VECTOR® 40 #2486
W x H x L: 3.6 x 2.5 x 19 mm, expansion 3.0 mm

10/50 pieces

CE

CA® TIP 1

Large spherical tip allows insertion of retention points in CA® splints, for increase of friction. Stainless steel, sterilisable.


○ CA® TIP 1 (hemisphere Ø 1.4 mm) #5281
 CE



CA® TIP 2

Small spherical tip allows forming of activation points in CA® splints. Stainless steel, sterilisable.



○ CA® TIP 2 (hemisphere Ø 1,0 mm) #5291
 CE



CA® TIP 3

Pliers with oval plateau, enable creation of activation surfaces in CA® splints. Stainless steel, sterilisable.

○ CA® TIP 3 #5292
 CE



CA® TIP 4

Pliers for torque control with curved plateau, enable creation of labial /lingual activation surfaces in CA® splints. Stainless steel, sterilisable.

○ CA® TIP 4 #5295
 CE

CA® POWER GRIP

CA® POWER GRIP silicone forms enable the production of individual attachments for the CA® CLEAR ALIGNER technology. The set contains nine different moulds that can be filled with light-curing composite for direct application on the prepared tooth.

CA® POWER GRIP MFM ①

MFM (multifunctional forced movement): is used for extrusion and mesial/distal rotation. For use with CA® splints and elastics.

CA® POWER GRIP MDI ②

MDI (mesial/distal inclination): enables setting upright frontal teeth and for mesial/distal inclination.

CA® POWER GRIP DMI ③

DMI (distal mesial inclination): enables setting upright frontal teeth and for distal/mesial tooth inclinations.

CA® POWER GRIP FALA ④

FALA (force application LABIAL): used for protrusion, retrusion and torque control up to 3°. For creation of activation surfaces on the teeth in the gingival area with labial placement – for teeth: 11, 13, 21, 23, 33 and 43.

CA® POWER GRIP FALI ⑤

FALI (force application LINGUAL): used for protrusion, retrusion and torque control up to 3°. For creation of activation surfaces on the teeth in the gingival area with labial, lingual and palatal placement – for teeth: 12, 22, 31, 32, 41 and 42.

CA® POWER GRIP FALI ⑥

BUB (Build-Up bicuspid): round pad for molars (3 mm) used for occlusal bite elevation (approx. 2 mm).

CA® POWER GRIP FALI ⑦


BUM (Build-Up molar): oval molar pad (3 × 5 mm) used for occlusal bite elevation (approx 2 mm).

CA® POWER GRIP FALI ⑧

BUAS (Build-Up anterior small, width 2 mm): used as intrusion plateau for the treatment of deep bite. Helps maintaining intrusion forces while the patient is not wearing the CA® CLEAR ALIGNER, allows molar extrusion.

CA® POWER GRIP FALI ⑨

BUAL (Build-Up anterior large, width 4 mm) also serves as intrusion plateau for the treatment of deep bite, width adjusted to front teeth with larger palatal surface.



- CA® POWER GRIP comes with 2 holders and 18 silicone moulds (9 versions, 2 pieces each), instructions and CD. #5338
- CA® POWER GRIP MFM ①* #5330
- CA® POWER GRIP MDI ②* #5331
- CA® POWER GRIP DMI ③* #5332
- CA® POWER GRIP FALA ④* #5333
- CA® POWER GRIP FALI ⑤* #5334
- CA® POWER GRIP BUB ⑥* #5335
- CA® POWER GRIP BUM ⑦* #5336
- CA® POWER GRIP BUAS ⑧* #5337
- CA® POWER GRIP BUAL ⑨* #5339
- CA® POWER GRIP holder one holder each for frontal and posterior region, stainless steel, sterilisable. #5340

*2 pieces

CE

CA® FLASH

For insulation and easy removal when fabricating pressure moulded CA® splints. CA® Flash is applied to ISOFLAN® foil prior to pressure moulding.

- CA® FLASH #5308

100 ml pump spray bottle, without propellant gas



CA® Cyano Veneer Fast / CA® Pipettes

Crystal-clear, fast-curing bonding agent for adhesion of CA® plastic buttons on CA® splints, can also be used for secure adhesion of screws for Herbst and IST® appliances. CA® pipettes for precise dosing of adhesive.

- CA® Cyano Veneer Fast, 5 g #5302

€

- CA® pipettes, 25 pieces #5300



CA® Plastic buttons

Transparent plastic buttons allow placing elastics in the treatment of class II and III.

- CA® plastic buttons #5301

4/40 pieces

€

CA® Coloured pencils

Coloured pencils for marking the CA® measuring points on the hard plaster models, with red and blue lead.

- CA® coloured pencils #5290

2 pieces





CA® Grinding set / CA® Carbide bur

Grinding set consists of six different instruments and dressing stone. CA® carbide bur made of fine grit carbide characterised by high cutting performance, vibration-free concentric running and very long endurance.

- CA® Grinding set #5263
Comes with set-up and diamond separating disc, DIMO® PRO with mandrel, HM-tricutter (fine), set-up trimmer, dressing stone and CA® carbide bur.
- CA® Carbide bur #5262

1 piece

CA® Book volume 1

CA® Book by Dr. Pablo Echarri on CA® splint therapy with detailed description of the CA® CLEAR ALIGNER treatment concept as well as diagnosis, indication and fabrication.

More than 430 pages with many coloured illustrations.

- CA® Book volume 1 German #5298
- CA® Book volume 1 English #5305
- CA® Book volume 1 Spanish #5306

1 piece



CA® Book volume 2 Advanced

Complements volume 1, for advanced users. With advanced therapeutic options of the CA® CLEAR ALIGNER system and survey of the state-of-the-art in CA® CLEAR ALIGNER technology, such as Progressive Stripping Therapy (PST), CA® POWER GRIP moulds No.1-9, VECTOR® 30, CA® Distalizer Kit acc. to Dr. Echarri, digital technology, extraction cases and combined therapy.

More than 250 pages with many coloured illustrations.

- CA® Book volume 2 German #5307
- CA® Book volume 2 English #5309
- CA® Book volume 2 Spanish #5318

1 piece

NEW!



SD box

Retainer box for storing various orthodontic and snoring appliances, such as positioner, IST® device, TAP® / TAP®-T splint, with mirror and safe closure. Air slits provide optimal ventilation.

Dimensions (LxWxH) 95 x 70 x 31 mm

- SD Box red #5437
- SD Box green #5438
- SD Box blue #5439

10 pieces



CA® Gauge

Set of stainless steel for control of enamel reduction.
Autoclavable up to max. 135° Celsius.



○ CA® gauge (0.1 - 0.5 mm)

#5316

1 piece

CE

CA® Stripping tools

Instruments for interproximal enamel reduction (IER), diamond coated on one or both sides, available in 12 versions. Coating on left, right or both sides. Additional CA® stripping tool (15 µm) with extremely fine and smooth diamond coating on both sides available for polishing.

For safe use of CA® stripping tools without EVA hand piece sterilisable stainless steel holders (in straight or angled version) are best suited.



Diamond particle size	Colour	diamond coating left	diamond coating right	double-sided diamond coating
● 90 µm*	green	#5342	#5343	#5344
● 60 µm	grey	#5345	#5346	#5347
● 40 µm	red	#5348	#5349	#5350
○ 25 µm	white	#5351	#5352	#5353
● 15 µm	yellow	-	-	#5354
○ CA® stainless steel holder, straight				#5355
○ CA® stainless steel holder, angled				#5356

CE

ECHARRI PST Set / PST Set Standard

Set acc. to Dr. Pablo Echarri for use in the progressive stripping therapy (PST) containing CA® stripping tools with double-sided diamond coating in five resp. three diamond particle sizes. Ensures a very smooth enamel finish using an EVA handle. Autoclavable.



- ECHARRI PST Set #5341
comes in aluminium case containing 1 × CA® stripping tool in 90 µm* (green), 60 µm (grey), 40 µm (red), 25 µm (white) and 15 µm (yellow), double-sided diamond coating and CA® gauge
- PST set standard #5357
comes in aluminium case containing 1 × CA® stripping tool in 40 µm (red), 25 µm (white) and 15 µm (yellow), double-sided diamond coating


CE

* not to be used for enamel removal, only for reduction of fillings and crowns.

Digital Orthodontics

Having the future in sight.

Clinical and dental processes in orthodontics are undergoing radical change. Actually, a striking shift of paradigms is taking place: Terms like scanning, 3D planning and printing of working models, splints or individual functional trays are nowadays most commonly used in the orthodontic language. We are among the pioneers of this process of digitalisation and are looking forward to the promising developments of the future. It is our ambition to bring this future to life in your practice or laboratory and to assist you in implementing new technologies and efficient working processes.





NEW!

Asiga MAX

The IMPRIMO® system now includes the Asiga MAX 3D printer, offering intuitive operation with fast and precise printing results. The built-in LCD panel with touch screen function provides real-time information on the building process and enables interactive control. Multiple print jobs can be stored on the internal solid-state drive and managed via touch screen or remotely via web browser. Printing jobs can be transferred over Ethernet or Wifi. In the absence of networking infrastructure, the Asiga MAX can operate as a wireless hot-spot for an independent production solution. Thanks to the IMPRIMO® 385 nm material range and reliable DLP technology, the Asiga MAX can be used in all fields of dentistry. Printing jobs can be created using the Asiga Composer software, included license-free with every system.

- *Asiga MAX* #6501
- *Tray* #6505



IMPRIMO® LC Model

Light-curing resin at the base of methacrylate for fabrication of high-precision dental models with smooth surfaces. Compatible with 3D printing systems with 385 nm light source.

○ IMPRIMO® LC Model

#6502



IMPRIMO® LC IMPRESSION

Light-curing resin at the base of methacrylate for fabrication of individual impression trays and bases for bite registrations, neon orange. Compatible with 3D printing systems with 385 nm light source.

○ IMPRIMO® LC IMPRESSION

#6506




IMPRIMO® LC Splint

Light-curing resin at the base of methacrylate, for fabrication of high-precision transparent occlusal splints and drilling templates. When cured, the resin is characterised by relative flexibility and provides high wearing comfort, reduced fracture susceptibility and excellent polishing properties. Contains no diluting monomers, compatible with 3D printing systems with 385 nm light source.

○ IMPRIMO® LC Splint

#6503



NEW!

OnyxCeph3™ CA® SMART 3D

The independent software solution for digital planning of aligners and virtual bracket positioning with numerous interfaces to the practice software. An automatic control of transpositions during set-up is defined in the CA® CLEAR ALIGNER treatment philosophy and correspondingly reflected in the software. The control of the complete treatment process and value creation during all involved treatment steps remain in the hands of the user. Available as clinical or laboratory version as license for up to five workstations in the network. One day software training included.

System requirements:

Network server and workstation with Intel 2.3 GHz (or higher), 8 GB RAM, 64 bit system, 200 GB HD space, operating system Windows 7 or higher/server OS, 1000 Mbps Ethernet, graphics card NVIDIA GeForce GTX I GTS or equivalent.

The purchase of the CA® SMART 3D software is subject to certification. For further information on training courses etc. please visit www.scheu-dental.com/en/scheu-academy/seminars.

Function OnyxCeph3™ CA® SMART 3D	clinic	laboratory
Scanning		
Importing scan data	✓	✓
Starting the scan process	✓	✓
Interface to specified 3D scanners	✓	✓
Set-up planning/data processing		
Superposition of 3D data sets	✓	✓
Possibility for corrections to already established set-ups	✓	✓
Creation of base moulds	✓	✓
Model labelling	✓	✓
Processing of 2D data sets	✓	✓
Processing of 3D data sets		✓
Virtual segmentation of teeth		✓
Control of set-ups using the CA® philosophy reflected in the software		✓
Digital bracket positioning		
Possibility for corrections to already positioned brackets	✓	✓
3D evaluations		✓
Bracket positioning on facial-axis point		✓
Exporting		
Generation of 3D reports	✓	✓
Exporting unchanged stl-data	✓	✓
Exporting customised set-ups as stl data		✓
Administration		
Interface to administration software, e.g. Computer konkret	✓	✓
Archiving model data	✓	✓



OnyxCeph3™ CA® SMART 3D – clinic

Meets the demands of orthodontists for efficient support in daily practice routines. Enables features such as importing scan data or starting of scan processes, superposition of 3D data sets, creation of base moulds, labelling, exporting and archiving of model data.

- OnyxCeph3™ CA® SMART 3D - clinic #6512
License for up to five workstations in the network

CE



OnyxCeph3™ CA® SMART 3D – laboratory

Meets the demands of lab owners for precise planning, transpositions and control of digital aligner and multi-band solutions. Enables processing of 2D and 3D data sets, virtual segmentation of teeth, digital bracket positioning and export of printing data.

- OnyxCeph3™ CA® SMART 3D - laboratory #6513
License for up to five workstations in the network

CE

Dental



Sleep Medicine



More quality of life.

We assure restful sleep. Trust in the market leader in dental sleep medicine. The IST[®] and TAP[®] appliances set the standard for successful treatment of snoring and sleep apnoea. We recommend our SCHEU-ACADEMY courses on indication and fabrication of these devices.



IST® Kit

Intra-oral anti-snoring device, pressure-moulded. Connection between upper and lower jaw is achieved by continuously adjustable HF telescope. Lower jaw and tongue are advanced by protrusion.

Easy adjustment of protrusion, adjustment can be securely fixed by counter nut.

○ **IST® Kit HF 9**

#5461

incl. 2 telescopes, 4 DURAN® (2.0 mm), 4 ISO-FOLAN®, DURASPLINT® resin, 4 fastening elements with O-ring, 4 SCHEU-anchor, Alastik chain, silicone tube, security key, SD box for patient, occlusion trainer, 2 CETRON® cleansing powder bags of 15 g, instructions and patient information.

1 piece

CE 0044



TELESCOPES HF 5/HF 9/HF 15


Telescopes with limit stop are available in 3 sizes, (HF 5, HF 9 and HF 15). Subsequent activation is possible.

- Telescope HF 5 #5448
- Telescope HF 9 #5450
- Telescope HF 15 #5452

1 pair

CE 0044

NEW: with counter nut



	length (min.)	length (max.)	stroke	activation length
HF 5	21 mm	29 mm	5 mm	3 mm
HF 9	26 mm	40 mm	9 mm	5 mm
HF 15	32 mm	52 mm	15 mm	5 mm




Positioning aid and other accessories

Provides exact angle position for inserting fastening elements.

- Positioning aid for IST® (12°/O°), 1 piece #5469
- Positioning aid for IST® (12°/O°), 4 pieces #5468
(with O-ring and stable retention, incl. silicone tube)
- O-rings for IST®, 10 pieces #5466
- Fastening element for IST® and Herbst*, 4 pieces #5457
(with hexagon socket screw incl. silicone protection tube)
- Hexagon socket screw, 4 pieces #5458

*Herbst is a registered trademark of DENTAURUM GmbH, Ispringen.

CE 0044



GEORGE GAUGE™ Bite Registration

For simple and precise registration of occlusion in snoring therapy and orthodontics. The bite registration is pre-adjusted using 2 or 5 mm bite forks.

- GEORGE GAUGE™ Bite Registration set #5471
comes with 3 bite forks (1 for each size)
- GEORGE GAUGE™ bite fork, 10 pieces #5472
2mm, grey, small (S)
- GEORGE GAUGE™ bite fork, 10 pieces #5473
5mm, white, small (S)
- GEORGE GAUGE™ bite fork, 10 pieces #5462
5mm, white, large (L)

1 piece

CE



CLEAR-BLOKKER®

Light-curing, high-transparent gel for building-up pressure moulded DURAN® splints, reinforced with DURASPLINT® LC. Gel is best suited for removal of air inclusions and insertion of canine guidance. Can also be used for fixation of IST® components on DURAN® splints. Polymerisation in wavelength range 360–420 nm, for example with LC Mini Light or LC-6 Light Oven (see p. 111). Prior to processing, LC primer has to be applied to DURAN® splints.

- CLEAR-BLOKKER® #5261
4 x 1.5 g, 2 cannula (Ø 0.4 mm) and 8 cannula (Ø 0.9 mm)

Single references of cannula can be found on p. 94.

CE

SD Cyano Veneer Fast / SD pipettes

Fast-curing bonding agent, fluid and crystal clear, suitable for fixation of components when fabricating TAP® or IST® splints. Can also be used for secure adhesion of screws in sleep medicine and orthodontics.

SD pipettes ensure precise dosing of SD Cyano Veneer Fast by micro tip.

- SD Cyano Veneer Fast, 5 g #3701

CE

- SD pipettes, 25 pieces #3703



SD Metal Primer

Bonding agent for perfect adhesion of acrylics and any kind of metal. Reduces crevice corrosion. For insertion of TAP® components, IST® fastening elements, polymerised retention and clasp elements in orthodontics and prosthetic acrylics as well as in veneer technique.

- SD Metal Primer #9038

5 ml

CE 0044



SD Box

Retainer box for storing various orthodontic and anti-snoring appliances, such as positioner, IST® splint, TAP®/TAP®-T splint, with mirror and safe closure. Air slits provide optimum ventilation.

Measurements (W x H x L): 95 x 70 x 31 mm

- SD Box *blue* #5439
- SD Box *red* #5437
- SD Box *green* #5438

10 pieces



OPM

For simple fabrication of a pressure moulded OPM protrusion splint in snoring therapy.

Activation is done in 1mm steps by replacing the plastic bars on both sides.

- OPM Kit #5495
incl. 6 pairs of plastic bars (25 - 30 mm long), 4 fastening elements with O-rings, 1 bag CETRON® cleansing powder of 15 g, occlusion trainer, silicone protection tube
- OPM plastic bars, 6 pairs #5494
(25 - 30 mm)

CE



Occlusion trainer

Condylar reposition device for patients' training after wearing protrusion splints during the night.

- Occlusion trainer #5470

10 pieces

CE

TAP® Standard

Thornton Adjustable Positioner (TAP®), enables a controlled and effective treatment of snoring and sleep apnoea. The TAP® splint has been developed by the American dentist and engineer Dr. Keith Thornton.

- TAP® Standard kit #9000
kit incl. TAP® threaded rod, hexagon wrench, occlusion trainer, 5 bags CETRON® cleansing powder of 15 g, SD box blue and patient information
- cover plate incl. 2 screws and hexagon key #9015
- positioning aid incl. socket wrench #9016

CE



TAP®-T

Developed on the basis of TAP® Standard in delicate design. Nickel-free components for optimum bio-compatibility. Upper protrusion adjustment device can be continuously activated up to 7 mm. Reactivation is possible.

Also available as pre-assembled set.

- TAP®-T kit #9030
kit comes with TAP® case, upper protrusion adjustment device made of titanium, lower guiding device made of nickel-free CoCr-alloy, mounting aid, activation key, occlusion trainer, 5 bags of 15 g CETRON® cleansing powder, CETRON® cleansing box and patient information
- TAP®-T kit, pre-assembled #9030.2
kit comes with TAP® case, upper protrusion adjustment device made of titanium, pre-assembled with sticky wax, lower guiding device made of nickel-free CoCr-alloy, mounting aid, activation key, occlusion trainer, 5 bags of 15 g CETRON® cleansing powder, CETRON® cleansing box and patient information
- Mounting aid, 5 pieces #9039

CE

TAP® products are exclusively delivered to certified participants of the TAP® technical course.

TAP® -T Reverse

TAP®-T Reverse incorporates a reduced upper protrusion adjustment device made of titanium. It has been specifically designed for patients with deep bite. Lower protrusion is adjusted by activation key (palatal), as this reduced adjustment device can be inserted behind the upper incisors in the DURASOFT® splint. Reactivation is possible.

Also available as pre-assembled set.

- TAP®-T Reverse kit #9050
kit comes with TAP® case, small upper protrusion adjustment device made of titanium, lower guiding device made of nickel-free CoCr-alloy, mounting aid, activation key, occlusion trainer, 5 bags of 15 g CETRON® cleansing powder, CETRON® cleansing box and patient information
- TAP®-T Reverse set , pre-assembled #9050.2
set comes with TAP® case, upper protrusion adjustment device made of titanium, pre-assembled with sticky wax, lower guiding device made of nickel-free CoCr-alloy, mounting aid, activation key, occlusion trainer, 5 bags of 15 g CETRON® cleansing powder, patient information
- Mounting aid, 5 pieces #9039



TAP®-T / TAP®-T Reverse starter set

Simplified mounting thanks to combination of TAP®-T or TAP®-T Reverse upper protrusion adjustment device and TAP®-T lower guiding device, pre-assembled with sticky wax.

- TAP®-T starter set #9031
set comes with TAP®-T positioning aid, pre-assembled with sticky wax, activation key, SD box blue, patient information
- TAP®-T Reverse starter set #9051
set comes with TAP®-T positioning aid, pre-assembled with sticky wax, activation key, SD box blue, patient information

TAP® products are exclusively delivered to certified participants of the TAP® technical course.

2 mm additional mobility!

TAP® T XL Lower guiding device

Increased version (2 mm additional mobility) for more wearing comfort. Can be used with TAP®-T and TAP®-T Reverse upper protrusion adjustment devices.

○ TAP® T XL Lower guiding device, incl. matrix #9029



TAP® -T Matrix*

Reusable auxiliary element made of heat-resistant silicone for efficient fabrication of the sticky wax spacer. Makes sure the TAP®-T titan hook runs free in the guiding device in the lower splint.

○ TAP® -T Matrix ,2 pieces #9037

* Not for use with TAP® -T XL lower guiding device.



TAP® Splint case for patients

Zippered case for storing and carrying TAP® and TAP®-T splints.

○ TAP® Splint case for patients #9082

TAP® products are exclusively delivered to certified participants of the TAP® technical course.

NEW!**TAP® PAP**

Comfortable nasal pillow mask designed for patients with low CPAP compliance. Since the TAP® PAP mask is anchored to the upper teeth with a special ThermAcryl™ mouthpiece, it is very stable, will not move or leak nor leave annoying marks on the nose and face.

- TAP® PAP comes with nasal pillow in small, medium and large

#9073

CE

myTAP™

Diagnostic testing appliance for getting acquainted with the simple and patient-friendly treatment of snoring and sleep apnoea. Metal-free and easy adjustable trays enable protrusion up to 20 mm, assuring considerable lateral mobility.

- myTAP™ upper and lower jaw splint, incl. nut, storage box and instructions

#9070

CE

NEW!**AM Aligner**


Thermoplastic occlusion trainer for condylar reposition after wearing protrusion splints during the night. The AM Aligner is heated and adapted to the patient's upper jaw by the clinician.

- AM Aligner

#9072

CE





50 years of competence
and experience.

Tried and tested. Our pressure moulding material has been tested in terms of cell biology according to the international standards DIN EN ISO 10993 and 7405 and has been released for clinical purposes. The bar code on each blank ensures quick and simple handling.

Pressure Moulding Technique





BIOCRYL® C

Break-resistant acrylic material of pure PMMA, monomer-free, without allergenic irritants, bonds to acrylic.

Application

Orthodontic plates and retainers

○ 2.0 × 125 mm*	clear	#3147
● 2.0 × 125 mm**	red	#3149
● 2.0 × 125 mm**	blue	#3150
● 2.0 × 125 mm**	yellow	#3151
○ 2.0 × 125 mm*	rose-transparent	#3152
○ 3.0 × 125 mm*	clear	#3153

*10 /100 pieces ** 10 pieces

CE

BIOCRYL® M

Hard-elastic, monomer-free material, without allergenic irritants, bonds to acrylic.

Application

Orthodontic plates and retainers

○ 2.0 × 125 mm	Multicolor Set*	#3228
● 2.0 × 125 mm	Polka-Dots	#3229
● 2.0 × 125 mm	Zebra	#3230
● 2.0 × 125 mm	Rainbow	#3231
● 2.0 × 125 mm	Tiger	#3232
● 2.0 × 125 mm	Camouflage	#3233

10 pieces *5 x 2 pieces

CE



NEW!

BIOCRYL® Form X

Visco-elastic material, does **not** bond to acrylic, blue.

Application

Radiopaque material for the exact visualisation of the clinical situation in computer tomography.

● 1.5 × 125 mm	blue	#3148
----------------	------	-------

CE

DURAN®

Hard-elastic transparent material, bonds to DURASPLINT® and DURASPLINT® LC. Apply LC primer prior to using DURASPLINT® LC. No pre-drying required.

Application

Miniplast splints, occlusal splints, interim and fluorisation splints, drilling/X-ray templates, ant-snoring appliances, invisible retainers, temporary crowns and bridges.

- | | |
|------------------|-------|
| ○ 0.4 × 125 mm* | #3470 |
| ○ 0.5 × 125 mm | #3411 |
| ○ 0.625 × 125 mm | #3409 |
| ○ 0.75 × 125 mm | #3413 |
| ○ 1.0 × 125 mm | #3415 |
| ○ 1.5 × 125 mm | #3417 |
| ○ 2.0 × 125 mm | #3419 |
| ○ 3.0 × 125 mm | #3421 |

10 /100 pieces * 10 pieces

CE

Now available in 0.4 mm material thickness!



DURAN®+

Hard-elastic transparent material, bonds to DURASPLINT® and DURASPLINT® LC. Apply LC primer prior to using DURASPLINT® LC. **With integrated insulation foil.** No pre-drying required.

Application

Miniplast splints, occlusal splints, interim and fluorisation splints, X-ray/drilling templates, anti-snoring devices, temporary crowns and bridges, invisible retainer.

- | | |
|------------------|-------|
| ○ 0.4 × 125 mm* | #3472 |
| ○ 0.5 × 125 mm | #3426 |
| ○ 0.625 × 125 mm | #3428 |
| ○ 0.75 × 125 mm | #3430 |
| ○ 1.0 × 125 mm | #3432 |
| ○ 1.5 × 125 mm | #3434 |
| ○ 2.0 × 125 mm | #3436 |

10 /100 pieces *10 pieces

CE

Now available in 0.4 mm material thickness!





DURAN®+ A2 pd

New material in proven DURAN® quality, now available in aesthetic colour (VITA A2). **With integrated insulation foil.** Already pre-dried and vacuum sealed for instant use.

Application

Tooth coloured temporaries, splints for planning and demonstration for aligner treatments

- 0.5 × 125 mm #3422
- 0.625 × 125 mm #3424
- 0.75 × 125 mm #3423

10 pieces

CE

NEW!

DURAN®+ white pd

New material in proven DURAN® quality, now available in white (Vita OM1). **With integrated insulation foil.** Already pre-dried and vacuum sealed for instant use.

Application

Cosmetic splints

- 0.5 × 125 mm #3473
- 0.625 × 125 mm #3474
- 0.75 × 125 mm #3475

10 pieces

CE



DURASOFT® pd

Unbreakable, abrasion-resistant, transparent compound material – hard on one side and soft on the other. Hard side bonds to acrylic. Already pre-dried and vacuum sealed for instant use.

Application

Miniplast splints, bruxism splints, occlusal splints

- 1.2 × 125 mm* *hard/soft 0.7 mm / 0.5 mm* #3342
- 1.8 × 125 mm* *hard/soft 0.9 mm / 0.9 mm* #3388
- 2.5 × 125 mm* *hard/soft 1.6 mm / 0.9 mm* #3387
- 3.0 × 125 mm** *hard/soft 2.1 mm / 0.9 mm* #3341
- 4.0 × 125 mm** *hard/soft 3.1 mm / 0.9 mm* #3345

*10 /100 pieces ** 10 pieces

CE



DURASOFT® blue pd

The proven break and abrasion resistant material in blue-transparent colour. Already pre-dried and vacuum sealed for instant use.

Application

Specifically for simple fabrication of coloured intra-oral anti-snoring devices. Blue colour increases aesthetics.

- 1.8 × 125 mm #3347
- 2.5 × 125 mm #3348

10 pieces

CE

DURASOFT® seal

Clear pressure moulding material made of thermoplastic polyurethane for use with DURASOFT® pd. Forms a robust inner layer and provides high quality hard-soft splints with excellent resistance against discolouration. Also available in blue and magenta.


Application

Discolouration resistant inner layer for DURASOFT® pd splints. Specifically designed for long-term use splints such as intra-oral anti-snoring devices.

- 0.1 × 125 mm clear #3338
- 0.1 × 125 mm blue #3339
- 0.1 × 125 mm magenta #3337

10 pieces

CE



BIOPLAST® XTREME Kit*

Two different degrees of hardness in one pressure moulding blank thanks to the integrated reinforcement in the delicate anterior region. Absorbs hard and flat impacts in contact sports. V-shaped harder insert for variable size of front tooth protection.

Kits come with 1 BIOPLAST® XTREME 5.0 x 125 mm, 1 bag of 15 g CETRON® cleansing powder, 1 BIOPLAST® insulating agent, retainer box.

- BIOPLAST® XTREME Kit, red #3283
- BIOPLAST® XTREME Kit, blue #3284
- BIOPLAST® XTREME Kit, black #3285

* only for MINISTAR® from 4/92 and BIOSTAR® from 4/99



BIOPLAST® XTREME PRO Kit*

For professional sports where hard and flat impacts often occur. Additional hard-elastic layer for integrated reinforcement in the anterior region.

Kits come with 1 BIOPLAST® XTREME 4.0 x 125 mm, 1 BIOPLAST® clear 2.0 x 125 mm, 1 DURAN® 0.75 x 125 mm, 1 bag of 15 g CETRON® cleansing powder, 1 BIOPLAST® insulating agent, retainer box.

- BIOPLAST® XTREME Kit, red #3292
- BIOPLAST® XTREME Kit, blue #3293
- BIOPLAST® XTREME Kit, black #3294

* only for MINISTAR® from 4/92 and BIOSTAR® from 4/99



BIOPLAST® XTREME*

BIOPLAST® XTREME material is ideally suited for simple and safe preparation of individual mouth guard. Hard V-shaped reinforcement layer ensures additional protection in the anterior region, allowing precise positioning in the desired extension. Soft layer can be used for lamination with BIOPLAST® or for insertion of the counter bite.

- 4.0 x 125 mm red #3295
- 4.0 x 125 mm blue #3296
- 4.0 x 125 mm black #3297
- 4.0 x 125 mm clear #3298
- 5.0 x 125 mm red #3286
- 5.0 x 125 mm blue #3287
- 5.0 x 125 mm black #3288
- 5.0 x 125 mm clear #3289

10 pieces * only for MINISTAR® from 4/92 and BIOSTAR® from 4/99





BIOPLAST® XTREME Deco*

The proven XTREME material for the fabrication of individual mouth guard is now available in many individual designs. A new manufacturing technology enables photo-realistic, multi-colour designing. Material does **not** bond to acrylic.

○ 5.0 × 125 mm	<i>chinese</i>	#3312
○ 5.0 × 125 mm	<i>dragon</i>	#3313
○ 5.0 × 125 mm	<i>kiss</i>	#3314
○ 5.0 × 125 mm	<i>aqua</i>	#3315
○ 5.0 × 125 mm	<i>stone</i>	#3316
○ 5.0 × 125 mm	<i>lion</i>	#3317

5 pieces

CE

* only for MINISTAR® from 4/92 and BIOSTAR® from 4/99

BIOPLAST® Color

Soft, pliable rubber like material, does **not** bond to acrylic.

Application

Coloured mouth guards, positioners

● 3.0 × 125 mm	<i>blue-transparent</i>	#3356
● 3.0 × 125 mm	<i>red-transparent</i>	#3355
○ 3.0 × 125 mm	<i>white</i>	#3445
● 3.0 × 125 mm	<i>black</i>	#3448
● 3.0 × 125 mm	<i>neon red</i>	#3352
● 3.0 × 125 mm	<i>neon green</i>	#3353
● 3.0 × 125 mm	<i>azure</i>	#3354
○ 3.0 × 125 mm	<i>Color set</i>	#3271

10 pieces

CE



BIOPLAST® Multicolor

Soft, pliable rubber like material, does **not** bond to acrylic.

Application

Coloured mouth guards, positioners

● 3.0 × 125 mm	<i>Germany</i>	#3390
● 3.0 × 125 mm	<i>France</i>	#3391
● 3.0 × 125 mm	<i>Italy</i>	#3392
● 3.0 × 125 mm	<i>Spain</i>	#3393
● 3.0 × 125 mm	<i>Sweden</i>	#3394
● 3.0 × 125 mm	<i>Brazil</i>	#3395
● 3.0 × 125 mm	<i>Transylvania</i>	#3396
● 3.0 × 125 mm	<i>Atlantis</i>	#3397
● 3.0 × 125 mm	<i>White dots</i>	#3398
● 3.0 × 125 mm	<i>Green dots</i>	#3399
○ 3.0 × 125 mm	<i>Multicolor-Set</i>	#3400

10 pieces

CE





BIOPLAST®

Soft pliable rubber like material, does **not** bond to acrylic, clear-transparent.

Application

Duplicating moulds, mouth guard, bruxism splints, positioner, bleaching splints, OSAMU-Retainer®

- 1.0 × 125 mm #3270
- 1.5 × 125 mm #3183
- 2.0 × 125 mm #3185
- 3.0 × 125 mm #3187
- 4.0 × 125 mm** #3188
- 5.0 × 125 mm*/** #3374

10/100 pieces * 10 pieces
 **only for MINISTAR® from 4/92 and BIOSTAR® from 4/99

CE



BIOPLAST® bleach

Soft pliable rubber like material, best suited for fabrication of bleaching splints. Material does **not** bond to acrylic .

- 1.0 × 125 mm

#3190

10/100 pieces

CE



NEW: BRUX CHECKER® white



BRUX CHECKER® acc. to Prof. Sato

Foil with one-side colour coating for analysis and therapy planning as well as for control of prosthetic restorations.

- 0.1 × 125 mm, red
- 0.1 × 125 mm, white

#3209

#3210

10 pieces

CE

Manual "Diagnosis of occlusion in bruxism" (Kanagawa Dental College, Japan) can be found at www.scheu-dental.com/en/downloads/other

IMPRELON®

Hard material, not elastic, bonds to acrylic.

Application

Material in clear and opaque for impression and medical trays, material in white for chin caps and bite registration.

○ 2.0 × 125 mm	clear	#3164
○ 3.0 × 125 mm*	clear	#3168
○ 3.0 × 125 mm	white	#3274
○ 3.0 × 125 mm	opaque	#3361

10 pieces * 10/100 pieces

CE



IMPRELON® S pd

Extremely break and abrasion resistant material. New formulation! Already pre-dried and vacuum sealed for instant use.

Application

OSAMU-Retainer®, miniplast splints, retention splints

○ 0.5 × 125 mm	clear	#3320
○ 0.75 × 125 mm*	clear	#3322
○ 1.0 × 125 mm*	clear	#3324
○ 1.5 × 125 mm	clear	#3326
○ 2.0 × 125 mm	clear	#3328
○ 3.0 × 125 mm	clear	#3330

10 pieces * 10/100 pieces

CE

OSAMU-Retainer® Kit

Transparent retainer acc. to Dr. Yoshii. Pressure-moulded using BIOPLAST® 1.5 mm and IMPRELON® S pd 0.75 mm and OSAMU Bond for optimal compound of material. OSAMU-Retainer® Kit for fabrication of 10 retainers. All components also available separately.

○ OSAMU-Retainer® Kit	#3242
10 IMPRELON® S pd (0.75 x 125 mm), 10 BIO-PLAST® (1.5 x 125 mm), 10 ml OSAMU-bond	

1 piece

CE



COPYPLAST®

Visco-elastic material, opaque, does **not** bond to acrylic.

Application:

Duplicating moulds, bleaching splints, bracket transfer trays, copings, moulds for temporary crowns and bridges

- 0.5 × 125 mm #3177
- 0.75 × 125 mm #3169
- 1.0 × 125 mm #3179
- 1.5 × 125 mm #3175
- 2.0 × 125 mm #3181

10/100 pieces

CE



COPYPLAST® C

Transparent, visco-elastic material, unbreakable and abrasion resistant, does **not** bond to acrylic.

Application:

Long term and invisible retainers, orthodontic splints

- 1.0 × 125 mm #3403

10/100 pieces

CE

HARDCAST®

Hard-elastic material, does **not** bond to acrylic, combustible leaving no residue, opaque.

Application:

Copings, spacer foil for bleaching splints

- 0.4 × 125 mm #3194
- 0.6 × 125 mm #3195
- 0.8 × 125 mm #3200

10 pieces

CE





Spacer Foil

For fabrication of copings, will be removed before embedding the model, clear.

○ 0.1 × 125 mm

#3202

100 pieces



ISOFOLAN®

Ensures optimum insulation for plaster against acrylic in the pressure moulding technique. Best suited as spacer foil when pressurising hard-elastic material for splints.

○ 0.1 × 125 mm

#3207

100 pieces



Skin Packaging

Protection packaging for sending out high quality models, consisting of cardboard and coating foil.

○ Skin Packaging

#3206

○ Coating foil (0.15 x 125 mm)

#3203

Protection cover for study and orthodontic models

100 pieces

Indirect bonding technology

Indirect bonding for precise transfer of the bracket position with considerable reduction of treatment time.

Thanks to the specifically designed SD products described on these pages, the transfer matrix can be quickly and easily fabricated in the laboratory.



SD Bond

Water soluble, thixotropic adhesive for temporary fixing and positioning of brackets on the model. Optimum curing times and complete water solubility.

○ SD Bond, 5 cartridges of 3 g
incl. cannula (ø 1,2 mm)

#3602

SD Silicone

Silicone with excellent pouring and adhesive characteristics for fabricating bracket transfer trays. For blocking-out and precise, stable fixing of brackets on the models during pressure-moulding, can be removed leaving no residue. Comes in two-component syringe, with 15 mixing cannula.

○ SD silicone, double cartridge of 2× 10ml
with 15 mixing cannula

#3603



SD Scalpel

Scalpel for precise orientation cuts in the bracket transfer foil in the area of the bracket slots, serving as basis for the future extension of the tray and for finishing.

○ SD scalpel with 6 blades

#3604





SD Cutter

Smooth knife for precise reduction of the pressure-moulded bracket transfer tray below the bracket slots. Blade can be fixed for precise cuts, doesn't require further finishing.

○ SD Cutter

#3605

SD Separating disc, serrated

Discs for precise and smooth cuts in the bracket transfer foil in the area of bracket slots, serving as orientation for the future extension of the tray and for finishing. Can also be used for finishing DURAN®.

○ SD Separating disc, serrated
3 discs, with mandrel

#3606



SD Bracket pincers

Pliers with delicate tip for inserting and removing brackets from pressure-moulded bracket transfer tray. Made of stainless steel, sterilisable.

○ SD Bracket pincers

#3607



SD Bracket remover

To release brackets from bracket transfer tray after pressure-moulding and for easier removal of the bracket transfer tray from patient's mouth. Made of stainless steel, sterilisable.

○ SD Bracket remover

#3608





- BLUE-BLOKKER® #5276
4 x 1,5 g, 2 cannula (Ø 0.6 mm) and
8 cannula (Ø 0.9 mm)
- Micro cannula* (Ø 0,4 mm) #5277
- Micro cannula* (Ø 0,6 mm) #5264
- Micro cannula* (Ø 0,9 mm) #5278

*5 pieces

BLUE-BLOKKER®

Light-curing, blue-transparent blocking-out and modelling putty with excellent flow characteristics. For fabrication of space maintainers for bleaching splints and medical trays as well as for sealing and blocking set-up models. Combustible, leaving no residue, can therefore also be used in prosthetics when fabricating fixed or removable prostheses. Telescopic models are stable and can be removed without tension. Curing in the wavelength range 380-470 nm.



SIL-KITT

Permanently plastic silicone, used specifically for blocking-out undercuts in pressure moulded splints, anti-snoring and orthodontic appliances.

- SIL-KITT red #3442
- SIL-KITT transparent #3443

150 g



Blocking-out putty

Permanently plastic putty for blocking-out large-scale undercuts.

○ *Blocking-out putty*

#3220

250 g



BIOPLAST® Insulation agent

Insulation for plaster models when using BIOPLAST® material for fabricating positioners and mouth guards. BIOPLAST® insulating agent provides absolute transparency and smooth inner surfaces. Also available in practical applicator bottle of 100 ml for quick and non-drip insulation of plaster models.

○ *BIOPLAST® Insulating agent, 100 ml*

#3166

Applicator bottle with 3 tips

○ *BIOPLAST® insulating, 500ml*

#3189



DURASPLINT®

High-transparent cold-cure resin specifically designed for fabricating occlusal splints with adjusted occlusion. Very simple to use, little shrinkage, optimised trimming.

- DURASPLINT® Kit #3257
incl. 200 g polymer and 100 ml monomer, mixing and dosing cup, 20ml orthodontic insulating agent incl. brush, pipette, glass spatula and detailed instructions
- DURASPLINT® polymer, 200/400 g #3258
- DURASPLINT® monomer, 100/250 ml #3259

CE



BIOCRYL®-RESIN

Special orthodontic acrylic for attaching clasp arms and relining expansion screws when pressure-moulding retention and expansion plates using BIOCRYL® C and BIOCRYL® M. Material polymerises during the cooling phase in the pressure moulding unit. No time-consuming polymerisation in pressure pot.

- BIOCRYL®-RESIN polymer, 400/1000 g #3221
- BIOCRYL®-RESIN monomer, 250/500 ml #3218

CE



DURASPLINT® LC

Light-curing, transparent resin, specifically designed for the adjustment of DURAN® splints with a minimal thickness of 1.5 mm. Doesn't require curing in autoclave. For polymerisation, we recommend using the LC 6 light oven (see p.111) with pre-set polymerisation times of 5 or 10 min. Curing in the wavelength of 350 – 450 nm. Prior to processing, LC primer has to be applied to DURAN® splints.

No yellow discolouration of material after light-curing.

- DURASPLINT® LC Kit #3256
3 bars (38 g in total) DURASPLINT® LC, 1.5 g CLEAR BLOKKER® with 2 cannula, 10 ml LC primer with brush, 10g LC insulation, 1 LC bur A, B and C each, 1 SD polishing paste A and B each, 5 SD polishing brushes, SD felt discs and SD high-gloss polishing mops with mandrel
- DURASPLINT® LC, 3 bars, 38 g #3261
W x L x H 32 x 120 x 8 mm

CE



CLEAR-BLOKKER®

Light-curing, high-transparent gel for building-up pressure moulded DURAN® splints, reinforced with DURASPLINT® LC. Gel is best suited for removal of air inclusions and insertion of canine guidance. Can also be used for fixation of IST® components on DURAN® splints. Polymerisation in wavelength range 360–420 nm, for example with LC Mini Light or LC-6 Light Oven (see p. 111). Prior to processing, LC primer has to be applied to DURAN® splints.

- CLEAR-BLOKKER® #5261
4 x 1.5 g, 2 cannula (Ø 0.4 mm) and 8 cannula (Ø 0.9 mm)

Single references of cannula can be found on p. 94.

CE



LC Primer

Light-curing primer for pre-treating DURAN® splints before adjusting them with DURASPLINT® LC or CLEAR BLOKKER®. For polymerisation, we recommend using the LC 6 Light Oven (see p.111) with pre-set polymerisation times of 5 or 10 min. Curing in the wavelength of 350 – 450 nm.

○ LC Primer, with brush

#3263

10 ml

CE



LC Insulation

Facilitates modelling DURASPLINT® LC on DURAN® splints. Ensures clean and neat working results with smooth surfaces, considerably simplifying finishing.

○ LC Insulation

#3266

10 g



SD Composite primer

Bonding agent for fabrication of temporary crowns and bridges in the dental pressure moulding technique. Primer ensures a secure adhesion between composite resin (also cartridge dispensing systems) and DURAN®. Brush for easy application of primer is integrated in bottle cover.

○ SD Composite primer

#3705

25 ml

CE 0044



OSAMU-Bond

Bonding agent for optimal compound of BIOPLAST® and IMPRELON® S pd in the pressure moulding technique.

○ OSAMU Bond

#3246

10 ml

CE

NEW: LC Bur C



LC Burs

Specifically designed for finishing cold cure resin or light-cured adjustments on splints. Allows neat cutting without clogging and leaving any residue. LC bur A comes in conical form for specific finishing or occlusal grinding. LC bur B is pear-shaped, ensuring very smooth grinding surface patterns when it comes to larger surfaces and LC bur C comes in bud form for quick and uniform removal of roughness.

- LC Bur A #3267
- LC Bur B #3268
- LC Bur C #3265

SD Polishing paste

Specifically designed for finishing (SD polishing paste A) and for high-gloss polishing (SD polishing paste B) of splints, adjusted with DURASPLINT® or DURASPLINT® LC. Can also be used for polishing orthodontic appliances.

- SD polishing paste A, bar of 60g #3710
- SD polishing paste B, bar of 100g #3711



SD Polishing brush

High-quality goat hair brush with micro fibre insert, for use in combination with SD polishing paste A and B for polishing adjusted splints and orthodontic appliances.

SD Felt wheel

For use in combination with SD polishing pastes when polishing built-up splints and orthodontic appliances.

SD High-gloss polishing mop

Extra soft mop structure for finishing and high-gloss polishing of built-up splints and orthodontic appliances.



- SD polishing brush with mandrel #3712
- SD felt wheel with mandrel #3713
- SD high-gloss polishing mop with mandrel #3714

5 pieces

Impression tray handles

Stable, handy plastic handles which can be fixed to tray with cold cure resin and individually formed by heating.

○ *Impression tray handles* #3216
 10/100 pieces €



Foam discs

Foam discs can be used as space holder and insulation for pressure moulded impression trays, providing retention to impression material.

○ *Foam discs (Ø 103 mm)* #3444
 10 pieces



SD Foil scissors A

SD foil scissors A are delicately made of a high-quality auto-clavable stainless steel alloy. For better handling, scissors are laterally curved. Micro serration prevents slipping.

○ *SD foil scissors A* #3460
 1 piece €

Pellet Cover

Used for covering stainless steel pellets while pressurising splints. Cover is reusable and avoids adhesion of pellets to the heated soft material.

- Pellet cover, large #3006
- Pellet cover, small #3007

5 pieces



Stainless steel pellets

Stainless steel pellets, highly polished by galvanic treatment, for embedding models in pressure moulding machines. Thanks to the cylindrical pellet shape, models are securely positioned and can't be pressed into the pellets.

- Stainless steel pellets #3222

1000 g

Pellet brush

With integrated magnet for quick removal of pellets and lifting the model platform from the model cup.

- Pellet brush #3136

SD Foil scissors B

SD foil scissors B are robust and durable for cutting splints or inserting cuts into splints after pressurising.

- SD foil scissors B #3461

1 piece



Finishing Set

Set containing eight grinding instruments for cutting, trimming and finishing of soft-elastic and hard-elastic material. With color coding for easy identification.

○ **Finishing Set**

#3378

Set comes with one each cutting bur, finishing bur, HM carbide bur, HM tricutter fine, BIOPERM trimmer, POLYFIX, OSAMU polisher, DIMO®PRO SLIM and dressing stone.



Cutting bur

Spiral cutting bur for quick, neat cutting of hard pressure moulding material.

HM Carbide cutter

Triangular carbide cutting bur for quick, neat cutting of hard pressure moulding material.

HM Tricutter fine

Fine triangular bur for interdental finishing of splints.

Finishing bur

Square cutting bur made of special soft alloy that won't damage clasps and wires.

CA® Carbide bur

CA® carbide bur made of fine grit carbide is characterised by high cutting performance, vibration-free concentric running and very long endurance.

- Cutting bur (without) #3214
- HM Carbide cutter (yellow) #3369
- HM Tricutter fine (green) #3370
- Finishing bur (red) #3377
- CA® Carbide bur* #5262

6 pieces * 1 piece

BIOPERM Trimmer

For trimming and initial polishing of soft-elastic and compound material.

POLYFIX

For trimming and initial polishing of soft-elastic and compound material.

OSAMU Polisher

Finishing and polishing of soft-elastic and compound material.

- *BIOPERM trimmer, blue, 2 pieces* #3226
- *POLYFIX on mandrel, brown, 1 piece* #3371
- *POLYFIX without mandrel, brown, 6 pieces* #3372
- *OSAMU polisher, grey, 2 pieces* #3247



DIMO®

Stable and flexible trimming wheels with prolonged processing time, three types of grains, incl. mandrel.

Application:

Trimming and initial polishing of elastomer, sandwich material, mouth guard and orthodontic appliances.

- *DIMO® brown, coarse** #3382
- *DIMO® green, medium** #3383
- *DIMO® maroon, fine** #3384
- *DIMO® mandrel* #3380

10 pieces *10/50 pieces



DIMO®PRO / DIMO®PRO SLIM

Trimming wheel made of special elastomer with fine fibrous structure, extremely long working time, incl. mandrel.

Application

DIMO®PRO: For finishing large-scale marginal areas and for final and high-gloss polishing of splints, positioners, mouth guard and orthodontic appliances. Provides highly-polished rims and uniform surfaces of hard/soft compound material.

DIMO®PRO SLIM: for finishing splints, specifically in the area of ligaments and interdentality.

- *DIMO®PRO disc 24 × 6 mm* #3381
- *DIMO®PRO SLIM* disc 24 × 4 mm* #3376

10/50 pieces * tapered rim

Laboratory equipment

We are setting benchmarks!

We are setting benchmarks, instead of just meeting standards.

From O-ring to the microprocessor control – our specialists in the equipment manufacturing are well versed. All parts are manufactured with maximum precision and attention to detail, which is why our pressure moulding units are delivered with a three-year warranty.







NEW!

Asiga MAX

The IMPRIMO® system now includes the Asiga MAX 3D printer, offering intuitive operation with fast and precise printing results. The built-in LCD panel with touch screen function provides real-time information on the building process and enables interactive control. Multiple print jobs can be stored on the internal solid-state drive and managed via touch screen or remotely via web browser. Printing jobs can be transferred over Ethernet or Wifi. In the absence of networking infrastructure, the Asiga MAX can operate as a wireless hot-spot for an independent production solution. Thanks to the IMPRIMO® 385 nm material range and reliable DLP technology, the Asiga MAX can be used in all fields of dentistry. Printing jobs can be created using the Asiga Composer software, included license-free with every system.

- *Asiga MAX* #6501
- *Tray* #6505



IMPRIMO® LC Model

Light-curing resin at the base of methacrylate for the fabrication of high-precision dental models with smooth surfaces. Compatible with 3D printing system with 385 nm light source.

○ IMPRIMO® LC Model

#6502



IMPRIMO® LC Impression

Light-curing resin at the base of methacrylate for the fabrication of individual impression trays and bases for bite registrations, neon orange. Compatible with 3D printing system with 385 nm light source.

○ IMPRIMO® LC Impression

#6506

CE



IMPRIMO® LC Splint

Light-curing resin at the base of methacrylate for the fabrication of high-precision transparent occlusal splints and drilling templates. When cured, the resin is characterised by relative flexibility and provides high wearing comfort, reduced fracture susceptibility and excellent polishing properties. Contains no diluting monomers, compatible with 3D printing system with 385 nm light source.

○ IMPRIMO® LC Splint

#6503

CE

BIOSTAR®

Versatile pressure moulding machine for all applications in dental pressure moulding and orthodontics, equipped with scan technology and 6.0 bar working pressure, adjustable by pressure reducing valve on the unit. The unique thermally controlled infra-red heater reaches working temperature in one second. All parameters such as temperature, heating and cooling times are programmed by reading the scanner or entering the corresponding code. **The coloured 3.2" LCD display with a resolution of 320 x 240 pixel leads you through the single working steps of the pressure moulding process, providing further information on the scanned material, its indication as a well as a library with illustrated application hints for fabrication. Rear port for easy updates.** Large model cup facilitates embedding of mounted models. Pellet container ensures optimum dispensing and storing of pellets.



Power: 230/115/100 V, 750 W

Working pressure: 0.5 - 6.0 bar

Dimensions (W x H x D): 520 x 240 x 440 mm

Weight: 14 kgs

3 years warranty, 10 years availability of spare parts

- *BIOSTAR® 230 V* #3001
- *BIOSTAR® 115 V* #3011
- *BIOSTAR® 100 V* #3021
- *BIOSTAR® protection cover made of transparent plastic, avoids soiling* #3449



New generation

MINISTAR S®

The compact pressure moulding unit MINISTAR S® equipped with scanning technology: scanner reads the coded material, allowing programming of all important parameters such as heating time, temperature and cooling time. In other words: precise pressure moulding results are reached in record time thanks to the extremely fast, patented heater and working pressure of 4.0 bar. **The new generation is equipped with a user-friendly membrane keypad and a large display. All operating parameters are monitored and indicated on the display.**

Power: 230/115/100 V, 750 W

Working pressure: 0.5 - 4.0 bar

Dimensions (W × H × D): 500 × 210 × 400 mm

Weight: 9 kgs

3 years warranty, 10 years availability of spare parts

- | | |
|--|-------|
| <input type="radio"/> MINISTAR S® 230 V | #3501 |
| <input type="radio"/> MINISTAR S® 115 V | #3511 |
| <input type="radio"/> MINISTAR S® 100 V | #3521 |
| <input type="radio"/> MINISTAR S® protection cover made of transparent plastic, avoids soiling | #3447 |



Sampler of pressure moulding material

Sampler of all available pressure moulding materials with information on thickness, diameters and colours. Helpful information on working procedures and chemical and physical parameters is given as well.

- | | |
|---|-------|
| <input type="radio"/> Sampler of pressure moulding material | #3137 |
|---|-------|

CE



TWINSTAR® P

Pressure moulding machine for basic dental applications in the pressure moulding technique such as temporaries, splints, mouth guard, bite registrations and individual impression trays.

Characteristics

Medium-wave, time-controlled heater, programmed cooling phase, pre-adjusted 2 bar working pressure, range of material 0.1-3.0 mm, heating of the material side which is directly placed over the model, large model cup and ergonomic pellet container, simple and functional "one-switch-operation".

Power: 230/115/100 V, 340 W

Working pressure: 0.1 - 2.0 bar

Dimensions (W x H x D): 480 x 210 x 430 mm

Weight: 7.5 kgs

3 years warranty, 10 years availability of spare parts

- TWINSTAR® P 230 V
- TWINSTAR® P 115 V
- TWINSTAR® P 100 V

- #3140
- #3141
- #3138



TWINSTAR® V

Vacuum machine, doesn't require compressed air connection. For all basic dental applications in the thermo-forming technique such as temporaries, splints, mouth guard, bite registration, and individual impression trays.

Characteristics

Medium-wave, time-controlled infra-red heater, programmed cooling phase, range of material 0.1-3.0 mm, heating of the material side which is placed directly over the model, large model cup and ergonomic pellet container, easy and functional "one-switch-operation".

Power: 230/115/100 V, 340 W

Working vacuum: 0.8 bar

Dimensions (W x H x D): 440 x 210 x 430 mm

Weight: 7.5 kgs

3 years warranty, 10 years availability of spare parts

- TWINSTAR® V 230 V
- TWINSTAR® V 115 V
- TWINSTAR® V 100 V

- #3142
- #3143
- #3139

LC 6 Light Oven

High-capacity light polymerisation unit for up to 4 models and fixator. Its inner dimensions of 170 x 160 x 120 mm (W x D x H) offer enough space for articulated and mounted models. Six fluorescent light tubes (3 UVA and 3 blue light tubes) and an aluminium reflector guarantee homogeneous illumination. Two polymerisation times of 5 and 10 minutes are already pre-set in accordance with the light-curing DURASPLINT® LC material. Continuous operation can also be set. Spectral range of the light curing device between 350 and 450 nm. Auto start function for easier application.

Power: 230/115 V

Dimensions (W x H x D): 240 X 180 x 270 mm :

Weight: 3 kgs

3 years warranty, 10 years availability of spare parts



○ LC 6 Light Oven

#5610

○ LC 6 Light Oven

#5620

LC Mini Light

LED curing light with 14 LEDs designed for fixation and initial polymerisation of DURASPLINT® LC and CLEAR-BLOKKER®. Ideal for fixing fastening parts of IST® or OPM devices in combination with the LC primer. Final light polymerisation in stationary light curing unit such as LC 6 Light Oven. Use with customary AAA batteries.



Wavelength: 400-405 nm

weight (without batteries): 96 g

length (without shade): 120 mm

○ LC Mini Light with shade

#5606



Heating Support

Can be attached directly to the heater when fabricating positioners and mouth guard, providing the optimum distance between model and heater when heating the occlusal surfaces and for compound of upper and lower splint.

Measurements: (W x H x D): 190 x 160 x 70 mm

○ Heating support for BIOSTAR® and MINISTAR S® #3452



Model Platform

Model platform with 85 mm diameter provides increased surface for dental arcade models made of resin and hard plaster preventing the model from getting pressed into the stainless steel pellets.

○ Model platform, Ø 85 mm

#3134

Split-Fixator

Stable occlusal holder with sliding arms. By unscrewing the centring tips, models can be removed without any damage and reinserted at any time in exactly the same position. Precise dovetailed guidance avoids twisting of fixator arms.

○ *Split-fixator* #5380



Hinge fixator

Thanks to the hinge in the upper part of the fixator, grinding of functional orthodontic appliances and fabrication of bimaxillary appliances becomes easier. Mounting splints on fixator arms ensure easy removal and precise repositioning of models.

○ *Hinge fixator* #5374
○ *peaked screws, 10 pieces* #5376

Mounting splints

For hinge fixator

○ *Mounting splints* #5375
10 pieces



Model holder

With flexible clamps and movable joints, keeps small and large models fixed in any required position.

○ Model holder

#5270

1 piece



Parallelometer

Indispensable auxiliary instrument when fabricating precise adjusted splints. The parallelometer enables determination of insertion direction and marking of tooth equator. The extensions of DURAN® splints, built-up with DURASPLINT® LC, should not exceed the tooth equator to avoid cracks or fissures when inserting or removing the splint.

○ Parallelometer

#5362

incl. model holder, graphite leads and measurement accessories

1 piece



Indication Chart on Applications in the Pressure Moulding Technique

Indication	Product	Recommended Thickness
Anti-snoring splint	DURAN®	2.0 mm
	DURAN®+	2.0 mm
	DURASOFT® pd with DURASOFT® seal	1.8 / 2.5 mm
Bite taking / Bite registration	IMPRELON® white, clear	2.0 / 3.0 mm
Bleaching splint	BIOPLAST® bleach	1.0 mm
	COPYPLAST®	1.0 mm
Bracket transfer mask	COPYPLAST®	0.5 / 0.75 mm
Copings	HARDCAST®	0.4 / 0.6 / 0.8 mm
	COPYPLAST®	0.5 / 0.75 mm
Cosmetic splint	DURAN®+ A2	0.5 mm - 0.75 mm
	DURAN®+ white pd	0.5 mm - 0.75 mm
Etching mask	COPYPLAST®	1.0 mm
Fluorisation splint	DURAN®	1.0 / 1.5 mm
	DURAN®+	1.0 / 1.5 mm
Foil for therapy planning / analysis of bruxism	BRUX CHECKER®	0.1 mm
ImmediateDentures / temporary dentures	BIOCRYL® C (rose or clear)	2.0 / 3.0 mm
Individual impression tray, functional tray	IMPRELON® (clear)	2.0 / 3.0 mm
	IMPRELON® (opaque)	3.0 mm
Insulation Foil	ISOFOLAN®	0.1 mm
Invisible Retainer / Functional orthodontic splint	DURASOFT® pd with DURASOFT® seal	1.2 mm
	DURAN®	0.75 / 1.0 mm
	DURAN®+	0.75 / 1.0 mm
	IMPRELON®S pd	0.75 - 1.0 mm
	COPYPLAST® C	1.0 mm
Medical tray	DURAN®	1.0 mm
	DURAN®+	1.0 mm
Model duplication	BIOPLAST® (plaster)	mind. 3.0 mm
	COPYPLAST® (resin)	2.0 mm
Model protection / Protective coating / Skin-packaging	Coating Foil	0.15 mm
Moulds for temporary crowns and bridges	COPYPLAST®	0.5 - 2.0 mm
Mouth Guard, sports protection	BIOPLAST®	1.0 - 5.0 mm
	BIOPLAST® XTREME / XTREME PRO	4.0 / 5.0 mm
Occlusal splint / TMJ splint / Bruxism splint / Miniplast splint	DURAN®	0.75 - 1.5 mm
	DURAN®+	1.5 mm
	DURASOFT® pd with DURASOFT® seal	1.8 mm
Orthodontic plate / Retainer	BIOCRYL® C	2.0 - 3.0 mm
	BIOCRYL® M	2.0 mm
OSAMU-Retainer® (combination BIOPLAST® and IMPRELON®S pd)	BIOPLAST®	1.5 mm
	IMPRELON®S pd	0.75 mm
Positioner	BIOPLAST®	3.0 - 5.0 mm
Retainer / Orthodontic plates	BIOCRYL® C / M	2.0 mm (for Lower jaw 3.0)
	BIOCRYL®	3.0 mm
Retention splint / Long-term retainer	COPYPLAST® C	1.0 mm
	DURAN®	0.75 / 1.0 mm
	DURAN®+	0.75 / 1.0 mm
Spacer maintainer foil	ISOFOLAN®	0.1 mm
	Spacer Foil	0.1 mm
Spacer maintainer for bleaching splints	HARDCAST®	0.4 / 0.6 / 0.8 mm
Template for implants / Drillings / X-ray	DURAN®	2.0 / 3.0 mm
	DURAN®+	2.0 / 3.0 mm
Temporary Splint / Strip crown for temporary crowns and bridges	DURAN®	0.5 mm - 1.0 mm
	DURAN® +A2 pd	0.5 / 0.75 mm

Appendix

Material Qualities	Bonds to Acrylic	Recommended Acrylic	Recommended Insulation
hard-elastic	yes	DURASPLINT®	ISOFOLAN® foil
hard-elastic	yes	DURASPLINT®	
hard / soft	yes	DURASPLINT®	ISOFOLAN® foil
hard / not elastic	yes	DURASPLINT®	ISOFOLAN® foil
soft-elastic	no		BIOPLAST® Insulating agent / soak model into water
viscoplastic	no		Soak model shortly into water
viscoplastic	no		Soak model shortly into water
viscoplastic	no		Spacer maintainer foil
viscoplastic	no		Soak model shortly into water
hard-elastic	yes	veneering material	
hard-elastic	yes	veneering material	
viscoplastic	no		Soak model shortly into water
hard-elastic	yes		BIOPLAST® (space maintainer)
hard-elastic	yes		
impact-resistant	yes	BIOCRYL® RESIN	ISOFOLAN® foil
hard / not elastic	yes	DURASPLINT®	ISOFOLAN® foil or Foam Discs
hard / not elastic	yes	DURASPLINT®	ISOFOLAN® foil or Foam Discs
viscoplastic insulation foil	no		
hard / soft	yes	DURASPLINT®	ISOFOLAN® foil
hard-elastic	yes	DURASPLINT®	ISOFOLAN® foil
hard-elastic	yes	DURASPLINT®	
hard-elastic	no		ISOFOLAN® foil
viscoplastic	no		Soak model shortly into water
hard-elastic	yes	DURASPLINT®	ISOFOLAN® foil
hard-elastic	yes	DURASPLINT®	
soft-elastic	no		BIOPLAST® Insulating agent / soak model into water
viscoplastic	no		Soak model shortly into water
bonds to Skin cardboard			
viscoplastic	no		Soak model shortly into water
soft-elastic	no		BIOPLAST® insulating agent
soft-elastic	no		BIOPLAST® insulating agent
hard-elastic	yes	DURASPLINT® / DURASPLINT® LC	ISOFOLAN® foil
hard-elastic	yes	DURASPLINT® / DURASPLINT® LC	
hard / soft	yes	DURASPLINT®	ISOFOLAN® foil
impact-resistant	yes	BIOCRYL® RESIN	ISOFOLAN® foil
hard-elastic	yes	BIOCRYL® RESIN	ISOFOLAN® foil
soft-elastic	no		BIOPLAST® insulating agent
hard-elastic	no		
soft-elastic	no		BIOPLAST® insulating agent
impact-resistant / hard elastic	yes	BIOCRYL® RESIN	ISOFOLAN® foil
impact-resistant	yes	BIOCRYL® RESIN	ISOFOLAN® foil
viscoplastic	no		Soak model shortly into water
hard-elastic	yes	DURASPLINT®	ISOFOLAN® foil
hard-elastic	yes	DURASPLINT®	
viscoplastic insulation foil	no		
hard insulation foil	yes		
viscoplastic	no		Spacer maintainer foil
hard-elastic	yes	DURASPLINT®	ISOFOLAN® foil
hard-elastic	yes	DURASPLINT®	
hard-elastic	yes	veneering material	ISOFOLAN® foil
hard-elastic	yes	veneering material	

Notes

Opening and expansion parameters of VECTOR® expansion screws

Expansion screw	REF	Expansion per 1/1 turn	Expansion per ¼ turn	Opening dimension/ total expansion
VECTOR® 30	2484	0,35 mm	0,09 mm	3 mm
VECTOR® 40	2486	0,7 mm	0,18 mm	3 mm
VECTOR® 50	2489	0,7 mm	0,18 mm	3 mm
VECTOR® 80	2490	0,8 mm	0,2 mm	7,5 mm
VECTOR® 90	2491	0,7 mm	0,18 mm	4 mm
VECTOR® 100	2492	0,7 mm	0,18 mm	5 mm
VECTOR® 140	2494	0,8 mm	0,2 mm	5 mm
VECTOR® 160	2496	0,8 mm	0,2 mm	6 mm
VECTOR® 190	2497	0,8 mm	0,2 mm	5 mm
VECTOR® 200	2498	0,8 mm	0,2 mm	7,5 mm
VECTOR® 350	2509	0,35 mm	0,09 mm	4 mm
VECTOR® 400	2510	0,35 mm	0,09 mm	3 mm
VECTOR® 420	2512	0,35 mm	0,09 mm	5 mm
VECTOR® 440	2514	0,35 mm	0,09 mm	5 mm
VECTOR® 440R	2513	0,35 mm	0,09 mm	5 mm
VECTOR® 450	2515	0,35 mm	0,09 mm	4 mm
VECTOR® 460	2516	0,4 mm	0,1 mm	5 mm
VECTOR® 500	2520	0,35 mm	0,09 mm	3 mm
VECTOR® 520	2522	0,35 mm	0,09 mm	5 mm
VECTOR® 720	2530	0,8 / 0,4 mm	0,2/0,1 mm	7 / 5,5 mm
VECTOR® 730	2532	0,4 mm	0,1 mm	5,5 / 2 x 3,5 mm
VECTOR® 740	2540	0,7 / 0,35 mm	0,175 / 0,0875 mm	5 / 2,5 mm
VECTOR® 750	2542	0,35 mm	0,09 mm	2,5 / 2 x 2,5 mm
VECTOR® 900	2536	0,7 mm	0,18 mm	5 mm
VECTOR® 600	2524	0,9 mm	0,23 mm	9 mm
VECTOR® 620	2526	0,9 mm	0,23 mm	11 mm
VECTOR® 800	2533	0,8 mm	0,2 mm	9 mm
VECTOR® 820	2534	0,8 mm	0,2 mm	11 mm

Material list of alloys used by SCHEU-DENTAL (Analysis in Weight-%)

Material	DIN-no.	C	Si	Mn	Cr	Mo	Ni	P	S	Others	Fe	Product Name Product Group
stainless steel	1.4021	0.16-0.25	≤ 1.0	≤ 155	12.0-14.0			≤ 0.040	≤ 0.015			CARAT®
stainless steel	1.4034	0.43-0.5	≤ 1.0	≤ 1.0	12.5-14.5			≤ 0.040	≤ 0.015			CARAT®
stainless steel	1.4116	0.45-0.55	≤ 1.0	≤ 1.0	14.0-15.0	0.5-0.8				V 0.1-0.2		CA-pliers
stainless steel	1.4301	≤ 0.07	≤ 1.0	≤ 2.0	17.5-19.5		8.5-10.5	≤ 0.045	≤ 0.015		Rest	SCHEU-anchor, Clasps/Bars
stainless steel	1.4303	≤ 0.06	≤ 1.0	≤ 2.0	17.0-19.0		11.0-13.0	≤ 0.045	≤ 0.015		Rest	Profi-Bar
stainless steel	1.4305	≤ 0.10	≤ 1.0	≤ 2.0	17.0-19.0		8.0-10.0	≤ 0.045	0.15-0.35		Rest	Adamik-ancor
stainless steel	1.4310	0.05-0.15	≤ 2.0	≤ 2.0	16.0-19.0	≤ 0.8	6.0-9.0	≤ 0.045	≤ 0.03		Rest	CHROMIUM Prefomed Wires
stainless steel	1.4456	≤ 0.1	≤ 1.0	16.0-20.0	16.0-20.0	1.8-2.5	≤ 0.2	≤ 0.050	≤ 0.05	V ≤ 0.2 N 0.7-1.0	Rest	MENZANIUM®

Material	DIN-no.	Ni	Fe	MN	Cr	Mo	Ti	Si	C	Al	O	H	N	V	Cu	Nb	Product Name Product Group
Titan Grade 5	1.4021		0.25	≤ 0.1			Rest		0.08	5.5-6.75	0.2	0.015	0.05	3.5-4.5			CARAT®
CATAMOLD® P.A.N.A.C.E.A.		≤ 0.1	Rest	10-12.0	16.5-17.5	3.0-3.5		≤ 1.0	0.2				0.75-0.9				CARAT®
CATAMOLD® 17-4 PH		3.5	Rest	≤ 0.1	15.0-17.5			≤ 1.0							3.5	0.15-0.45	CA-pliers

Measuring

mm	Inch	°C	°F	DIN-no.	Material-Name	AISI	UNS
0.10	0.004	10	50	1.4021	X20Cr13	420	S42000
0.15	0.006	25	77	1.4034	X46Cr13	420C	
0.40	0.016	50	122	1.4116	X50CrMoV 15		
0.5	0.020	75	167	1.4035	X46CrS13	420F	
0.6	0.024	100	212	1.4301	X5CrNi18-10	304	S30400
0.625	0.025	125	257	1.4303	X4CrNi18-12	305	S30500
0.75	0.030	150	302	1.4305	X8CrNiS18-9	303	S30300
0.8	0.031	180	356	1.4310	X10CrNi18-8	301	S30100
1.00	0.0394	190	374	1.4456	X8CrMnMoN18-18-2		
1.15	0.045	200	392				
1.30	0.051	210	410				
1.50	0.059	220	428				
2.00	0.079	250	482				
3.00	0.118	280	536				
4.00	0.157	300	572				
5.00	0.197	350	662				
25.40	1.00	400	752				
		bar	at			N/mm²	p.s.i.
		1	≈ 1.02			≈ 0.1	14.504
		2	≈ 2.04			≈ 0.2	29.008
		3	≈ 3.06			≈ 0.3	43.511
		4	≈ 4.08			≈ 0.4	58.015
		5	≈ 5.10			≈ 0.5	72.519
		6	≈ 6.12			≈ 0.6	87.023
		7	≈ 7.14			≈ 0.7	101.526

Product Index

REF	Product	Page	REF	Product	Page	REF	Product	Page
1			2334	Face-Bow	43	3137	Sampler of Pressure Moulding Material	109
1002	O-Clasp	16	2335	Face-Bow	43	3138	TWINSTAR P®	110
1003	O-Clasp	16	2336	Face-Bow	43	3139	TWINSTAR V®	110
1007	OK-Clasp	16	2337	Face-Bow	43	3140	TWINSTAR P®	110
1009	N-Clasp	16	2338	Face-Bow	43	3141	TWINSTAR P®	110
1010	N-Clasp	16	2419	Coffin Spring	22	3142	TWINSTAR V®	110
1011	N-Clasp	16	2434	Activator Tubes	44	3143	TWINSTAR V®	110
1012	N BI-Form-Clasp	16	2436	Activator Tubes	44	3147	BIOCRYL® C	82
1013	NK-Clasp	17	2439	Face-Bow	43	3148	BIOCRYL® Form X	82
1017	J-Clasp	18	2450	ORTHOFIT®-Anatomic	42	3149	BIOCRYL® C	82
1018	J-Clasp	18	2451	ORTHOFIT®-Anatomic	42	3150	BIOCRYL® C	82
1019	J-Clasp	18	2453	ORTHOFIT®-Line	42	3151	BIOCRYL® C	82
1020	J-Clasp	18	2454	ORTHOFIT®-Line	42	3152	BIOCRYL® C	82
1021	J-Clasp	18	2455	ORTHOFIT®-Line	42	3153	BIOCRYL® C	82
1022	J-Clasp	18	2456	ORTHOFIT®-Line	42	3164	IMPRELON®	89
1023	JM-Clasp	18	2457	ORTHOFIT®-Line	42	3166	BIOPLAST® Insulating Agent	95
1026	T-Clasp	17	2470	Activator Tubes	44	3168	IMPRELON®	89
1027	T-Clasp	17	2471	Activator Tubes	44	3169	COPYPLAST®	90
1029	TK-Clasp	17	2472	Activator Tubes	44	3175	COPYPLAST®	90
1030	TK-Clasp	17	2477	VECTOR® Piston Screws	36	3177	COPYPLAST®	90
1038	Jaw Fracture Splint	19	2478	VECTOR® Piston Screws	36	3179	COPYPLAST®	90
1039	PROFI Bar	19	2479	VECTOR® Piston Screws	36	3181	COPYPLAST®	90
1040	PROFI Bar	19	2480	Mounting Screw	36	3183	BIOPLAST®	88
1047	J-Clasp	18	2481	Serrated Nut	36	3185	BIOPLAST®	88
1048	J-Clasp oblique	18	2482	Reamer	36	3187	BIOPLAST®	88
1049	Jaw Fracture Splint	19	2483	Screwdriver F	36	3188	BIOPLAST®	88
1051	Lingual Bar	19	2484	VECTOR® 30	30/58	3189	BIOPLAST® Insulating Agent	95
1052	Lingual Bar	19	2486	VECTOR® 40	30/58	3190	BIOPLAST® bleach	88
1053	Lingual Bar	19	2489	VECTOR® 50	30	3194	HARDCAST®	90
2			2490	VECTOR® 80	31	3195	HARDCAST®	90
2040	CARAT® Crampon Pliers	46	2491	VECTOR® 90	31	3200	HARDCAST®	90
2041	CARAT® Flat Pliers	46	2492	VECTOR® 100	31	3202	Spacer Foil	91
2042	CARAT® Young Pliers	46	2494	VECTOR® 140	31	3203	Coating Foil	91
2043	CARAT® Angle Pliers	47	2496	VECTOR® 160	31	3206	SKIN Packaging	91
2044	CARAT® Aderer Pliers large	47	2497	VECTOR® 190	31	3207	ISOFOLAN®	91
2046	CARAT® Aderer Pliers medium	47	2498	VECTOR® 200	31	3209	BRUX CHECKER®	88
2048	SCHEU-Anchor	18 / 22	2509	VECTOR® 350	32	3201	BRUX CHECKER®	88
2049	SCHEU-Anchor	18 / 22	2510	VECTOR® 400	32	3211	Silicone Tube	28
2050	SCHEU-Anchor	18 / 22	2512	VECTOR® 420	32	3214	Cutting Bur	102
2051	SCHEU-Anchor	18 / 22	2513	VECTOR® 440 R	33	3216	Impression Tray Handles	100
2052	SCHEU-Anchor	18 / 22	2514	VECTOR® 440	32	3218	BIOCRYL®-RESIN Monomer	96
2053	SCHEU-Anchor	18 / 22	2515	VECTOR® 450	33	3220	Blocking-out Putty	95
2055	Arrow Clasp	22	2516	VECTOR® 460	33	3221	BIOCRYL®-RESIN Polymer	96
2056	Arrow Clasp	22	2520	VECTOR® 500	33	3222	Stainless Steel Pellets	101
2083	Screwdriver M	36	2522	VECTOR® 520	33	3226	BIOPERM Trimmer	39/102
2100	Security expansion key	37	2524	VECTOR® 600	34	3228	BIOCRYL® M	82
2101	Key for expansion screws	37	2526	VECTOR® 620	34	3229	BIOCRYL® M	82
2104	Adams-Clasp	22	2530	VECTOR® 720	34	3230	BIOCRYL® M	82
2105	Adams-Clasp	22	2532	VECTOR® 730	34	3231	BIOCRYL® M	82
2106	Adams-Clasp	22	2533	VECTOR® 800	35	3232	BIOCRYL® M	82
2107	Adams-Clasp	22	2534	VECTOR® 820	35	3233	BIOCRYL® M	82
2108	Adams-Clasp	22	2536	VECTOR® 900	36	3236	CETRON® Cleansing Set	49
2109	Adams-Clasp	22	2540	VECTOR® 740	35	3238	CETRON® Care Spray	49
2118	U-Bar acc. to Karwetzky	22	2542	VECTOR® 750	35	3239	CETRON® Cleansing Powder	49
2119	Face-Bow	43	3			3242	OSAMU-Retainer® Kit	89
2130	Key for expansion screws	37	3001	BIOSTAR®	108	3246	OSAMU-Bond	98
2288	Triangle Clasp	22	3006	Pellet Cover	101	3247	OSAMU-Polisher	103
2322	Triangle Clasp	22	3007	Pellet Cover	101	3256	DURASPLINT® LC Kit	97
2329	Frontal Guiding Bar acc. to R. Hinz	22	3011	BIOSTAR®	108	3257	DURASPLINT® Kit	96
2330	Frontal Guiding Bar acc. to R. Hinz	22	3021	BIOSTAR®	108	3258	DURASPLINT® Polymer	96
2331	Positioning Support	22	3134	Model Platform	111	3259	DURASPLINT® Monomer	96
2332	Angle Tester acc. to Adamik	22 / 36	3136	Pellet Brush	101	3261	DURASPLINT® LC	97

REF	Product	Page	REF	Product	Page	REF	Product	Page
3263	LC Primer	98	3390	BIOPLAST® MULTICOLOR	87	3703	SD Pipettes	74
3265	LC Bur C	99	3391	BIOPLAST® MULTICOLOR	87	3705	SD Composite Primer	98
3266	LC Insulation	98	3392	BIOPLAST® MULTICOLOR	87	3710	SD Polishing Paste A	99
3267	LC Bur A	99	3393	BIOPLAST® MULTICOLOR	87	3711	SD Polishing Paste B	99
3268	LC Bur B	99	3394	BIOPLAST® MULTICOLOR	87	3712	SD Polishing Brush	99
3270	BIOPLAST®	88	3395	BIOPLAST® MULTICOLOR	87	3713	SD Felt Wheel	99
3271	BIOPLAST® COLOR	87	3396	BIOPLAST® MULTICOLOR	87	3714	SD High-Gloss Polishing Mop	99
3274	IMPRELON®	89	3397	BIOPLAST® MULTICOLOR	87	5		
3283	BIOPLAST® XTREME Kit	86	3398	BIOPLAST® MULTICOLOR	87	5261	CLEAR-BLOKKER®	74/97
3284	BIOPLAST® XTREME Kit	86	3399	BIOPLAST® MULTICOLOR	87	5262	CA® Carbide Bur	62/102
3285	BIOPLAST® XTREME Kit	86	3400	BIOPLAST® MULTICOLOR-Set	87	5263	CA® Grinding Set	62
3286	BIOPLAST® XTREME	86	3403	COPYPLAST® C	90	5264	Cannula	40/94
3287	BIOPLAST® XTREME	86	3405	CA® Foil soft	56	5266	Set-Up Plaster	38
3288	BIOPLAST® XTREME	86	3406	CA® Foil medium	56	5268	Set-Up Wax Sticks	39
3289	BIOPLAST® XTREME	86	3407	CA® Foil hard	56	5270	Model Holder	113
3292	BIOPLAST® XTREME PRO Kit	86	3409	DURAN®	83	5274	CA® CHECKER	54
3293	BIOPLAST® XTREME PRO Kit	86	3411	DURAN®	83	5276	BLUE-BLOKKER®	40/94
3294	BIOPLAST® XTREME PRO Kit	86	3413	DURAN®	83	5277	Cannula	40/94
3295	BIOPLAST® XTREME	86	3415	DURAN®	83	5278	Cannula	40/94
3296	BIOPLAST® XTREME	86	3417	DURAN®	83	5281	CA® TIP 1	59
3297	BIOPLAST® XTREME	86	3419	DURAN®	83	5282	CA® Splint Case for patients	57
3298	BIOPLAST® XTREME	86	3421	DURAN®	83	5283	CA® Splint Bags	57
3312	BIOPLAST® XTREME DECO	87	3422	DURAN® + A2 pd	84	5288	Diamond Separating Disc	40
3313	BIOPLAST® XTREME DECO	87	3423	DURAN® + A2 pd	84	5290	CA® Coloured Pencils	61
3314	BIOPLAST® XTREME DECO	87	3424	DURAN® + A2 pd	84	5291	CA® TIP 2	59
3315	BIOPLAST® XTREME DECO	87	3426	DURAN®+	83	5292	CA® TIP 3	59
3316	BIOPLAST® XTREME DECO	87	3428	DURAN®+	83	5293	Set-Up Separating Disc	38
3317	BIOPLAST® XTREME DECO	87	3430	DURAN®+	83	5295	CA® TIP 4	59
3320	IMPRELON® S pd	89	3432	DURAN®+	83	5298	CA® Book volume 1	62
3322	IMPRELON® S pd	89	3434	DURAN®+	83	5300	CA® Pipettes	61
3324	IMPRELON® S pd	89	3436	DURAN®+	83	5301	CA® Plastic Buttons	61
3326	IMPRELON® S pd	89	3438	CA® Retention foil pd	56	5302	CA® Cyano Veneer Fast	61
3328	IMPRELON® S pd	89	3440	CA® Retention foil pd	56	5303	CA® Travel Case	57
3330	IMPRELON® S pd	89	3441	CA® Foil extra soft	56	5305	CA® Book volume 1	62
3337	DURASOFT® seal	85	3442	SIL-KITT	94	5306	CA® Book volume 1	62
3338	DURASOFT® seal	85	3443	SIL-KITT	94	5307	CA® Book volume 2	62
3339	DURASOFT® seal	85	3444	Foam Discs	100	5308	CA® FLASH	61
3341	DURASOFT® pd	84	3445	BIOPLAST® COLOR	87	5309	CA® Book volume 2	62
3342	DURASOFT® pd	84	3447	MINISTAR® Protection Cover	109	5311	CA® SMART 2D	54
3345	DURASOFT® pd	84	3448	BIOPLAST® COLOR	87	5313	Extension Set (4.0) for CA® Photo Support	55
3347	DURASOFT® blue pd	85	3449	BIOSTAR® Protection Cover	108	5314	Fixing Block for height adjustment of CA® Photo Support	55
3348	DURASOFT® blue pd	85	3452	Heating Support	111	5315	CA® Photo Support	55
3352	BIOPLAST® COLOR	87	3453	CA® Foil A2 soft pd	56	5316	CA® Gauge	63
3353	BIOPLAST® COLOR	87	3454	CA® Foil A2 medium pd	56	5317	CA® Arch CHROMIUM	58
3354	BIOPLAST® COLOR	87	3455	CA® Foil A2 hard pd	56	5318	CA® Book volume 2	62
3355	BIOPLAST® COLOR	87	3460	SD Foil Scissors A	100	5319	CA® Distalizer Kit	57
3356	BIOPLAST® COLOR	87	3461	SD Foil Scissors B	101	5330	CA® POWER GRIP	60
3361	IMPRELON®	89	3470	DURAN®	83	5331	CA® POWER GRIP	60
3369	HM Carbide Cutter	102	3472	DURAN®+	83	5332	CA® POWER GRIP	60
3370	HM Tricutter fine	102	3473	DURAN®+ white pd	84	5333	CA® POWER GRIP	60
3371	POLYFIX on mandrel	39/102	3474	DURAN®+ white pd	84	5334	CA® POWER GRIP	60
3372	POLYFIX	39/102	3475	DURAN®+ white pd	84	5335	CA® POWER GRIP	60
3374	BIOPLAST®	88	3501	MINISTAR S®	109	5336	CA® POWER GRIP	60
3376	DIMO®PRO SLIM	103	3511	MINISTAR S®	109	5337	CA® POWER GRIP	60
3377	Finishing Bur	102	3521	MINISTAR S®	109	5338	CA® POWER GRIP Set	60
3378	Finishing Set	102	3602	SD Bond	92	5339	CA® POWER GRIP	60
3380	DIMO® mandrel	103	3603	SD Silicone	92	5340	CA® POWER GRIP Holder	60
3381	DIMO®PRO	103	3604	SD Scalpel	92	5341	ECHARRI PST Set	63
3382	DIMO®	103	3605	SD Cutter	93	5342	CA® Stripping Tools	63
3383	DIMO®	103	3606	SD Discs	93	5343	CA® Stripping Tools	63
3384	DIMO®	103	3607	SD Bracket Pincers	93	5344	CA® Stripping Tools	63
3387	DURASOFT® pd	84	3608	SD Bracket Remover	93	5345	CA® Stripping Tools	63
3388	DURASOFT® pd	84	3701	SD Cyano Veneer Fast	74	5346	CA® Stripping Tools	63

REF	Product	Page
5347	CA® Stripping Tools	63
5348	CA® Stripping Tools	63
5349	CA® Stripping Tools	63
5350	CA® Stripping Tools	63
5351	CA® Stripping Tools	63
5352	CA® Stripping Tools	63
5353	CA® Stripping Tools	63
5354	CA® Stripping Tools	63
5355	CA® Stripping Tool stainless steel holder	63
5356	CA® Stripping Tool stainless steel holder	63
5357	PST Set Standard	63
5359	Splint Remover	57
5361	CA® Model Holder for Photo Support	55
5362	Parallelometer	113
5368	Set-Up Wax Bars	39
5371	Set-Up Model Trimmer	38
5372	Set-Up Saw	39
5373	Saw Blades	39
5374	Hinge Fixator	112
5375	Mounting Splints	112
5376	Peaked Screws	112
5379	Saw Blades	39
5380	Split-Fixator	112
5435	Retainer boxes	29
5437	SD Box red	62
5438	SD Box green	62
5439	SD Box blue	75
5448	Telescope HF 5	72
5450	Telescope HF 9	72
5452	Telescope HF 15	72
5455	Hinge System acc. to Herbst	41
5456	Hinge System acc. to Herbst	41
5457	Fixing elements	73
5458	Hexagon Socket Screw	73
5461	IST®Kit HF 9	72
5462	GEORGE GAUGE™ bite fork	73
5466	O-Rings	73
5468	Fixing Elements	73
5469	Positioning Support	73
5470	Occlusion Trainer	75
5471	GEORGE GAUGE™ Set	73
5472	GEORGE GAUGE™ bite fork	73
5473	GEORGE GAUGE™ bite fork	73
5478	HERBST hinges	41
5479	Hinge system acc. to Herbst	41
5480	Hinge system acc. to Herbst	41
5481	Distance Rings	41
5482	Distance Rings	41
5483	Distance Rings	41
5484	Distance Rings	41
5485	Distance Rings	41
5489	Fixing Elements for Hinge System	41
5490	Hinge System 44 kit	41
5494	OPM Plastic Bars	75
5495	OPM Kit	75
5606	LC Mini Light	111
5610	LC-6 Light Oven	111
5620	LC-6 Light Oven	111
6		
6501	Asiga MAX	66/106
6502	IMPRIMO® LC Model	67/107
6506	IMPRIMO® LC Impression	67/107
6503	IMPRIMO® LC Splint	67/107

REF	Product	Page
6505	Tray	66/106
6510	3D Model Insulation	28
6512	OnyxCeph3™ CA® SMART 3D - clinic	69
6513	OnyxCeph3™ CA® SMART 3D - laboratory	69
7		
7243	RESIMIX®	28
7244	RESIMIX®	28
7284	RESIMIX®	28
8		
8010	Quad Helix	45
8011	Quad Helix	45
8020	Quad Helix	45
8021	Quad Helix	45
8030	Quad Helix	45
8031	Quad Helix	45
8080	Lingual-/Palatal sheath	45
8081	Lingual-/Palatal sheath	45
8130	STEADY-RESIN	27
8131	STEADY-RESIN	27
8134	STEADY-RESIN	26
8135	STEADY-RESIN	26
8136	STEADY-RESIN	26
8137	STEADY-RESIN	26
8138	Spray Nozzles	29
8139	Spray Nozzles	29
8140	STEADY-RESIN	26
8141	Spray Bottle	29
8142	Ortho Designs	29
8143	Ortho Designs	29
8144	Mica	27
8145	Mica	27
8146	Mica	27
8147	Mica	27
8148	Mica	27
8149	Nozzle for bottles	26
8150	Mica	27
8151	Mica	27
8152	Mica	27
8293	CHROMIUM Wire Rods	25
8294	CHROMIUM Wire Rods	25
8297	CHROMIUM Wire Rods	25
8299	CHROMIUM Wire Rods	25
8306	CHROMIUM Wire Rods	25
8307	CHROMIUM Wire Rods	25
8308	CHROMIUM Wire Rods	25
8310	CHROMIUM Wire Rods	25
8311	CHROMIUM Wire Rods	25
8312	CHROMIUM Wire Rods	25
8313	CHROMIUM Wire Rods	25
8314	CHROMIUM Wire Rods	25
8323	CHROMIUM Coil Wire	25
8324	CHROMIUM Coil Wire	25
8325	CHROMIUM Coil Wire	25
8326	CHROMIUM Coil Wire	25
8327	CHROMIUM Coil Wire	25
8328	CHROMIUM Coil Wire	25
8329	CHROMIUM Coil Wire	25
8330	CHROMIUM Coil Wire	25
8331	CHROMIUM Coil Wire	25
8332	CHROMIUM Coil Wire	25
8333	CHROMIUM Coil Wire	25
8334	CHROMIUM Coil Wire	25
8335	CHROMIUM Coil Wire	25

REF	Product	Page
8336	CHROMIUM Coil Wire	25
8338	CHROMIUM Coil Wire	25
8339	CHROMIUM Coil Wire	25
8350	STEADY-RESIN	27
8351	STEADY-RESIN	27
8352	STEADY-RESIN	27
8353	STEADY-RESIN	27
8354	STEADY-RESIN	27
8355	STEADY-RESIN	27
8356	STEADY-RESIN	27
8357	STEADY-RESIN	27
8358	STEADY-RESIN	27
8359	STEADY-RESIN	27
8360	STEADY-RESIN	27
8361	STEADY-RESIN	27
8364	Orthodontic Insulating Agent	28
8382	CHROMIUM Coil Wire	25
8383	CHROMIUM Coil Wire	25
8431	Symmetry Measurement Gauge	44
8432	ECHARRI Ceph Template	48
8433	ECHARRI Arch Template	48
8434	ECHARRI VTO Template	48
8435	ECHARRI Template Kit	49
8450	Ortho arch bender	47
8453	MENZANIUM® Coil Wire	24
8454	MENZANIUM® Coil Wire	24
8455	MENZANIUM® Coil Wire	24
8456	MENZANIUM® Coil Wire	24
8459	MENZANIUM® Coil Wire	24
8460	MENZANIUM® Coil Wire	24
8461	MENZANIUM® Coil Wire	24
8462	MENZANIUM® Coil Wire	24
8463	MENZANIUM® Coil Wire	24
8464	MENZANIUM® Coil Wire	24
8466	MENZANIUM® Coil Wire	24
8467	MENZANIUM® Wire Rods	24
8468	MENZANIUM® Wire Rods	24
8469	MENZANIUM® Wire Rods	24
8470	MENZANIUM® Wire Rods	24
8471	MENZANIUM® Wire Rods	24
8472	MENZANIUM® Wire Rods	24
8478	MENZANIUM® Wire Rods	24
8480	MENZANIUM® Wire Rods	24
8482	MENZANIUM® Wire Rods	24
9		
9000	TAP® Standard Kit	76
9015	Cover Plate	76
9016	Mounting Device	76
9029	TAP®-T XL Lower guiding device	78
9030	TAP®-T Set	76
9030.2	TAP®-T Set pre-assembled	76
9031	TAP®-T Starter-Kit	77
9035	Activation key	37
9037	TAP®-T Matrix	78
9038	SD Metal Primer	74
9039	Mounting Device	76/77
9050	TAP®-T Reverse Kit	77
9050.2	TAP®-T Reverse Kit pre-assembled	77
9051	TAP®-T Reverse Starter-Kit	77
9070	myTAP™	79
9072	AM Aligner	79
9073	TAP® PAP	79
9082	TAP® Splint case for patients	78



Welcome to SCHEU-DENTAL

SCHEU-DENTAL GmbH
 Am Burgberg 20
 58642 Iserlohn
 Germany



Motorway A45, at Hagener Kreuz
 → A46 Direction Iserlohn

Exit 46
 → Iserlohn – Oestrich

Please turn left after 200 m
 → Direction Oestrich

Please turn right at the roundabout
 → Gewerbegebiet Burgberg

Please turn right after 600 m
 → Schleddestraße

Please turn right immediately again
 → Am Burgberg



ALL GOOD THINGS COME IN THREES.

CA DIGITAL **YOUR PARTNER FOR DIGITAL ORTHODONTIC TREATMENT**

CA DIGITAL's **CA® CLEAR ALIGNER** is simply three times as good:
3 innovative material thickness grades for up to 1 mm tooth movement per step, **3 steps per impression** for more precise results - all of this online in just **3 clicks**. Thanks to innovative software solutions, integrated workflows and advanced treatment options, CA DIGITAL is the right partner for all practices that intend taking the path to digitization. CA DIGITAL - your digital future is just a click away: www.ca-digit.com





More information:

SCHEU-DENTAL GmbH
www.scheu-dental.com

phone +49 2374 9288-0
fax +49 2374 9288-90

Am Burgberg 20
58642 Iserlohn · Germany

