

# Splint Remover

Using the Splint Remover, the patient can easily remove pressure-moulded aligners or splints. Please proceed as follows:



1. Hook the Splint Remover into the rear edge of the aligner or splint. Then begin with loosening the splint at one side with a slight downwards movement
2. Proceed in the same way at the other side.
3. Repeat these steps (alternating between right and left side) until the aligner/splint can be easily removed in order to make sure not to damage it.

Non-sterile.


Ref: 5359.1 / 5 pcs.



SCHEU-DENTAL GmbH  
www.scheu-dental.com

phone +49 2374 9288-0  
fax +49 2374 9288-90

Am Burgberg 20  
58642 Iserlohn · Germany

 Facebook  
scheu.dental

**SCHEU**  
Dental Technology

DE/GB 5.000/02/19 G REF PM0175.01

