



-45%

-50%

-60%

BLACK DAYS

promotions

S L O V A K I A

.....

Valid from 20th until 24th
November 2023

INDEX

Personal Protection Equipment

ULTRA™ EARLOOP BLACK FACEMASK
with Secure Fit™ Mask Technology P.1

Diagnostic

DRESSING PLIERS P.2

MIRRORS P.2

PROBES P.2

EXPLORERS P.2

Hand Scaling

CURETTES AND SCALERS
EVEREDGE™ 2.0 P.3

Restorative

COMPOSCULP INSTRUMENT KIT P.3

PAPER FORCEPS P.3

EXCAVATORS P.3

CROWN REMOVERS P.4

PRIMARY CROWNS P.4

Endodontic

MEASURING DEVICES P.5

ENDODONTIC MIRRORS P.5

LOCKING PLIERS P.5

POST & SILVER POINT
REMOVAL FORCEPS P.5

ENDODONTIC EXPLORERS P.5

RUBBER DAM FRAME P.5

RUBBER DAM PUNCH P.5

Surgical

SCALPEL BLADES AND HANDLES P.6

CHEEK RETRACTORS P.6

PERIOSTEALS P.6

SCISSORS P.7

SURGICAL CURETTES P.8

TISSUE PLIERS P.8

NEEDLE HOLDERS P.9

PERSONAL PROTECTION EQUIPMENT CROSSTEX

-60%

ULTRA™ EARLOOP BLACK FACEMASK with Secure Fit™ Mask Technology

FEATURE & BENEFIT

Secure Fit™ Technology

The innovative Crosstex™ Ultra™ Earloop Facemask with patented Secure Fit™ Technology creates a snug, custom fit that adapts to most size or shape face to significantly minimizes gapping and provides high level of protection.

Fluid resistant outer layer

Spunbond outer layer material for maximum protection against fluids traveling through the mask.

Fluid resistant white inner layer

Unique layer barrier for maximum protection against fluids.

Extra-long aluminium nose piece

Significantly reduces gapping for maximum protection.

Not made with natural rubber latex

Prevents latex allergy

Fluid resistance of 160 mmHg

Ideal for procedures where high amounts of fluid, spray and/or aerosols are produced.



GCFBKSF

Ultra™ Earloop Facemask Secure Fit™ Technology, Black, (EN14683 Type IIR Enhanced), 50 pcs/box

Retail 20,40 €

Promo 8,16 €

-50%

DRESSING PLIERS



DP18
#18

Retail 33,93 €
Promo 16,96 €

DP2
#2 College

Retail 30,96 €
Promo 15,48 €

MIRRORS



MIR5DSHDE
#5 HD Double-Sided Mirror, EU Simple Stem, single

Retail 23,48 €
Promo 11,74 €

MHE6
#6 European Simple Stem Mirror Handle

Retail 30,90 €
Promo 15,45 €

PROBES

1-2-3-5-7-8-9-10 mm

3-6-8-11 mm



PQW
Qulix SE Williams Probe, handle #30

Retail 36,59 €
Promo 18,29 €



PCP116
Qulix SE Probe, handle #6

Retail 48,33 €
Promo 24,16 €

EXPLORERS



EXD11/126
#11/12 ODU DE Explorer handle #6

Retail 39,43 €
Promo 19,71 €

HAND SCALING HU-FRIEDY

RESTORATIVE HU-FRIEDY

CURETTES AND SCALERS EVEREDGE™ 2.0



-45%

Valid for the whole portfolio of EverEdge™ 2.0 Scalers and Curettes handle #9

SG5/695E2

Gracey 5/6 Curette, Standard, handle #9

Retail 77,37 €

Promo 42,56 €

SG13/1498E2

Gracey 13/14 Curette, Standard, handle #9

Retail 77,37 €


Promo 42,56 €

COMPOSCULP INSTRUMENT KIT

SMOOTH HANDLES

-45%

CompoSculp instruments were designed by Dr. Didier Dietschi with the intention of aiding clinicians in achieving superior esthetic restorations. The CompoSculp kit is comprised of 5 complimentary composite instruments for both anterior and posterior restorations.



Contains:

Dietschi CompoSculp 1/2, #CCH Handle	PFIDD1/2
Dietschi CompoSculp 3/4, #CCH Handle	PFIDD3/4
Dietschi CompoSculp 5/6, #CCH Handle	PFIDD5/6
Dietschi CompoSculp 7/8, #CCH Handle	PFIDD7/8
Dietschi CompoSculp 9/10, #CCH Handle	PFIDD9/10
IMS™ Infinity Series™ 1/4 DIN Cassette, 5 Instruments, Orange	IME14DIN53

PFIDDDIN

CompoSculp Instrument Kit, Smooth Handles #CCH

Retail 444,22 € each

Promo 244,32 € each

PAPER FORCEPS

-60%



APF2E

Miller Articulating Paper Forceps

Retail 57,97 €

Promo 23,19 €

EXCAVATORS



-50%

EXC18H

18 Heavy Excavator, handle #41

Retail 53,09 €

Promo 26,54 €

RESTORATIVE HU-FRIEDY

CROWN REMOVERS



-40%

CRU

Trial Crown
Remove

Retail 328,78 €

Promo 197,27 €

CRL

Trial Crown
Remover, Lower

Retail 328,78 €

Promo 197,27 €

CRCH1

Christensen Crown
Remover Straight,
Anterior

Retail 119,21 €

Promo 71,53 €

PRIMARY CROWNS

Stainless Steel Pedo Crowns come in 48 sizes for use on the first and second primary molars. These crowns feature a more narrow mesial-distal maxillary first primary molar width for quicker adjacent crown placement.

-45%



Upper Right Primary 1st #5, 5 pcs/box	SSC-URD5-CE
Upper Right Primary 2nd #3, 5 pcs/box	SSC-URE3-CE
Upper Left Primary 1st #3, 5 pcs/box	SSC-ULD3-CE
Upper Left Primary 1st #4 5 pcs/box	SSC-URD4-CE

Retail 30,65 € each

Promo 16,86 € each


ENDODONTIC HU-FRIEDY

DID YOU KNOW? All hinged instruments – forceps, rongeurs, scissors, pliers, hemostats, orthodontic pliers, etc. – should be kept lubricated. Regular use of proper lubricants, like HuFriedyGroup Instrument Lubricant Spray (ILS), will prevent rust, corrosion and stiff joints to ensure smooth operation. **All hinged instruments should be sterilised in the open position.**



MEASURING DEVICES

-45%



CLR6

Gracey 13/14 Curette, Standard, handle #9

Retail 45,24 €

Promo 24,88 €

ENDODONTIC MIRRORS

-45%



MM5

Micro Mirror Round, 5 mm

Retail 124,40 €

Promo 68,42 €

-45%



MM8

Micro Mirror, 3 mm x 9 mm

Retail 153,14 €

Promo 84,23 €

LOCKING PLIERS

-50%



EPL1

1 Endo Locking Pliers, 160-5 mm

Retail 57,04 €

Promo 28,52 €

POST & SILVER POINT REMOVAL FORCEPS

-50%



RF45

Steiglitz (Archer) 45°, 130-5 mm

Retail 155,12 €

Promo 77,56 €

ENDODONTIC EXPLORERS

-50%



EXDG16

#DG16 DE, Endodontic Explorer, handle #31

Retail 45,48 €

Promo 22,74 €

RUBBER DAM FRAME

-50%

RDAF6

Rubber Dam Frame, Adult - 15 cm

Retail 48,39 €

Promo 24,19 €

RUBBER DAM PUNCH

-50%



RDPA

Rubber Dam Punch, Ainsworth, 17,5 cm

Retail 302,33 €

Promo 151,16 €

SCALPEL BLADES AND HANDLES

-60%



10-255-11
#11 Sterile Scalpel
Blades, 100 pcs/box,
Carbon Steel

Retail 57,04 €

Promo 22,82 €

-45%



10-130-05E
#5 Scalpel Handle

Retail 104,32 €

Promo 57,38 €

-60%



10-256-12
#12D Sterile Scalpel
Blades, 100 pcs/box,
Carbon Steel

Retail 68,17 €

Promo 27,27 €

CHEEK RETRACTORS

-45%



CRO2
Orringer Retractor
Medium

Retail 62,29 €

Promo 34,26 €

CRW2
Weider Cheek & Tongue
Retractor, Large

CRW1
Weider Cheek & Tongue
Retractor, Small

Retail 107,10 € each

Promo 58,90 € each


CRM2
Cawood-Minnesota
Retractor

Retail 47,90 €

Promo 26,34 €

PERIOSTEALS

-45%



PWILL
Williger Periosteal,
Large

Retail 99,37 €

Promo 54,66 €

PP25
#25 Palti Periosteal
Large

Retail 125,39 €

Promo 68,97 €

DID YOU KNOW? All hinged instruments – forceps, rongeurs, scissors, pliers, hemostats, orthodontic pliers, etc. – should be kept lubricated. Regular use of proper lubricants, like HuFriedyGroup Instrument Lubricant Spray (ILS), will prevent rust, corrosion and stiff joints to ensure smooth operation. **All hinged instruments should be sterilised in the open position.**



SCISSORS

SUPER-CUT

Cutting blade honed/ serrated to „razor“ sharpness with opposing blade finished with fine serrations for reduced tissue slippage.

PERMA SHARP™

Featuring tungsten carbide inserts bonded to surgical-grade stainless steel on both cutting edges for smooth, controlled cutting and maximum durability.

-60%



S14

#14 LaGrange Scissors,
11,5 cm, 1 serrated blade,
Standard

Retail 88,62 €

Promo 35,45 €

-60%



S16

#16 Iris Scissors, 12,5
cm, 1 serrated blade,
Standard

Retail 66,37 €

Promo 26,55 €

-45%



S5082

Iris Perma Sharp™
Scissors, Straight,
11,5 cm

Retail 238,36 €

Promo 131,10 €

-45%



SCORT

Scissors Cortellini,
Curved, 12 cm

Retail 79,97 €

Promo 43,98 €

-45%



SCORTSC

Scissors Cortellini,
Super-Cut, Curved,
12 cm, 1 serrated blade

Retail 144,12 €

Promo 79,26 €

SURGICAL CURETTES

-45%



CL88

#88 Lucas DE Surgical Curette, handle #522

Retail 116,68 €

Promo 64,17 €

-45%



CM46

#4 Molt SE Curette, handle #6

Retail 101,23 €

Promo 55,68 €

TISSUE PLIERS

-45%



TP42

#42 Adson 1x2 Tissue Pliers, 12 cm

Retail 40,48 €

Promo 22,26 €



TP5043

Adson-Brown Perma Sharp™ Tissue Plier, 12 cm

Retail 243,68 €

Promo 134,02 €



TPG1

Gerald Tissue Pliers, Straight/Serrated, 17,78 cm

Retail 62,91 €

Promo 34,60 €



TPDAPV

Anatomic Microsurgical Tissue Plier, 18 cm

Retail 472,40 €

Promo 259,82 €



TPDSSPV

Microsurgical Tissue Plier, Straight, 18 cm

Retail 602,12 €

Promo 331,16 €

DID YOU KNOW? All hinged instruments – forceps, rongeurs, scissors, pliers, hemostats, orthodontic pliers, etc. – should be kept lubricated. Regular use of proper lubricants, like HuFriedyGroup Instrument Lubricant Spray (ILS), will prevent rust, corrosion and stiff joints to ensure smooth operation. **All hinged instruments should be sterilised in the open position.**



NEEDLE HOLDERS

1-0, 2-0, 3-0, 4-0



NH5072

Olsen-Hegar Perma Sharp™ Needle Holder, 14 cm

Retail 249,36 €

Promo 137,15 €



NHLSCHLEE

Schlee Microsurgical Needle Holder

Retail 947,27 €

Promo 521,00 €

5-0, 6-0, 7-0, 8-0, 9-0



SPNHDPVN

Swiss Perio Microsurgical Needle Holder, 18 cm

Retail 978,97 €

Promo 538,44 €

6-0, 7-0, 8-0, 9-0, 10-0




NHTMBH

Titanium Diamond Dusted Castroviejo Microsurgical Needle Holder, 18 cm

Retail 1.548,40 €

Promo 851,62 €

6-0, 7-0, 8-0



NHDPV

Diamond Dusted Castroviejo Microsurgical Needle Holder, Straight, 18 cm

Retail 981,14 €

Promo 539,63 €

1-0, 2-0, 3-0, 4-0



NHMH

Mayo-Hegar Needle Holder, 16 cm

Retail 64,83 €

Promo 35,66 €

6-0, 7-0, 8-0, 9-0, 10-0, 11-0



NHM-5027

Castroviejo Microsurgical Needle Holder, Curved, 13 cm

Retail 517,64 €

Promo 284,70 €

1-0, 2-0, 3-0, 4-0




NH5074

Mathieu-Kocher Perma Sharp™ Needle Holder, 14 cm

Retail 253,19 €

Promo 139,26 €



For more information on HuFriedyGroup products or to place an order,
call at 00800 48 37 43 39 or visit **HuFriedyGroup.eu**

As some of the products included in this catalogue might not be available in some countries, we invite you to verify the availability of the product with your local authorised dealer or to contact our customer care at the free phone number 00800 48 37 43 39.

Hu-Friedy Mfg. Co., LLC • European Headquarters • Lyoner Str. 9 • 60528 Frankfurt am Main, Germany • HuFriedyGroup.eu
All company and product names are trademarks of Hu-Friedy Mfg. Co. LLC, its affiliates or related companies, unless otherwise noted.
©2023 Hu-Friedy Mfg. Co., LLC All rights reserved. HFL-072/1123