



€166 **€ 116**



€63,60
€ 31,80

€104
€ 83



HU-FRIEDY NATIONAL PROMOTIONS

SLOVAKIA
AUTUMN 2019

Valid from 1. October until 31. December 2019


VAT included

For more information about Hu-Friedy's Total Solutions visit HU-FRIEDY.EU

©2019 Hu-Friedy Mfg. Co., LLC. All rights reserved.

How the best perform

Hu-Friedy



INDEX

INSTRUMENT MANAGEMENT SYSTEM

- 01
- 1 • IMS Infinity Cassettes
 - 2 • Enzymax and Sporeclear

DIAGNOSTIC

- 03
- 3 • Probes, Expos and Explorers
 - 4 • Mirrors and Mirror Handles

ENDODONTIC

- 05
- 5 • Endo Accessories

PERIODONTAL

- 06
- 6 • Titanium Scalers
 - 7 • EverEdge 2.0 Scalers and Curettes

POWER

- 09
- 9 • Ultrasonic Tips

ORTHODONTIC

- 10
- 10 • Cutters and Pliers

SURGICAL


- 11
- 11 • Needle Holders
 - 12 • Scissors and Surgical Curettes
 - 13 • Elevators and Periosteals
 - 14 • Tunneling Knives and Tissue Pliers
 - 15 • Scalpel Blades and Handles

RESTORATIVE

- 16
- 16 • XTS Instruments




Visit our new **IMS INSTRUMENT MANAGEMENT SYSTEM** site
Organisation, Efficiency and Safety in a simple click!
 Use the QRCode



IMS INFINITY CASSETTES

INFECTION PREVENTION WITH EFFICIENT SAFETY

Optimal infection control and prevention are facilitated by the use of tools that allow for improved efficiency and safety. IMS allows clinicians to move instruments from cleaning through chairside without touching or damaging the instruments.



1/4 DIN CASSETTE

IME14DIN5_
 IMS Infinity Cassette
 1/4 DIN, 5 Instruments
 Dimensions: 72 x 184 x 34 mm.

Sugg. Retail Price

€ 140

Promo Price

€ 98

1/2 DIN CASSETTE

IME12DIN1_
 IMS Infinity Cassette 1/2 DIN, 8
 Instruments with an accessory area.
 Ideal for diagnostics and prophylaxis.
 Dimensions: 152 x 184 x 34 mm

Sugg. Retail Price

€ 158

Promo Price

€ 110

DIN CASSETTE

IMEDIN18_
 IMS Infinity Cassette DIN,
 18 Instruments with accessory area.
 Dimensions: 290 x 184 x 34 mm

Sugg. Retail Price

€ 212

Promo Price

€ 149

DIN SURGICAL CASSETTE

IMEDIN18_S
 IMS Infinity Cassette DIN,
 18 Instruments Surgery + 4 clips
 (2 x 11000) for hinged instruments.
 Dimensions: 290 x 184 x 34 mm

Sugg. Retail Price

€ 228

Promo Price

€ 160

SIGNATURE SERIES TUBS

Together, procedure tubs and cassettes create the most efficient and organised system for managing instruments and consumable products. The time savings achieved by using IMS (Instrument Management System) cassettes and procedure tubs will allow the staff to spend time on more value added activities in the dental office.

IMS-140_

Includes:

- Signature Series Tub
- Tub Liner
- Tub Divider (package of 8)
- Single Tub Cup w/Cover
- Double Tub Cup w/Cover
- Divided Slide Tray
- Locking Tub Cover



Sugg. Retail Price

€ 88


Promo Price

€ 70

CHOOSE THE COLOR.

Just add the color coding number to the last part of the partcode (replace "_" with the number)

Red	0	Lavender	6
Grey	1	Black	7
White	2	Blue	8
Orange	3	Green	9
Purple	4	Ocean Blue	11
Yellow	5		



SPORECLEAR SURFACES DISINFECTANT

It offers a fast, simple and effective answer to creating a safe work environment protecting dental staff and their patients from harmful pathogens. **Sporeclear** is **non-alcohol based, free of aldehydes, chlorine, phenols, peroxide** and **nonylphenol ethoxylates**, (NPEs), making them **compatible with all washable surfaces**.

Fast solutions for effective disinfection
SPORICIDAL - TUBERCULOCIDAL - BACTERICIDAL - VIRUCIDAL - FUNGICIDAL



MSCD354160
Wipes, 160 wipes canister

Sugg. Retail Price
€ 13

Promo Price
€ 11



MSCD354161
160 wipes refill pack

Sugg. Retail Price
€ 11

Promo Price
€ 9



MSCD356750
Liquid Spray, 750 ml bottle
(ready-to-use)

Sugg. Retail Price
€ 16

Promo Price
€ 12

ENZYMAX DUAL ENZYME ULTRASONIC DETERGENT/PRESOAK

A superior level of instrument cleaning powered by a dual-enzyme design that works quickly to remove and prevent biofilms while thoroughly cleaning instruments. Extended instrument life enabled by Steelgard agents and corrosion inhibitors.



IMS-1224
IMS Enzymax Liquid,
Quart Bottle,
Makes 121 lts.

Sugg. Retail Price
€ 60

Promo Price
€ 39



IMS-1226
IMS Enzymax Liquid,
Gallon,
Makes 485 lts.

Sugg. Retail Price
€ 121

Promo Price
€ 79



ENZYMAX PRECLEANER

Enzymax Spray Gel: Dual enzyme gel with corrosion inhibitor acts as a **ready-to-use instrument pre-cleaner**. Maintains moisture to help prevent debris from drying on instruments before cleaning and **helps dissolve blood, tissue, and other debris**. Can also be used as a linen spot remover.

IMS-1229
IMS Enzymax Spray Gel,
709 ml,
up to 450 Sprays

Sugg. Retail Price
€ 28

Promo Price
€ 18

PROBES



PCPUNC156
#CP-15 UNC Color-Coded, handle #6

Sugg. Retail Price
€ 40

Promo Price
€ 28

EXPROS



XP23/86
#23/CP-8 Color-Coded, handle #6



XP23/11.56
#23/CP-11.5B Color-Coded

Sugg. Retail Price
€ 47

Promo Price
€ 33

EXPLORERS



EXS236
#23 Shepherd's Hook, handle #6

Sugg. Retail Price
€ 23

Promo Price
€ 16



EXS96
#9, handle #6

Sugg. Retail Price
€ 23

Promo Price
€ 16



EXD57
#5, handle #7

Sugg. Retail Price
€ 34

Promo Price
€ 24





MIRRORS AND MIRROR HANDLES



MIR3HDE
#3 European
Single-Sided, Single
Ø 20 mm

Sugg. Retail Price

€ 15

Promo Price

€ 11



MIR4HDE
#4 European
Single-Sided, Single
Ø 22 mm

Sugg. Retail Price

€ 15

Promo Price

€ 11



MIR5HDE
#5 European
Single-Sided, Single
Ø 24,1 mm

Sugg. Retail Price

€ 15

Promo Price

€ 11



MHE6
Mirror Handle #6,
European

Sugg. Retail Price

€ 26

Promo Price

€ 17



ENDO ACCESSORIES




RF45
45° Steiglitz (Archer) Post and Silver
Point Removal Forceps 13,5 cm (5-1/4")

Sugg. Retail Price

€ 128

Promo Price

€ 77



RF90
45° Steiglitz (Archer) Post and Silver
Point Removal Forceps 13,5 cm (5-1/4")

Sugg. Retail Price

€ 128

Promo Price

€ 77



DID YOU KNOW?

All hinged instruments – forceps, rongeurs, scissors, pliers, hemostats, orthodontic pliers, etc. – should be kept lubricated. Regular use of proper lubricants, like Hu-Friedy Penetrating Oil (IPS) and Instrument Lubricant Spray (ILS), will prevent rust, corrosion and stiff joints to ensure smooth operation. All hinged instruments should be sterilised in the open position.

IPS
IMS Penetrating
Oil, 8 oz (236 ml)

Retail **Promo**

€ 31 € 22

ILS
IMS Instrument
Lubricant Spray,
8 oz (236 ml)

Retail **Promo**

€ 31 € 22





TITANIUM IMPLANT SCALERS



IMEDIMPLM
IMS Infinity Cassette with
5 Titanium Instruments

Sugg. Retail Price

€ 440

Promo Price

€ 352



IMPLG1/2T
Langer 1/2

Sugg. Retail Price

€ 75

Promo Price

€ 60



IMPM1/2T
Mini Five 1/2

Sugg. Retail Price

€ 75

Promo Price

€ 60



IMP204SDT
204SD

Sugg. Retail Price

€ 75

Promo Price

€ 60



IMPM11/12T
Mini Five 11/12

Sugg. Retail Price

€ 75

Promo Price

€ 60



IMPM13/14T
Mini Five 13/14

Sugg. Retail Price

€ 75

Promo Price

€ 60

EVEREDGE® 2.0 THE MOST ADVANCED SCALER EVER

SAVE 50% OFF BY HU-FRIEDY!

RETURN AT LEAST 5 OLD SCALERS AND CURETTES FROM ANY OTHER BRAND AND **GET A 50% DISCOUNT ON THE SAME QUANTITY OF NEW HU-FRIEDY EVEREDGE 2.0 SCALER AND CURETTE.**

Contact your local representative for more info and instructions on how return your old instruments and **see the next page to choose your new instrument between the complete EverEdge 2.0 range ▶**

*OFFER VALID ON THE COMPLETE EVEREDGE 2.0 SCALER AND CURETTES RANGE WITH #9 AND #RESIN8 HANDLES.





EVEREDGE® 2.0

THE MOST ADVANCED SCALER EVER

STANDARD GRACEY CURETTES

Description	Handle 9	Resin 8
1/2 Gracey Standard	SG1/291E2	SG1/2C8E2
3/4 Gracey Standard	SG3/499E2	SG3/4C8E2
5/6 Gracey Standard	SG5/695E2	SG5/6C8E2
7/8 Gracey Standard	SG7/897E2	SG7/8C8E2
9/10 Gracey Standard	SG9/1096E2	SG9/10C8E2
11/12 Gracey Standard	SG11/1293E2	SG11/12C8E2
13/14 Gracey Standard	SG13/1498E2	SG13/14C8E2
15/16 Gracey Standard	SG15/1693E2	SG15/16C8E2
17/18 Gracey Standard	SG17/1898E2	SG17/18C8E2
11/14 Gracey Standard	SG11/149E2	SG11/14C8E2
12/13 Gracey Standard	SG12/139E2	SG12/13C8E2
1/2 Gracey Standard, Rigid	SG1/2R9E2	SG1/2RC8E2
3/4 Gracey Standard, Rigid	SG3/4R9E2	SG3/4RC8E2
5/6 Gracey Standard, Rigid	SG5/6R9E2	SG5/6RC8E2
7/8 Gracey Standard, Rigid	SG7/8R9E2	SG7/8RC8E2
9/10 Gracey Standard, Rigid	SG9/10R9E2	-
11/12 Gracey Standard, Rigid	SG11/12R9E2	SG11/12RC8E2
13/14 Gracey Standard, Rigid	SG13/14R9E2	SG13/14RC8E2
15/16 Gracey Standard, Rigid	SG15/16R9E2	SG15/16RC8E2
17/18 Gracey Standard, Rigid	SG17/18R9E2	SG17/18RC8E2

AFTER FIVE GRACEY CURETTES

Description	Handle 9	Resin 8
1/2 After Five Gracey	SRPG1/291E2	SRPG1/2C8E2
3/4 After Five Gracey	SRPG3/499E2	SRPG3/4C8E2
5/6 After Five Gracey	SRPG5/695E2	SRPG5/6C8E2
7/8 After Five Gracey	SRPG7/897E2	SRPG7/8C8E2
11/12 After Five Gracey	SRP11/1293E2	SRP11/12C8E2
13/14 After Five Gracey	SRP13/1498E2	SRP13/14C8E2
15/16 After Five Gracey	SRP15/1693E2	SRP15/16C8E2
11/14 After Five Gracey	SRPG11/149E2	SRPG11/14C8E2
12/13 After Five Gracey	SRPG12/139E2	SRPG12/13C8E2
1/2 After Five Gracey, Rigid	SRP1/2R9E2	SRP1/2RC8E2
3/4 After Five Gracey, Rigid	SRP3/4R9E2	SRP3/4RC8E2
5/6 After Five Gracey, Rigid	SRP5/6R9E2	SRP5/6RC8E2
7/8 After Five Gracey, Rigid	SRP7/8R9E2	SRP7/8RC8E2
11/12 After Five Gracey, Rigid	SRP11/12R9E2	SRP11/12RC8E2
13/14 After Five Gracey, Rigid	SRP13/14R9E2	SRP13/14RC8E2
15/16 After Five Gracey, Rigid	SRP15/16R9E2	SRP15/16RC8E2

MINI FIVE GRACEY CURETTES

Description	Handle 9	Resin 8
1/2 Mini Five Gracey	SAS1/291E2	SAS1/2C8E2
3/4 Mini Five Gracey	SAS3/49E2	SAS3/4C8E2
5/6 Mini Five Gracey	SAS5/695E2	SAS5/6C8E2
7/8 Mini Five Gracey	SAS7/897E2	SAS7/8C8E2

MINI FIVE GRACEY CURETTES

Description	Handle 9	Resin 8
11/12 Mini Five Gracey	SAS11/1293E2	SAS11/12C8E2
13/14 Mini Five Gracey	SAS13/1498E2	SAS13/14C8E2
15/16 Mini Five Gracey	SAS15/1693E2	SAS15/16C8E2
17/18 Mini Five Gracey	SAS17/1898E2	-
11/14 Mini Five Gracey	SAS11/149E2	SAS11/14C8E2
12/13 Mini Five Gracey	SAS12/139E2	SAS12/13C8E2
1/2 Mini Five Gracey, Rigid	SAS1/2R9E2	SAS1/2RC8E2
3/4 Mini Five Gracey, Rigid	SAS3/4R9E2	SAS3/4RC8E2
5/6 Mini Five Gracey, Rigid	SAS5/6R9E2	SAS5/6RC8E2
7/8 Mini Five Gracey, Rigid	SAS7/8R9E2	SAS7/8RC8E2
11/12 Mini Five Gracey, Rigid	SAS11/12R9E2	SAS11/12RC8E2
13/14 Mini Five Gracey, Rigid	SAS13/14R9E2	SAS13/14RC8E2
15/16 Mini Five Gracey, Rigid	SAS15/16R9E2	SAS15/16RC8E2

MICRO MINI FIVE GRACEY CURETTES

Description	Handle 9	Resin 8
1/2 Micro Mini Five Gracey	SMS1/291E2	-
7/8 Micro Mini Five Gracey	SMS7/897E2	-
11/12 Micro Mini Five Gracey	SMS11/1293E2	-
13/14 Micro Mini Five Gracey	SMS13/1498E2	-

BIOGENT CURETTES

Description	Handle 9	Resin 8
1/2 Biogent	SG1/291E2	-
7/8 Biogent	SG7/897E2	-
11/12 Biogent	SG11/1293E2	-
13/14 Biogent	SG13/1498E2	-

UNIVERSAL CURETTES

Description	Handle 9	Resin 8
1/2 Barnhart	SBH1/29E2	SBH1/2C8E2
5/6 Barnhart	SBH5/69E2	SBH5/6C8E2
1/2 Barnhart, USC/Suter	SBH1/2USC9E2	-
5/6 Barnhart, USC/Suter	SBH5/6USC9E2	-
1/2 Barnhart, Rigid	SBH1/2R9E2	SBH1/2RC8E2
5/6 Barnhart, Rigid	SBH5/6R9E2	SBH5/6RC8E2
2R/2L Columbia	SC2R/2L9E2	SC2R/2LC8E2
4R/4L Columbia	SC4R/4L9E2	SC4R/4LC8E2
13/14 Columbia	SC13/149E2	SC13/14C8E2
13/14 McCall	SMT13/149E2	SMT13/14C8E2
YG7/B Younger-Good	SYG7/89E2	SYG7/8C8E2
1/2 Langer	SL1/29E2	SL1/2C8E2
3/4 Langer	SL3/49E2	SL3/4C8E2
5/6 Langer	SL5/69E2	SL5/6C8E2
17/18 Langer	SL17/189E2	SL17/18C8E2

UNIVERSAL CURETTES

Description	Handle 9	Resin 8
1/2 After Five Langer	SL1/2AF9E2	-
3/4 After Five Langer	SL3/4AF9E2	-
5/6 After Five Langer	SL5/6AF9E2	-
1/2 Mini Five Langer	SL1/2MF9E2	-
3/4 Mini Five Langer	SL3/4MF9E2	-
5/6 Mini Five Langer	SL5/6MF9E2	-

SICKLE SCALERS

Description	Handle 9	Resin 8
204S Scaler	S204S9E2	S204SC8E2
204SD Scaler	S204SD9E2	S204SDC8E2
H5/33 Scaler	SH5/339E2	SH5/33C8E2
H6/H7 Scaler	SH6/79E2	SH6/7C8E2
204 IUFW Scaler	SIUFW2049E2	SIUFW204C8E2
30/33 Jacquette	SJ30/339E2	SJ30/33C8E2
31/32 Jacquette	SJ31/329E2	SJ31/32C8E2
34/35 Jacquette	SJ34/359E2	SJ34/35C8E2
0/00 Morse	SM0/009E2	SM0/00C8E2
13S/14S McCall	SM13/14S9E2	SM13/14SC8E2
17/18 McCall	SM17/189E2	SM17/18C8E2
17/18 McCall, Ptd Tip	SM17/18PT9E2	-
23 Scaler	SM239E2	SM23C8E2
23A Scaler	SM23A9E2	SM23AC8E2
128 Nebraska	SM128D9E2	SM128DC8E2
135 Scaler	SM1359E2	SM135C8E2
137 Scaler	SM1379E2	SM137C8E2
U15/30 Scaler	SU15/309E2	SU15/30C8E2
U15/33 Scaler	SU15/339E2	SU15/33C8E2
1 Nevi Anterior	SCNEV19E2	SCNEV19C8E2
2 Nevi Posterior	SCNEV29E2	SCNEV29C8E2
3 Nevi Scaler Posterior	SCNEV39E2	SCNEV39C8E2
4 Nevi Scaler Posterior	SCNEV49E2	SCNEV49C8E2
1 Nevi H5, Anterior	SNEV1/H59E2	SNEV1/H5C8E2

	Handle 9	€ 63,60	Promo Price € 31,80
	Resin 8	€ 60	

ULTRASONIC TIPS

SERIE S

BUY 2 AND GET 1 FOR FREE

Compatible with Satelec/Acteon* generators.

TOP-SELLER

POWER SETTING MEDIUM

S10 Universal.
Standard diameter tip used to remove heavy supragingival deposits.

USS10

Sugg. Retail Price € 121

POWER SETTING MEDIUM

100 Universal.
Probe-sized tip used to remove light deposits in shallow to moderate pockets and for periodontal maintenance.

US100

Sugg. Retail Price € 137

POWER SETTING MEDIUM

P10 Universal.
Standard diameter tip used to remove light deposits and plaque. Tip features smooth edges for safe removal of soft and hard deposits.

USP10

Sugg. Retail Price € 121

POWER SETTING LOW

3 Perio
Thin-bladed curette-like tip.

USP3

Sugg. Retail Price € 143

ULTRASONIC TIPS

SERIE E

BUY 2 AND GET 1 FOR FREE

Compatible with all EMS* generators

TOP-SELLER

POWER SETTING MEDIUM

S10 Universal.
Standard diameter tip used to remove heavy supragingival deposits.

UES10

Sugg. Retail Price € 127

POWER SETTING MEDIUM

P10 Universal.
Standard diameter tip used to remove moderate to heavy deposits located supragingivally and in shallow pockets.

UEP

Sugg. Retail Price € 127

POWER SETTING LOW

100 Universal.
Probe-sized tip used to remove light deposits in shallow to moderate pockets and for periodontal maintenance.

UE100

Sugg. Retail Price € 152

POWER SETTING LOW


3 Perio
Thin-bladed curette-like tip.

UEP3

Sugg. Retail Price € 146



* Our offered items are compatible with but not original products of device manufacturing companies Satelec and/or EMS and are not produced, licensed or recommended from them. All used brand names and definitions are registered trademarks and brands of Acteon, Inc. and Electro Medical Systems respectively. They are intended only to show compatibility of items with products of different manufacturers.



ORTHODONTIC LINE

VALID FOR ALL PART CODES OF ORTHO STANDARD PLIERS & CUTTERS



DID YOU KNOW?

All hinged instruments – forceps, rongeurs, scissors, pliers, hemostats, orthodontic pliers, etc. – should be kept lubricated. Regular use of proper lubricants, like Hu-Friedy Penetrating Oil (IPS) and Instrument Lubricant Spray (ILS), will prevent rust, corrosion and stiff joints to ensure smooth operation. All hinged instruments should be sterilised in the open position. [see page 3](#)



CUTTERS

With our diamond-honed cutting edges and orbit-formed joint, you will experience sharp and accurately-aligned cutting edges for every procedure, every time.

For each individual area of focus, we provide a specifically engineered cutter. To meet these stringent standards, all of our blades are hand-finished to maximise cutting performance for the appropriate wire size.

Universal Cut & Hold Distal End Cutter, Long Handle, wire size <= .012" to .021"x.025" (.30 - .53x.64 mm) **678-101L**



UTILITY PLIERS

Strong, yet lightweight and perfectly balanced, the dimensions of these instruments are perfectly proportioned, they fit your grasp with unsurpassed ease.

Designed to provide superior control with effortless versatility, Hu-Friedy's Utility Pliers assure comfort to both you and your patients.

Slim Weingart Plier, wire size <= .030" (.76 mm) **678-202**



WIRE FORMING PLIERS

These instruments perform as natural extensions of your professional artistry.

Our impeccably smooth finished tips allow you to bend the wire without risk of damage. Whether you are shaping stainless steel or nickel titanium, Hu-Friedy has a hand-finished, precision instrument to help bring your brilliant workmanship to life.


Tweed Loop Forming Plier, wire size <= .022" x .025" (.56 x .64 mm) **678-306**

20% OFF FROM RETAIL PRICES

Sugg. Retail Price € 200 **Promo Price € 160**

Sugg. Retail Price € 175 **Promo Price € 140**

Sugg. Retail Price € 175 **Promo Price € 140**



NEEDLE HOLDERS



NH5024SIM
Simion Straight Castroviejo Perma Sharp, 3-0, 4-0* 18 cm

Sugg. Retail Price € 481

Promo Price € 336



NH5038
Crile-Wood Perma Sharp, 3-0, 4-0, 5-0, 6-0*, 15 cm

Sugg. Retail Price € 166

Promo Price € 116



NHCV
Castroviejo Microsurgical, 5-0, 6-0, 7-0, 8-0*, 14 cm

Sugg. Retail Price € 299


Promo Price € 209



NHCW
Crile-Wood, grooved, 3-0 - 6-0*, 15 cm

Sugg. Retail Price € 49

Promo Price € 34



NHMH
Crile-Wood, grooved, 1-0, 2-0, 3-0, 4-0*, 16 cm

Sugg. Retail Price € 55

Promo Price € 38


* Recommended suture sizes.



SCISSORS

SUPER-CUT DESIGN

Featuring a razor sharp blade and serrated edge for reduced tissue slippage.



S14SC
#14 LaGrange Double-Curved Super-Cut, 11.5 cm, 1 serrated blade

Sugg. Retail Price	Promo Price
€ 124	€ 87



S16
#16 Goldman-Fox Curved Scissors, 12.5 cm

Sugg. Retail Price	Promo Price
€ 55	€ 38




SCORTSC
Cortellini, Curved Super-Cut, 12 cm

Sugg. Retail Price	Promo Price
€ 120	€ 84

DID YOU KNOW?

All hinged instruments – forceps, rongeurs, scissors, pliers, hemostats, orthodontic pliers, etc. – should be kept lubricated. Regular use of proper lubricants, like Hu-Friedy Penetrating Oil (IPS) and Instrument Lubricant Spray (ILS), will prevent rust, corrosion and stiff joints to ensure smooth operation. All hinged instruments should be sterilised in the open position. [see page 3](#)



SURGICAL CURETTES



CL84
#84 Lucas Spoon Shape, handle #522 Hexagon



CL85
#85 Lucas Spoon Shape, handle #522 Hexagon



CL86
#86 Lucas Spoon Shape, handle #522 Hexagon



CL87
#87 Lucas Spoon Shape, handle #522 Hexagon



CM2/4
#2/4 Molt, handle #522 Hexagon

Sugg. Retail Price (per instrument)	Promo Price (per instrument)
€ 94	€ 61

ELEVATORS




EHB13/14
#13/14 Heidbrink Root Tip Pick, handle #522 Hexagon



EHO2/3
#2/3 Howard Root Tip Pick, handle #522 Hexagon

Sugg. Retail Price (per instrument)	Promo Price (per instrument)
€ 97	€ 72



EL3C
3 mm Curved Luxating Elevator, handle #510



EL3S
3 mm Straight Luxating Elevator, handle #510



EL5C
5 mm Curved Luxating Elevator, handle #510



EL5S
5 mm Straight Luxating Elevator, handle #510

Sugg. Retail Price (per instrument)	Promo Price (per instrument)
€ 113	€ 85

PERIOSTEALS



P24G
#24G, handle #2

Sugg. Retail Price	Promo Price
€ 88	€ 70



PPR3S
#3 Prichard, handle #522 Hexagon

Sugg. Retail Price	Promo Price
€ 104	€ 83




PPBUSER6
Buser, handle #6


Sugg. Retail Price	Promo Price
€ 100	€ 80



P146
#14 Goldman-Fox, handle #6

Sugg. Retail Price	Promo Price
€ 89	€ 72





TUNNELING KNIVES



TKN1
#1, handle #6

Sugg. Retail Price € 116
Promo Price € 81



TKN2
#2, handle #6

Sugg. Retail Price € 116
Promo Price € 81



PPAEL
Allen, handle #6

Sugg. Retail Price € 120
Promo Price € 84

TISSUE PLIERS



TPG1
Serrated/Straight
Gerald Plain, 18 cm

Sugg. Retail Price € 52
Promo Price € 42



TPG2
1x2 Curved/Anatomical
Gerald, 18 cm

Sugg. Retail Price € 52
Promo Price € 42



TP41
#41 Plain Adson,
12 cm

Sugg. Retail Price € 33
Promo Price € 26



TP43
#43 7x7 Adson-Brown,
12 cm

Sugg. Retail Price € 64
Promo Price € 51

SCALPEL BLADES AND HANDLES



10-255-12
Scalpel Blades #12,
100pcs./pack., Expiry Date

Sugg. Retail Price € 73
Promo Price € 47



10-255-15
Scalpel Blades #15,
100pcs./pack., Expiry Date

Sugg. Retail Price € 73
Promo Price € 47



10-256-12
Scalpel Blades #12D,
100pcs./pack., Expiry Date

Sugg. Retail Price € 88
Promo Price € 57




10-256-15
Scalpel Blades #15C,
Box of 100, Carbon Steel

Sugg. Retail Price € 88
Promo Price € 57



10-130-05E
No. 5E European Style Round Scalpel Handle, Matte Finish

Sugg. Retail Price € 87
Promo Price € 61





XTS INSTRUMENTS

OFFER VALID FOR ALL XTS PART CODES



TNCVIPC
Interproximal Carver XTS Composite Instrument DE
Sugg. Retail Price € 55
Promo Price € 39



TNCIGFT2
#2 Goldstein Flexi-Thin XTS Composite Instrument DE
Sugg. Retail Price € 55
Promo Price € 39



#3 Extra-Flex Goldstein Flexi-Thin XTS Composite Instrument, Anterior Kit, DE
TNCIGFT3
Sugg. Retail Price € 55
Promo Price € 39



TNPFIAG
#A6 XTS Composite Instrument DE
Sugg. Retail Price € 55
Promo Price € 39



TNPFI8A
#8A XTS Composite Instrument DE
Sugg. Retail Price € 55
Promo Price € 39



TNCIGFT1
#1 Goldstein Flexi-Thin XTS Composite Instrument DE
Sugg. Retail Price € 55
Promo Price € 39



EXTENDING YOUR INSTRUMENTS LIFE

Our newly upgraded **technical care center**, based in our facilities in Tuttlingen, further increases our goal of **maximising the useful life of our instruments**. This post-sale service allows dental professionals to ensure their instruments always perform at their best and for the longest time possible.

TECHNICAL CARE CENTER

WHAT CAN WE DO FOR YOU?

Your dedicated craftsmanship provides you with a full set of services to extend instruments life performance.

SHARPENING



RE-SERRATING

JOINT REALIGNMENT

CARBIDE INSERT REPLACEMENT (UP TO THREE REPLACEMENTS)



HOW IT WORKS?



CONTACTS

Riccardo Boschetto

Junior Territory Manager Eastern Europe,
Via Mauro Macchi, 27 - 20124 Milan, Italy
Tel. +39 02 36 58 96 00 ext. 202
Mob. +39 340 885 96 43
Mail rboschetto@hu-friedy.com



DOWNLOAD OUR CATALOGUE

For all of us at Hu-Friedy, this is much more than a book that contains the highest quality dental instruments and products. It's a manifestation of the dedication, innovative spirit and hard work of a company that has served and supported the dental community for over 110 years.

FOLLOW US ON:



@HuFriedyEU



@hu_friedy_europe

For more information about
Hu-Friedy's Total Solutions visit

HU-FRIEDY.EU


PREVENTIVE




RESTORATIVE



ORTHODONTIC




INSTRUMENT
MANAGEMENT



ENDODONTIC



SURGICAL



POWER



Chicago | Milan | Frankfurt | Shanghai | Tokyo | Tuttlingen

Manufacturer: Hu-Friedy Mfg. Co., LLC | 3232 N. Rockwell St. | Chicago, IL 60618 | USA

©2019 Hu-Friedy Mfg. Co., LLC. All rights reserved. The contents of this promotions including, but not limited to, all text, photographs and graphics, are the sole and exclusive property of Hu-Friedy Mfg. Co., LLC and may not be copied, reproduced or distributed without prior written permission. Hu-Friedy Mfg. Co., LLC. reserves the right to change features of products, images and text at any time without notice. SporeClear is a registered trademark of Altos Medical, LLC. HFL-072/0919

How the best perform

