

LC-Primer

Manufacturer:

SCHEU-DENTAL GmbH
Am Burgberg 20
D-58642 Iserlohn
Tel: +49 2374 9288-0
Fax: +49 2374 9288-90
URL: www.scheu-dental.com
E-Mail: service@scheu-dental.com



Description:

LC- Primer is a light-curing primer for a chemical compound between PET-G pressure moulding material and DURASPLINT® LC and CLEAR-BLOKKER®. It can also be used for the compound of already polymerised DURASPLINT® LC with CLEAR-BLOKKER® in case of build-ups or corrections.

Composition:

Methyl methacrylate, acrylic acid ester, photoinitiator

Curing:

5 minutes in light-curing oven (LC-6 Light oven) equipped with 3 UV-A tubes and 3 VIS tubes (wavelength range 350 – 450 nm), alternating order of tubes.

Medical product:

CE

Packaging unit:

10 ml

Storage condition:

4 – 25°C, not to be exposed to direct sunlight.

