



Item **Description**
TN003 **Medium Blade with Large Plugger**

- Contouring large Class I, II, V
- Packing composite material into large preparations
- Convex plugger tip minimizes pullback of composite material

TN003



TN004 **Extra Thin Blades**

- Reach under margins to hone and sculpt
- Sculpting interproximally
- Packing cord in gingival crevice

TN004



TN005 **Long Extra Thin Blades**

- Allows additional access for hard-to-reach
- Posterior contouring and shaping
- Sculpting interproximally
- Packing cord in gingival crevice

TN005



TN006 **Small/Medium
Multi-Function Plugger**

- Placement of composite material in smaller Class II and III
- Packing material toward gingival seats
- Carving occlusal anatomy with acorn carver

TN006

