

Composite Instruments Nonstick titanium-nitride coated

Our multi-function instruments have double condenser configurations and carving blades on opposite ends. Now you can place, pack, shape and carve your way to perfectly contoured composite restorations without ever having to switch instruments!

- Large-diameter, lightweight handles for improved ergonomics.
- Provides non-stick surface for placement and shaping of composites.
- 40% harder than stainless steel for a non-porous, non-stick surface.
- Withstands diverse cleaning/sterilization methods, including ultrasonic cleaning.



TN009,
Louis Mackenzie,
UK



TN008

5-in-1

Item
TN008

Description
5-in-1 Multi-Function Instrument
2 condensers, ball burnisher, blade, acorn



TN009

5-in-1

Item
TN009

Description
5-in-1 Multi-Function Instrument
2 blades, 2 condensers, acorn