

# Composi-Tight® 3D Fusion™ Firm Matrix Bands

**NEW!**

New **Composi-Tight® 3D Fusion™** Firm matrix bands with **eZ-Place™** technology maintain their shape during placement making them *significantly easier* to insert!

eZ-Place™ technology combines specially processed stainless steel with a proprietary matrix band forming process to yield Fusion Firm's remarkable handling characteristics. It resists deformation nearly 70% better than typical dead-soft bands yet is pre-contoured and super thin at 0.0015". You're ensured the tight contacts and ideal anatomy you've come to expect from the Composi-Tight® 3D Fusion™ system.

**Fusion Firm:**

- Nearly 70% better at resisting deformation!
- Super-thin at only 0.0015"
- Five ideal sizes!

**Perfect for:**

- Tight interproximal spaces!
- More conservative preps!
- First-time users and veterans alike!

**Composi-Tight® 3D Fusion™ Firm matrix bands are now available in a Composi-Tight® 3D Fusion™ kit.**

FXR-HHF-10



0.0015" thick

## Composi-Tight® 3D Fusion™ Firm Bands

Description	Item	Qty.
Assorted Band Kit: 5 sizes	FXHB04	300
Assorted Band Mini Kit: 5 sizes	FXHB05	150

FXH Series	Composi-Tight® 3D Fusion™ Firm matrix bands with eZ-Place™ technology for easier placement.					
	NA	FXH100 4.3 mm	FXH150 4.3 mm w/ext 5.9 mm	FXH175 5.3 mm	FXH200 6.1 mm	FXH300 6.1 mm w/ext 7.7 mm
thickness 0.0015"						
height						

Item	Standard (-M)		Bulk	
	Qty.		Qty.	
FXH100	50		100	
FXH150	50		100	
FXH175	50		100	
FXH200	50		100	
FXH300	50		100	

Be sure to check out the website at [www.garrisdental.com](http://www.garrisdental.com) for info on other innovative Garrison products, demo videos and our 6-month money-back guarantee.

**Garrison**  
Dental Solutions

150 DeWitt Lane  
Spring Lake, MI USA  
49456

To order: Call toll free **888.437.0032**  
or online at [www.garrisdental.com](http://www.garrisdental.com)




Mention Code  
**IN120FXH**



# Garrison BAND COMPARISON CHART

Dental Solutions



FX Series  
SXR Series  
M – SM – AU Series  
B Series

	Pediatric	Bicuspid	Bicuspid with deep cervical preps	Bicuspid and small molars	Molars	Molars with deep cervical preps			
FXH Series	Composi-Tight® 3D Fusion™ Firm matrix bands with eZ-Place™ technology for easier placement.							Standard (-M)	Bulk
	thickness 0.0015"						Item	Qty. Price	Qty. Price
height	NA	FXH100 4.3 mm	FXH150 4.3 mm w/ext 5.9 mm	FXH175 5.3 mm	FXH200 6.1 mm	FXH300 6.1 mm w/ext 7.7 mm	FXH100	50	100
							FXH150	50	100
							FXH175	50	100
							FXH200	50	100
							FXH300	50	100
FX Series	Composi-Tight® 3D Fusion™ Full Curve Bands have non-stick coating and an occlusal placement tab.							Standard (-M)	Bulk
	thickness 0.0016"						Item	Qty. Price	Qty. Price
height	NA	FX100 4.0 mm w/ext 4.4 mm	FX150 4.1 mm w/ext 6.0 mm	FX175 5.0 mm w/ext 5.6 mm	FX200 5.4 mm w/ext 6.6 mm	FX300 6.0 mm w/ext 8.7 mm	FX100	50	100
							FX150	30	60
							FX175	50	100
							FX200	50	100
							FX300	30	60
NTG Series	NiTin™ full contour matrix bands for ideal anatomy.								Bulk
	thickness 0.0015"						Item		Qty. Pk.
height	NA	NTG100 4.0 mm w/ext 4.4 mm	NTG150 4.1 mm w/ext 6.0 mm	NTG175 5.0 mm w/ext 5.6 mm	NTG200 5.4 mm w/ext 6.6 mm	NTG300 6.0 mm w/ext 8.7 mm	NTG100		50
							NTG150		50
							NTG175		50
							NTG200		50
							NTG300		50
SXR Series	Slick Bands™ XR Matrices have non-stick coating and Grab-Tabs™ for easier placement.							Standard (-M)	Bulk
	thickness 0.0016"						Item	Qty. Pk.	Qty. Pk.
height	NA	SXR100 4.6 mm	SXR150 3.8 mm w/ext 5.4 mm	SXR175 5.4 mm	SXR200 6.4 mm	SXR300 6.4 mm w/ext 8.6 mm	SXR100	50	100
							SXR150	30	60
							SXR175	50	100
							SXR200	50	100
							SXR300	30	60
SM Series	Slick Bands™ Matrices have non-stick coating.							Standard (-M)	Bulk
	thickness 0.0016"						Item	Qty. Price	Qty. Price
height		SM050 3.2 mm w/ext 4.2 mm	SM100 4.6 mm	SM150 3.8 mm w/ext 5.4 mm	SM175 5.5 mm	SM200 6.4 mm	SM050	30	60
							SM100	50	100
							SM150	30	60
							SM175	50	100
							SM200	50	100
							SM300	30	60
M Series	Composi-Tight® bands are contoured in 3 dimensions to provide an accurate contact at the tooth's natural height of contour.							Standard (-M)	Bulk
	thickness 0.0015"						Item	Qty. Pk.	Qty. Pk.
height		M050 3.2 mm w/ext 4.2 mm	M100 4.6 mm	M150 3.8 mm w/ext 5.4 mm	M175 5.5 mm	M200 6.4 mm	M050	25	50
							M100	50	100
							M150	25	50
							M175	50	100
							M200	50	100
							M300	25	50
AU Series	The Gold bands are the same as the Composi-Tight® bands.							Standard (-M)	Bulk
	thickness 0.0015"						Item	Qty. Pk.	Qty. Pk.
height		AU050 3.2 mm w/ext 4.2 mm	AU100 4.6 mm	AU150 3.8 mm w/ext 5.4 mm	NA	AU200 6.4 mm	AU050	25	50
							AU100	50	100
							AU150	25	50
							AU200	50	100
							AU300	25	50
B Series	The original bands have a slightly shorter length.							Standard (-M)	Bulk
	thickness 0.0015"						Item	Qty. Pk.	Qty. Pk.
height		B050 3.2 mm w/ext 4.2 mm	B100 4.6 mm	B150 3.8 mm w/ext 5.4 mm	NA	B200 6.4 mm	B050	25	50
							B100	50	100
							B150	25	50
							B200	50	100
							B300	25	50
MC Series	Clear bands are transparent and translucent to allow for trans-enamel polymerization without interference from metal bands.							Standard (-M)	Bulk
	thickness 0.0020"						Item	Qty. Pk.	Qty. Pk.
height	NA	MC100 5.0 mm	NA	MC175 5.9 mm	MC200 6.8 mm	NA	MC100	50	100
							MC175	50	100
							MC200	50	100

To order standard packs, add a -M after the item #.